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

### ENGLISH

#### 1- Charging

- **Full charging time:** approximately 2.5 hours.
  - **Option A:** Plug the USB-C charger cord into the charging port at the bottom of the pipette handle and connect the other end of the charger to a wall power socket. Use this option to continue working while charging.
  - **Option B:** Charge the pipette via inductive charging with the charging base. Connect the USB-C charger cord into the base of the pipette stand, then place the pipette in the charging position for wireless power charging.
- Charging is complete when all bars are displayed simultaneously on the LCD screen.
- To maintain optimal battery health and prolong its lifespan, it is recommended to recharge the pipette's battery when the power indicator shows only one bar remaining.
- The pipette can be used while charging if needed.

#### 2- Battery and Charging Features

- **Charging temperature:** 10°C to 35°C.
- **Charger requirements:** Input AC 100–240V, Output 5V, 1000mA.
- The battery level is displayed on the LCD screen.
- The pipette uses a rechargeable lithium battery: 3.7V, 3000mAh.
- **When storing the device for extended periods, please disconnect the battery.**
- Built-in protection prevents overcharging, allowing you to recharge it at any time without risking damage to the battery.
- Stop using the battery if abnormal heat, odor, discoloration, deformation, or abnormal condition is detected during use, charge, or storage. Dispose of the battery immediately according to rules in your country.

#### 3- Pipette Control Features

This pipette offers intuitive controls for precise liquid handling:

##### 1. Operation Buttons

- The handle has two main operation buttons on the front of the handle:
  - **The top button is for liquid suction.**
  - **The bottom button is for liquid discharge.**

You can precisely control both the suction and discharge speeds by simply varying the pressure you apply to these buttons.

**2. Speed Regulation Button:** Located just below the LCD screen, the speed regulation button allows you to pre-set your desired liquid suction and dispensing speeds.

**3. Pipette Adapter:** The pipette adapter assembly includes both the nosecone and silicone pipette holder. It is compatible with both plastic and glass pipettes (including Pasteur pipettes) ranging in volume from 0.1 ml to 100 ml.

**4. Pipette Filter:** When removing the pipette adapter, you'll access the 0.22  $\mu$ m PTFE hydrophobic filter. This filter is a crucial component to help protect the pipette.

##### Important Note:

- **Never use the pipette without a filter.**
- If the filter gets wet from drawing in too much liquid, it will significantly slow down or completely block liquid flow. In such cases, you must replace the filter immediately.

#### 4- Operating Your Pipette

**1. Adjusting Speed :** You can set your desired liquid suction and dispensing speeds using the speed control button located below the LCD screen. Both suction and dispensing speeds have five levels. Use the " $\Delta$ " button to increase speed and the " $\nabla$ " button to decrease it. Your selected speed will be shown at the bottom of the LCD screen. You can also use gravity for liquid discharge.

**Note:** While operating, you can also adjust the dispensing speed by varying the pressure you apply to the operation button with your fingers.

**2. Inserting the Pipette:** When inserting the pipette, grip it as close to the top as possible. Then, carefully push it into the pipette holder until you feel resistance.

##### 3. Aspirating Liquid

- a. Select the appropriate speed.
- b. Press the liquid suction button. Draw in slightly more liquid than your target volume.
- c. Lift the pipette off the liquid surface in the container.
- d. Press the liquid dispensing button to release some liquid until the meniscus (the curved surface of the liquid) in the suction tube reaches your desired volume scale.

**4. Dispensing Liquid:** The pipette offers two ways to dispense liquid:

- **Gravity Discharge:** Gently press the liquid discharge button until you feel a slight click. At this point, the liquid will flow out naturally by gravity (this works when the motor speed is zero).
- **Accelerated Discharge:** Press the liquid dispensing button down further. This activates the micro-pump, which significantly speeds up the liquid discharge.

**Note:** Regardless of the method you choose, you can finely adjust the dispensing speed by varying the pressure you apply with your finger.

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