

Full range of
PROGENE[®]
SYRINGE FILTERS



Certified reliability ■ Assured performance ■ 100% quality control ■ ISO certification

Innovative Products for Today's Modern Laboratory

UltiDent
Scientific

T : : 514 335-3433

Toll free : : 1 800 463-4163

F : : 514 335-0992

ultident.com



- Choose from sterile or non-sterile formats
- Filters have a female Luer-Lock inlet and a male Luer-Slip outlet
- 100% of production is tested by Quality Control on a lot by lot basis
- 0.20 µm pore size is ideal for sterilization and fluid ultrafiltration; 0.45 µm for fluid purification
- Manufactured under ISO 9001:2000 and ISO 13485:2003 Certification in a Class 10,000 Clean Room
- **CA : Hydrophilic** - Low protein binding - Ideal for cell culture and biological research
- **PES : Hydrophilic** - Fast flow rate and better chemical resistance than CA - Ideal for cell culture and biological research
- **Nylon : Hydrophilic** - High protein binding - Ideal for filtration of organic and aqueous solutions
- **PTFE : Hydrophobic** - Unaffected by most organic solvents, concentrated acids, bases and cryogenic fluids

CAT. #	PORE SIZE	DIAMETER	PACKAGING	PRICE
CA - Cellulose Acetate Membrane, Sterile				
22-CA-2520	0.20 µm	25 mm	1/Blister Pack, 50/Box	54.95
22-CA-2545	0.45 µm	25 mm	1/Blister Pack, 50/Box	54.95
PES - Polyethersulfone Membrane, Sterile				
22-PES-1320	0.20 µm	13 mm	1/Blister Pack, 50/Box	64.95
22-PES-1345	0.45 µm	13 mm	1/Blister Pack, 50/Box	64.95
22-PES-2520	0.20 µm	25 mm	1/Blister Pack, 50/Box	64.95
22-PES-2545	0.45 µm	25 mm	1/Blister Pack, 50/Box	64.95
Nylon Membrane, Non-Sterile				
22-NY-1320	0.20 µm	13 mm	Bulk, 100/Box	99.95
22-NY-1345	0.45 µm	13 mm	Bulk, 100/Box	99.95
22-NY-2520	0.20 µm	25 mm	Bulk, 100/Box	99.95
22-NY-2545	0.45 µm	25 mm	Bulk, 100/Box	99.95
PTFE Membrane, Non-Sterile				
22-PTFE-1320	0.20 µm	13 mm	Bulk, 100/Box	109.95
22-PTFE-1345	0.45 µm	13 mm	Bulk, 100/Box	109.95
22-PTFE-2520	0.20 µm	25 mm	Bulk, 100/Box	109.95
22-PTFE-2545	0.45 µm	25 mm	Bulk, 100/Box	109.95

Innovative Products for Today's Modern Laboratory