



### PROGENE® Low Retention Pipette Tips

- Filter tips protect your samples from all possible contamination. They ensure sample integrity and do not impede air flow.
- 100% reliability Our low-retention filter tips are hydrophobic and do not contain any additives.
- Certified RNase, DNase, DNA and PCR inhibitor free.
- Certified non-pyrogenic.
- Sterilized by gamma irradiation.







Our Tip Refill System saves time & money!

#### LOW RETENTION TIP REFILL SYSTEM

PROGENE® low-retention pipette tips are engineered with a specially formulated resin, providing complete sample control. The inner surfaces of our low-retention tips are considerably more hydrophobic than standard industry tips enabling viscous fluids to be dispensed more completely. There is less fluid retention and less sample loss. PROGENE® tips prove that less is more when it comes to liquid handling!

#### **CLEAN LIQUID RELEASE**

PROGENE® tips are purposely designed for a much cleaner liquid release. Our thinner walls and more finely tapered tip construction eliminate the need to touch off. Users will maintain complete fluid control working with smaller volumes, and not compromise their pipetting accuracy or precision.

The PROGENE® Refill System provides handy lift-and-load simplicity. Our new system does not require any tools or special instructions to move from an empty insert to a fresh set of 96 tips.

The PROGENE® Tip Refill System reduces lab waste by 70% and require less shelf space than standard racks.

- Our tip rack is designed to be a time saver as it allows for onehanded operation. A simple press on the front of the lid releases it from the rack; otherwise the lid sits securely closed.
- The PROGENE® Tip Refill System saves precious time, storage space and plastic waste.
- PROGENE® tip racks are robust and designed to withstand multiple autoclave cycles.
- Our tip rack was designed to facilitate off set pipetting. The edge of the rack and insert are out of the way so users can easily load tips onto a multichannel pipette.

# Buy 4 Get 1 FREE!

\*equal or lesser value

## PROGENE® Pipette Tips

#### PROGENE® Research-Grade Pipette Tips

10μL P2/P10 Filter	Part Number	Description	Quantity	List Price	20% Promo
	87-A10	10µL bulk	1000/Bag	16.95	13.56
10μL P2/P10	87-A10R	10µL racked	10 hinged racks of 96	49.95	39.96
	87-A10XL	10µL bulk	1000/Bag	16.95	13.56
	87-A10XLR	10µL racked	10 hinged racks of 96	49.95	39.96
	87-A200C	200μL, clear, bulk	1000/Bag	13.95	11.16
10μL XL Filter	87-A200Y	200μL, yellow, bulk	1000/Bag	13.95	11.16
	87-A200CR	200μL, clear, racked	10 hinged racks of 96	49.95	39.96
10μL XL	87-A200YR	200μL, yellow, racked	10 hinged racks of 96	49.95	39.96
	87-A300	300µL bulk	1000/Bag	24.95	19.96
	87-A300R	300µL racked	10 hinged racks of 96	53.95	43.16
200µL	87-A1000B	1000µL bulk	1000/Bag	20.95	16.76
200μL	87-A1000BR	1000µL racked	10 hinged racks of 100	67.95	54.36
	87-A1250	1250µL bulk	1000/Bag	22.95	18.36
200µL Filter	87-A1250-R	1250µL racked	10 hinged racks of 96	62.95	50.36

#### PROGENE® Tip Refill System

Part Number	Description	Quantity	List Price	20% Promo
87-A10-RFL	10µL refill, Low-Retention	10 refill inserts of 96	44.95	35.96
87-A10XL-RFL	10μL refill, Low-Retention	10 refill inserts of 96	44.95	35.96
87-A200C-RFL	200μL refill, Low-Retention	10 refill inserts of 96	35.95	28.76
87-A300-RFL	300μL refill, Low-Retention	10 refill inserts of 96	49.95	39.96
87- A1250-RFL	1250µL refill	10 refill inserts of 96	61.95	49.56

#### PROGENE® Low Retention Filter Tips

Part Number	Description	Quantity	List Price	20% Promo
87-F10S-LR	10µL graduated	10 hinged racks of 96	77.95	62.36
87-F10XLS-LR	10µL graduated reach	10 hinged racks of 96	77.95	62.36
87-F30S-LR	20µL graduated	10 hinged racks of 96	77.95	62.36
87-F100S-LR	100µL graduated	10 hinged racks of 96	77.95	62.36
87-F200S-LR	200µL graduated	10 hinged racks of 96	77.95	62.36
87-F300S-LR	300µL graduated	10 hinged racks of 96	77.95	62.36
87-F1000S-LR	1000µL graduated	10 hinged racks of 100	82.95	66.36
87-F1250S-LR	1250µL graduated reach	10 hinged racks of 96	82.95	66.36

1000μL Filter

1000µL

100µL Filter

20µL Filter

300µL

300µL Filter

#### PROGENE® Low-Retention Gel Loading Tips

- Certified RNase, DNase, DNA, and PCR inhibitor-free.
- Racked format is packaged in multichannel-friendly 96-tip racks.
- Racks feature a patented lid-lock design for easy, one-hand operation.
- Low-retention resin reduces sample loss and enhances sample delivery.

Part Number	Description	Quantity	List Price	20% Promo
87-A200-GEL	200 μL round gel loaders, bulk	1000/Bag	58.95	47.16
87-A200-GEL-R	200 μL wround gel loaders, racked	2 racks of 96	21.95	17.56

1250µL Filter

1250µL




#### **PCR Plates**

# Buy 4 Get 1 FREE!

\*equal or lesser value

Our large selection of 96 and 384-well plates are designed to precisely fit the majority of PCR thermal cyclers and sequencers on the market today. Each is designed to meet and exceed the specifications and performance of plates supplied by the original equipment manufacturer.

- Well ridges reduce sample evaporation when sealing with film, foil or associated strip caps.
- High contrast black alphanumerics for easy well identification.
- Ultra-thin, uniform wells ensure optimal heat transfer for high reaction efficiency.
- Conveniently packaged in small quantities of 10 plates to minimize risk of contamination.
- Certified free from DNase, RNase, DNA, and PCR inhibitors.
- Made in the USA



Part Number	Description	Volume	Quantity	List Price	20% Promo
87-C96-NS	96-well plate, Non-skirted	200µL	10/Box, 50/Case	144.95	115.96
87-C96-LPNS	96-well plate, Non-skirted, Low-Profile	100µL	20/Box	67.95	54.36
87-C96-ABI-1	96-well plate, Half-skirted, ABI® Hybrid	200µL	10/Box, 50/Case	144.95	115.96
87-C96-ABI-2	96-well plate, Half-skirted, ABI® Standard	200µL	10/Box, 50/Case	144.95	115.96
87-C96-RR	96-well plate, Half-skirted, Raised-rim	200µL	10/Box, 50/Case	185.95	148.76
87-C96-FAST	96-well plate, Half-skirted, ABI® Fast PCR	100µL	10/Box, 50/Case	147.95	118.36
87-C96-SK	96-well plate, Full-Skirted, Low-profile	200µL	10/Box, 50/Case	142.95	114.36
87-C96-LC480-W	96-well plate, LightCycler® Low-profile,	200µL	10/Box	65.95	52.76
	white (Kit with qPCR film)				
87-C384-1N	384-well plate, 1-notch (A24), Full-Skirted	50μL	10/Box, 50/Case	259.95	207.96
87-C384-2N	384-well plate, 2-notch (A24/P24), Full-Skirted	50μL	10/Box, 50/Case	259.95	207,96
87-C384-LC480-W	384-well plate, LightCycler® Low-profile,	50μL	10/Box	84.95	67.96
	white (Kit with qPCR film)	0	6	2	



\*equal or lesser value

Our PCR tubes are engineered to precisely fit most thermal cyclers on the market today. Strip tubes with individually-attached caps are customizable, allowing users to easily split the strip into smaller sections. Strip tubes with white wells enhance fluorescence signal strength, making them ideal for qPCR analysis when combined with an optically-clear flat cap. We offer many standard and Low Profile versions that are all conveniently packaged in resealable bags.

- Secure, evaporation-resistant seal between cap and tube.
- Ultra-thin walls for optimum heat transfer and reaction efficiency.
- Optically-clear flat caps ideal for qPCR and standard PCR.
- Crystal-clear clarity easy to visualize content of each well.
- Available in full range of colors or white (enhances qPCR signals).
- Certified free from DNase, RNase, DNA, and PCR inhibitors.
- Made in the USA.





Part Number	Description	Quantity	List Price	20% Promo
Strip Tubes				
87-C200-8-TB	0.2 ml 8-strip tubes	125/Pack	73.95	59.16
87-C100-8-TB	0.1 ml 8-strip tubes	125/Pack	73.95	59.16
87-C200-12-TB	0.2 ml 12-strip tubes	80/Pack	104.95	83.96
Strip Caps				
87-C8-QPCR	8-strip flat caps	125/Pack	26.95	21.56
87-C8-DOME	8-strip domed caps	125/Pack	21.95	17.56
87-C12-QPCR	12-strip flat caps	80/Pack	24.95	19.96
87-C12-DOME	12-strip domed caps	80/Pack	24.95	19.96
Strip Tubes with Attached	Caps			
87-C200-8-IND	0.2 ml 8-strip tubes with individual-attached	120 strips of 8	88.95	71.16
	optically flat caps			
Rotor-Gene®				
Part Number	Description	Quantity	List Price	20% Promo
87-C100-4-RG	4-strip 0.1 ml Rotor-Gene® Tubes & Caps	250/Pack	138.95	111.96

### Individual PCR Tubes

Part Number	Description	Quantity	List Price	20% Promo
87-C200-D	0.2 ml domed cap	1000/Pack	52.95	42.36
87-C200-F	0.2 ml frosted flat-cap	1000/Pack	44.95	35.96
87-C500-F	0.5 ml frosted flat-cap	1000/Pack	40.95	32.76

### Sealing Films & Foils

- Sealing films feature removable end tabs for easy handling and application.
- Certified free from DNase, RNase, DNA, and PCR inhibitors.
- Made in the USA.

Part Number	Description	Quantity	List Price	20% Promo
87-CFILM-QPCR	Pressure-Sensitive qPCR sealing film	100 films	202.95	162.36
87-CFILM-AL	Pierceable Aluminum Foil	100 films	93.95	74.99
87-CFILM-TS	Adhesive routine PCR sealing film	100 films	76.95	61.56




# 96- and 384-Well Plates

LightCycler® compatible

- Rapid, even heat transfer
- Real-time compatible
- Certified RNase, DNase, DNA, and PCR inhibitor free
- Tested pyrogen-free
- Contact us for free samples!

NON-SKIRTED, STANDARD 0.2mL Rimmed wells		87-C96-NS	A12 corner notch, natural	10/Box, 50/Case
NON-SKIRTED, LOW-PROFILE 0.1mL Raised rims		87-C96-LPNS	A12 corner notch, natural	20/Box
HALF-SKIRTED, ABI-TYPE 0.2mL Fits most cyclers, including ABI Prism® and GeneAmp® Systems.		87-C96-ABI-1	A12 corner notch, natural	10/Box, 50/Case
HALF-SKIRTED, RAISED RIMS 0.2mL		87-C96-RR	A12 corner notch, natural	10/Box, 50/Case
HALF-SKIRTED, LOW-PROFILE 0.1mL LightCycler® compatible		87-C96-LC480-W	H12 corner notch, white	10/Box, 50/Case
HALF-SKIRTED, LOW-PROFILE RAISED RIMS 0.1mL Fits ABI Fast® cyclers	**************************************	87-C96-FAST	A1 corner notch, natural	10/Box, 50/Case
HALF-SKIRTED, STANDARD 0.2mL		87-C96-ABI-2	A12 corner notch, natural	10/Box, 50/Case
FULL-SKIRTED, LOW-PROFILE RAISED RIMS 0.1mL Stackable Rimmed wells		87-C96-SK	A12 corner notch, natural	10/Box, 50/Case
FULL-SKIRTED, RAISED RIMS, 384-WELL 0.2mL Conical flat-bottom wells		87-C384-2N	A24, P24 corner notch, natural	10/Box, 50/Case



# Microcentrifuge Tubes

\*equal or lesser value

1.5

Microcentrifuge tubes are available in 0.5ml, 1.5ml, and 2.0ml volumes and are solidly constructed for high speed centrifugation. All tubes feature graduations every  $100\mu L$  allowing for convenient sample estimation.

- Rated to 20,000 x g for high-speed centrifugation.
- Crystal-clear material allows visualization of hard-to-see pellets.
- Available with proprietary Low-Retention resin to ensure maximum sample recovery.
- Autoclavable.
- Certified free from DNase, RNase, DNA, and PCR inhibitors.
- Made in the USA.

Part Number	Description	Quantity	List Price	20% Promo
87-B600-C	0.6mL Boil-Proof Homopolymer, Clear	1000/Pack	26.95	21.56
87-B600-AS	0.6mL Boil-Proof Homopolymer , Assorted	1000/Pack	33.95	27.16
87-B600-C-LR	0.6mL Boil-Proof Homopolymer, Low retention, Clear	500/Pack	24.95	19.96
87-B600-CS	0.6mL Boil-Proof Homopolymer, Clear, Sterile	500/Pack	26.95	21.56
87-B150-C	1.5mL Boil-Proof Homopolymer, Clear	500/Pack	14.95	11.96
87-B150-AS	1.5mL Boil-Proof Homopolymer, Assorted	500/Pack	17.95	14.36
87-B150-C-LR	1.5mL Boil-Proof Homopolymer, Low retention, Clear	250/Pack	14.95	11.96
87-B150-CS	1.5mL Boil-Proof Homopolymer, Clear, Sterile	250/Pack	15.95	12.76
87-B175-C	1.7mL UltraClear Copolymer, Clear	500/Pack	21.95	17.56
87-B200-C	2.0mL Boil-Proof Homopolymer, Clear	500/Pack	17.95	14.36
87-B200-AS	2.0mL Boil-Proof Homopolymer, Assorted	500/Pack	19.95	15.96
87-B200-C-LR	2.0mL Boil-Proof Homopolymer, Low retention, Clear	500/Pack	18.95	15.16
87-B200-CS	2.0mL Boil-Proof Homopolymer, Clear, Sterile	250/Pack	15.95	12.76
48-B500-C	5.0mL Macro tube, Clear	200/Pack	53.95	43.16
48-B500-CS	5.0mL Macro tube, Clear, Sterile	200/Pack	59.95	47.96
48-B500-AS	5.0mL Macro tube, Assorted	200/Pack	56.95	45.56

Valid to March 31, 2016



4517 Dobrin, St. Laurent, QC (Canada) H4R 2L8 Tel: 514.335.3433 | 1.800.463.4163 Fax: 514.335.0992

ISO 9001:2008 & 13485:2003