


PROGENE® PCR® TUBES AND STRIPS

- Certified RNase, DNase, DNA and Pyrogen-Free
- All plates feature high contrast printed alphanumeric
- Ultra-thin and consistent wall thickness results in rapid thermal transfer
- Made in the USA



PROGENE®

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT.#	DESCRIPTION	QUANTITY	PRICE
Individual Tubes			
87-C200-F	0.2 ml - Flat-Cap	1000/Pack	35.95
Strip Tubes			
87-C200-8-TB	0.2 ml 8-Strip Tubes	125/Pack	59.95
87-C200-8-IND	0.2 ml 8-strip tubes with Optical Flat Caps	120/Pack	71.95
87-C100-4-RG	0.1 ml 4-Strip with Caps for RotorGene®	250/Pack	111.95
Strip Caps			
87-C8-QPCR	8-Strip Optical Flat Caps	125/Pack	21.95
87-C8-DOME	8-Strip Domed Caps	125/Pack	17.95

PROGENE® SEALING FILMS

- Made in USA
- Certified RNase, DNase, DNA and Pyrogen-Free
- All place feature high contrast printed alphanumeric
- Ultra-thin and consistent wall thickness results in rapid thermal transfer



87-CFILM-TS



87-CFILM-AL



87-CFILM-QPCR

PROGENE®

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT.#	DESCRIPTION	QUANTITY	PRICE
Sealing Films			
87-CFILM-QPCR	Pressure sensitive qPCR sealing film	100 films	162.95
87-CFILM-AL	Pierceable Aluminum foil	100 films	74.95
87-CFILM-TS	Adhesive routine PCR sealing film	100 films	61.95

BUY ANY 4 GET 1 FREE!


Valid To November 30th 2016

Quote # Q387-2

PROMOTION

PROGENE® 96 WELL PCR® PLATES

- Well ridges reduce sample evaporation when sealing with film, foil or associated strip caps.
- High contrast black alphanumeric for easy well identification
- Ultra thin, uniform wells ensure optimal heat transfer for high reaction efficiency
- Certified free from DNase, RNase, DNA and PCR inhibitors
- Made in USA




87-C96-NS



87-C96-SK

PROGENE®

CERTIFIED
RNase, DNase & Pyrogen-Free



87-C96-ABI-2



87-C96-FAST



87-C96-LC480-C

CAT.#	DESCRIPTION	VOLUME	QUANTITY	PRICE
PCR Plates				
87-C96-NS	Non-Skirted	200 µl	10/Pack, 50/Case	115.95
87-C96-FAST	Half-Skirted - ABI Fast PCR®	100 µl	10/Pack, 50/Case	118.95
87-C96-ABI-2	Half-Skirted - ABI Standard	200 µl	10/Pack, 50/Case	115.95
87-C96-SK	Full-Skirted - Low-Profile	200 µl	10/Pack, 50/Case	114.95
87-C96-LC480-C	Plate, LightCycler Low profile, Clear	200 µl	10/pack	52.95
87-C96-LC480-W	Plate, LightCycler Low profile, White	200 µl	10/pack	52.95