



SPRING SUMMER 2023

FEATURED PRODUCTS




Quote # 2023-Q3

VALID TO AUGUST 31ST 2023

CONTENTS

	A	Liquid Handling	3
	B	Tubes and Vials..	8
	C	PCR and Sealing Film	14
	D	Cell Culture and Microplates.....	20
	E	Histology, Pathology and Cuvettes	25
	F	Gloves and Apparel	29
	G	Laboratory Essentials and Racks.....	34
	H	Equipment	44

EXCLUSIVE PROMOTIONS



CELLTREAT™ LOW-FORCE TIP
SOLUTION PIPETTE TIPS (LFTS)



**BUY 3
GET 1 FREE**

details on page 5



AURELIA® SONIC
POWDER-FREE NITRILE GLOVES
aurelia

**BUY 6
GET 4 FREE**

details on page 29



PROGENE®
VACUUM FILTRATION SYSTEM
PROGENE

**PROMO
PRICE**

details on page 22

ULTIDENT[®] BRAND UNIVERSAL-FIT LOW-RETENTION FILTER TIPS, **STERILE**

- Manufactured from super clear virgin polypropylene
- Certified Non-Pyrogenic & DNase/RNase-free - Gamma sterile
- Free of Nucleic Acids, Pyrogens/Endotoxins, PCR Inhibitors and trace Metals
- Tips are super hydrophobic, significantly reducing sample-loss and increasing precision
- Provides high reproducibility
- Extended length (Reach) tips available in 10 µl and 1250 µl sizes

Filter barrier made of high quality polyolefin ensures sample integrity

Super hydrophobic technology for maximum sample recovery




**BUY 4,
GET 1
FREE**

equal or lower price



Compatible with all popular pipettes



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
50-6010	10 µl, Graduated	10 Hinged Racks of 96	54.95
50-6010-XL	10 µl, Reach, Graduated	10 Hinged Racks of 96	54.95
50-6011	20 µl, Graduated	10 Hinged Racks of 96	54.95
50-6012	100 µl, Graduated	10 Hinged Racks of 96	54.95
50-6013	200 µl, Graduated	10 Hinged Racks of 96	54.95
50-6014	300 µl, Graduated	10 Hinged Racks of 96	54.95
50-6015	1000 µl, Graduated	10 Hinged Racks of 96	69.95
50-6015-XL	1250 µl, Reach, Graduated	10 Hinged Racks of 96	74.95

PROGENE[®] UNIVERSAL-FIT TIPS

- Tips are fully autoclavable - Choose Research-Grade or Low-Retention formats
- All products made with colorants are guaranteed to be heavy metals-free
- Manufactured from virgin, medical-grade FDA-approved, bovine-free plastic resins in medically-clean facility
- All molds are diamond-polished to a glass-like finish ensuring optimum performance
- Bulk Formats: Packaged in resealable Zip-Lock[®] bags
- **Low-Retention Formats:** Minimal sample retention and maximum sample recovery



CERTIFIED
RNase, DNase & Pyrogen-Free


CAT. #	DESCRIPTION	PACKAGING	PRICE
Research-Grade			
87-A10	10 µl, Graduated	1000/Pk	13.95
87-A10XL	10 µl, Reach, Graduated	1000/Pk	13.95
87-A200C	200 µl, Graduated, Clear	1000/Pk	12.95
87-A200Y	200 µl, Graduated, Yellow	1000/Pk	12.95
87-A1000B	1000 µl, Graduated, Blue	1000/Pk	18.95
Low-Retention			
87-A10-LR	10 µl, Graduated	1000/Pk	17.95
87-A10XL-LR	10 µl, Reach, Graduated	1000/Pk	17.95
87-A200-LR	200 µl, Graduated	1000/Pk	16.95
87-A300-LR	300 µl, Reach, Graduated	1000/Pk	24.95
87-A1000C-LR	1000 µl, Graduated	1000/Pk	19.95
87-A1250	1250 µl, Reach, Graduated	1000/Pk	26.95



LIQUID HANDLING

VERTEX® UNIVERSAL-FIT LOW-RETENTION FILTER TIPS, **STERILE**

- Ultra-low retention for unmatched accuracy
- Ultra-clear resin for easy sample visualization
- Space saving, multichannel friendly-rack design
- Bold, Autoclavable front-of-rack labels with volumes and key information
- Graduated for easy volume estimation
- Unsurpassed Value



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

Vertex®

CAT. #	DESCRIPTION	PACKAGING	PRICE
4117NSFS	10 µL, Graduated	10 Hinged Racks of 96	64.95
4137NSFS	10 µL, Reach, Graduated	10 Hinged Racks of 96	64.95
4237NAFS	20 µL, Graduated	10 Hinged Racks of 96	64.95
4237NBFS	100 µL, Graduated	10 Hinged Racks of 96	64.95
4237NSFS	200 µL, Graduated	10 Hinged Racks of 96	64.95
4927NSFS	300 µL, Graduated	10 Hinged Racks of 96	64.95
4337NSFS	1000 µL, Graduated	10 Hinged Racks of 96	74.95
4347NSFS	1250 µL, Reach, Graduated	10 Hinged Racks of 96	74.95

Competitor's Tip



Vertex Tip



Less is more!

VERTEX® UNIVERSAL-FIT ECOPAC® NOSTICK® REFILL SYSTEM

- Designed to specifically fit Vertex Tip Boxes
- Color coded (red = 10µl, yellow = 200µl, blue = 1250µl)
- Small footprint to minimize benchtop and shelf space
- Multichannel friendly
- Fully autoclavable and recyclable



Patented, press-to-open rack technology enables easy onehanded operation while also remaining securely closed when not in use.

Designed to facilitate offset pipetting and accommodate multichannel pipettors.

Imprinted with tip size and lot number data, visible even when seated inside the rack.

Vertex®

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
4135N00S	10 µL, Reach, Graduated	10 Wafers of 96/Pk	37.95
4235N00S	200 µL, Graduated	10 Wafers of 96/Pk	33.95
4345N00S	1250 µL, Reach, Graduated	8 Wafers of 96/Pk	42.95
87-ARFL-10-300	Empty Box for 10 to 300 µl tips	Each	4.95
87-ARFL-1250	Empty Box for 1250 µl tips	Each	4.95


LIQUID HANDLING

CELLTREAT™ LOW-FORCE TIP SOLUTION PIPETTE TIPS (LFTS)

- LFTS style pipette tips are uniquely designed to pair with pipettes featuring LFTS technology
- Certified RNase, DNase & Pyrogen-free
- Tips are designed to provide a more comfortable pipetting experience
- Lids feature stacking supports, and are removable to save bench space
- 10 µl tips feature an extended length for easier access to the bottom of vessels

**BUY 3,
GET 1
FREE**


equal or lower price



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CELLTREAT®
SCIENTIFIC PRODUCTS

CAT. #	DESCRIPTION	PACKAGING	PRICE
Racked & Sterile			
229074	10 µl Reach	10 Hinged Racks of 96	71.95
229075	200 µl	10 Hinged Racks of 96	77.95
229076	1000 µl	10 Hinged Racks of 96	94.95
Racked & Sterile with Filter			
229071	10 µl Reach	10 Hinged Racks of 96	70.95
229072	200 µl	10 Hinged Racks of 96	80.95
229073	1000 µl	10 Hinged Racks of 96	109.95



BIOMEK® FX-TYPE ROBOTIC TIPS FOR AP96 HEADS


- Certified RNase, DNase, DNA, PCR inhibitors and Pyrogen-free
- No setup changes necessary - Use existing tip definitions
- Engineered for Beckman Coulter® Biomek® and FX Span-8 Automated Systems

SSIbio

CERTIFIED
RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
4612-00	20 µl Biomek FX-Type Tips, Racked	10 Racks of 96	55.95
4612-S0	20 µl Biomek FX-Type Tips, Racked & Sterile	10 Racks of 96	61.95
4612-SF	20 µl Biomek FX-Type Tips, Racked & Sterile with Filter	10 Racks of 96	137.95
4632-00	50 µl Biomek FX-Type Tips, Racked	10 Racks of 96	55.95
4632-S0	50 µl Biomek FX-Type Tips, Racked & Sterile	10 Racks of 96	61.95
4632-SF	50 µl Biomek FX-Type Tips, Racked & Sterile with Filter	10 Racks of 96	137.95
4622-00	250 µl Biomek FX-Type Tips, Racked	10 Racks of 96	55.95
4622-S0	250 µl Biomek FX-Type Tips, Racked & Sterile	10 Racks of 96	61.95
4622-SF	250 µl Biomek FX-Type Tips, Racked & Sterile with Filter	10 Racks of 96	137.95




LIQUID HANDLING

LABCON™ ZAP™ SUPERSLIK™ UNIVERSAL FIT FILTER BARRIER PIPETTE TIPS, **STERILE**

- Earth Friendly® - Made with Solar Power™
- Autoclavable at 122°C for a maximum of 20 minutes at 15 PSI
- Certified RNase, DNase, DNA, ATP, Heavy Metal and Pyrogen-Free
- These tips include an endotoxin-free 10 micron patented porous hydrophobic filter
- Low Carbon footprint - Packaging made from 100% renewable materials
- Tips made from pure resins
- Superslik™: Low retention




CAT. #	DESCRIPTION	PACKAGING	PRICE
1171-965	0.1-10 µl UltraFine™ Point	12 Racks of 96	119.95
1172-965	0.1-10 µl Reach with TubeGuard	10 Racks of 96	119.95
1173-965	0.1-10 µl Reach with UltraFine™ Point	10 Racks of 96	119.95
1175-965	40 µl UltraFine™ Point	10 Racks of 96	119.95
1174-965	100 µl UltraFine™ Point	10 Racks of 96	119.95
1179-965	200 µl UltraFine™ Point	10 Racks of 96	119.95
1178-965	300 µl Bevel Point™	10 Racks of 96	119.95
1177-965	1000 µl	6 Racks of 96	119.95



LABCON™ ECLIPSE™ UNIVERSAL-FIT PIPETTE TIP REFILLS

- Earth Friendly® - Made with Solar Power™
- Eclipse® pipette tip refills: Simply lift, place, and snap
- Autoclavable at 122°C for a maximum of 20 minutes at 15 psi
- Certified RNase, DNase, DNA, ATP, Heavy Metal and Pyrogen-Free
- **10 µl with UltraFine™ Point:** Graduated at 2 and 10 µl - 31 mm in length
- **10 µl with TubeGuard:** TubeGuard prevents pipette barrel from touching inside wall of tubes - 38 mm in length
- Low Carbon footprint - Packaging made from 100% renewable materials
- Tips made from pure resins



CAT. #	DESCRIPTION	PACKAGING	PRICE
Eclipse™ Refills - Research-Grade			
1038-260	0.1-10 µl UltraFine™ Point	1248/Refill (13x96)	68.95
1036-260	0.1-10 µl Reach with TubeGuard	1344/Refill (14x96)	57.95
1030-260	200 µl Bevel Point™, Yellow	960/Refill (10x96)	30.95
1093-260	200 µl Bevel Point™, Graduated	960/Refill (10x96)	32.95
1014-260	300 µl UltraFine™ Point	576/Refill (6x96)	21.95
1045-260	1250 µl UltraFine™ Point, Graduated	480/Refill (5x96)	35.95


Empty Racks for Eclipse™ Refills

1280-960	Empty Rack for 10, 200 & 300 µl Refills	10/Pk	21.95
1279-960	Empty Rack for 1250 µl Refills	6/Pk	21.95

VERTEX NoStick® MULTI-VOLUME XL TIPS

One tip does it all!

- Graduations at 10, 50, 100, 150 and 200 µl
- Eliminate the need to keep multiple tips in your inventory.
- Universal fit for broad compatibility with most major pipettor brands
- Extended length 91 mm
- Filter placement is optimized to accommodate 20 µl, 50 µl, 100 µl, and 200 µl pipettors



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
4287N00	300 µl XL Non-Filtered, Non-sterile	2 Hinged Racks of 96	21.95
4287NSFS	200 µl XL Filtered, Sterile	2 Hinged Racks of 96	34.95

ULTIDENTBRAND TRANSFER PIPETTES

- Molded from LDPE - Choose from non-sterile and sterile options
- Inert to biological fluids and most acids - Can be sealed and frozen down to -196 °C




UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Non-Sterile			
170-CTB700-1	5.0 ml Graduated	500/Pk	14.95
170-CTB700-212	3.0 ml Graduated, Small Bulb	500/Pk	14.95
170-CTB700-251	4.0 ml Narrow-Stem	500/Pk	14.95
170-CTB700-274	5.0 ml Graduated, Blood Bank	500/Pk	14.95
170-CTB700-3	7.0 ml Graduated	500/Pk	14.95
Sterile			
170-CTB700-1S	5.0 ml Graduated	20/Pk, 400/Cs	36.95
170-CTB700-1SPP	5.0 ml Graduated	1/Pk, 400/Cs	43.95
170-CTB700-3S	20/Pk, 400/Cs	36.95	
170-CTB700-3SPP	7.0 ml Graduated	1/Pk, 400/Cs	43.95
170-CTB700-233S	5.0 ml Fine-Tip	20/Pk, 400/Cs	53.95

PROGENE® LOW-RETENTION GEL LOADING PIPETTE TIPS

- Choose from round and flat tip options
- Racks feature a patented lid-lock design for easy, one-hand operation
- Low-binding resin reduces sample loss and enhances sample delivery
- Extended length (91 mm / 3.58 inch) easily reaches the bottom of a 12 x 75 mm tube



CERTIFIED

RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-A10-GEL-F37	10 µl 0.37 mm Flat, Racked	2 Hinged Racks of 96	84.95
87-A10-GEL-R	10 µl 0.57 mm Round, Racked	2 Hinged Racks of 96	34.95
87-A10-GEL-RS	10 µl 0.57 mm Round, Racked, Sterile	2 Hinged Racks of 96	39.95
87-A200-GEL	200 µl 0.57 mm Round, Bulk	1000/Bag	69.95
87-A200-GEL-R	200 µl 0.57 mm Round, Racked	2 Hinged Racks of 96	26.95
87-A200-GEL-RS	200 µl 0.57 mm Round, Racked, Sterile	2 Hinged Racks of 96	37.95
87-A200-GEL-F37	200 µl 0.37 mm Flat, Racked	2 Hinged Racks of 96	84.95
87-A200-GEL-F37S	200 µl 0.37 mm Flat, Racked, Sterile	2 Hinged Racks of 96	89.95

**BUY 3,
GET 1
FREE**

equal or lower price

PROGENE® ASPIRATING PIPETS, **STERILE**

- Certified RNase, DNase and Pyrogen-free
- Supplied without graduations or cotton plugs
- Ideal for aspiration of media and sample waste removal without needing to measure the volume



CERTIFIED

RNase, DNase & Pyrogen-Free
STERILE



CAT. #	DESCRIPTION	PACKAGING	PRICE
229251	1 ml Aspirating Pipet, Bulk Packed	25/Bag, 400/Cs	71.95
229252	2 ml Aspirating Pipet, Bulk Packed	25/Bag, 400/Cs	74.95
229255	5 ml Aspirating Pipet, Bulk Packed	25/Bag, 200/Cs	47.95
229261	1 ml Aspirating Pipet, Paper/Plastic	1/Envelope, 400/Cs	72.95
229262	2 ml Aspirating Pipet, Paper/Plastic	1/Envelope, 400/Cs	76.95
229265	5 ml Aspirating Pipet, Paper/Plastic	1/Envelope, 200/Cs	48.95
229247	10 ml Aspirating Pipet, Paper/Plastic	1/Envelope, 200/Cs	53.95

TUBES AND VIALS

PROGENE® PREMIUM BOIL-PROOF MICROTUBES

- Rated to 20,000 Xg and functional from -90°C to +121°C
- 0.6, 1.5, 1.7 and 2.0 ml tubes are graduated at every 100 µl



CERTIFIED

RNase, DNase & Pyrogen-Free

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-B600-C	0.6 ml Boil-Proof	1000/Pk	25.95
87-B600-C-LR	0.6 ml Boil-Proof, Low-Retention	500/Pk	25.95
87-B600-CS	0.6 ml Boil-Proof, Sterile	500/Pk	25.95
87-B600-AS	0.6 ml Boil-Proof, Asst'd Colors	1000/Pk	33.95
87-B150-C	1.5 ml Boil-Proof	500/Pk	13.95
87-B150-C-LR	1.5 ml Boil-Proof, Low-Retention	250/Pk	13.95
87-B150-CS	1.5 ml Boil-Proof, Sterile	250/Pk	13.95
87-B150-AS	1.5 ml Boil-Proof, Asst'd Colors	500/Pk	17.95
87-B150-NC	1.5 ml Microtubes, No-Cap	500/Pk	18.95
87-B175-C	1.7 ml Boil-Proof, Large Cap	500/Pk	15.95
87-B200-C	2.0 ml Boil-Proof	500/Pk	17.95
87-B200-C-LR	2.0 ml Boil-Proof, Low-Retention	250/Pk	17.95
87-B200-CS	2.0 ml Boil-Proof, Sterile	250/Pk	17.95
87-B200-AS	2.0 ml Boil-Proof, Asst'd Colors	500/Pk	21.95
87-B200-NC	2.0 ml Microtubes, No-Cap	500/Pk	21.95
48-B500-C	5 ml Snap-Cap	200/Pk	50.95
48-B500-CS	5 ml Snap-Cap, Sterile	50/Bag, 200/Pk	53.95
48-B500-AMS	5 ml Snap-Cap, Amber, Sterile	50/Bag, 200/Pk	62.95
48-C2530	5 ml Screw-Cap	500/Cs	119.95
48-C2540	5 ml Screw-Cap, Sterile	50/Rack, 500/Cs	139.95

PROGENE® PESTLES

- Homogenize cells or resuspend precipitates such as protein or DNA pellets

PROGENE



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-PES	Microtube Pestles	10/Pk	10.95

ULTIDENTBRAND MICROTUBES, **STERILE**

- Tested up to 28,000 Xg and functional from -90°C to +121°C
- Manufactured under clean room conditions, with SuperClear co-polymer polypropylene
- N.B. Only 48-C150-CS is subdivided into bags of 50



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C600-CS	0.6 ml Microtubes, Clear	500/Pk	23.95
48-C150-CS	1.5 ml Microtubes, Clear	50/Bag, 500/Pk	29.95
48-C150-AS	1.5 ml Microtubes, Assorted	500/Pk	23.95
48-C200-CS	2.0 ml Microtubes, Clear	500/Pk	23.95

ULTIDENTBRAND CRYOGENIC COLOR DOTS

- Temperature Range: -196°C to +110°C
- Available in 15 individual colors as well as packs of 15 assorted colors
- Specially designed to comfortably fit 0.2, 0.5 / 0.6 and 1.5 / 2.0 ml microtubes
- Waterproof, they are coated with a permanent adhesive that can withstand prolonged storage



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
53-LT-7-__	7.1 mm Dots for 0.2 ml PCR Tubes	800/Pk	17.95
53-LT-7-A	7.1 mm Dots for 0.2 ml PCR Tubes, Assorted	2400/Pk	51.95
53-LT-9-__	9 mm Dots	520/Pk	15.95
53-LT-9-A	9 mm Dots, Assorted	1560/Pk	41.95
53-LT-11-__	11 mm Dots	462/Pk	15.95
53-LT-11-A	11 mm Dots, Assorted	1155/Pk	33.95
53-LT-13-__	13 mm Dots	420/Pk	15.95
53-LT-13-A	13 mm Dots, Assorted	900/Pk	30.95

Color Options: WH (White), SI (Silver), BE (Beige), YE (Yellow), LE (Lemon),
GA (Green Apple), MI (Mint), PI (Pink), LA (Lavender) BL (Blue),
OR (Orange), BR (Moka), RE (Red), VI (Violet), GS (Green Seafoam)

TUBES AND VIALS

ULTIDENTBRAND 10 ML TUBES, **STERILE**

**PROMO
PRICE**



- Leak-proof screw-caps are installed on these tubes
- Molded from premium grade polypropylene
- 3000 RCF
- 16 x 100 mm

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
11-6293	10 ml PP Centrifuge Tubes with Caps, Sterile	100/Bag, 1000/Cs	164.95

PROGENE® 87-SERIES SCREW-CAP MICROTUBES

- Non-sterile tubes and O-ring caps sold separately
- Rated to 20,000 x g, Temperature range: -196°C to +121°C

PROGENE
CERTIFIED
RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

Tubes Only, Non-Sterile

87-S500	0.5 ml, Conical	500/Pk	29.95
87-S500-SK	0.5 ml, Skirted	500/Pk	29.95
87-S150	1.5 ml, Conical, Graduated	500/Pk	29.95
87-S150-SK	1.5 ml, Skirted	500/Pk	29.95
87-S200	2.0 ml, Conical, Graduated	500/Pk	29.95
87-S200-SK	2.0 ml, Skirted, Graduated	500/Pk	29.95

O-ring caps only, Non-Sterile

87-SC-C	O-Ring Caps, Clear	500/Pk	46.95
87-SC-A	O-Ring Caps, Assorted	500/Pk	49.95

ULTIDENTBRAND 12 x 75 MM TUBES

- Molded from premium grade polypropylene
- Leak-proof screw-caps are installed on these tubes
- Raided at 3000 rpm - Autoclavable and freezable

UltiDent
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
11-6148W	12 X 75 mm PP Screw-Cap Tubes, Non-sterile	250/Bag, 1000/Cs	122.95

PROGENE® 71-SERIES SCREW-CAP MICROTUBES, **STERILE**

- Flat caps with frosted writing surface
- Packaged in tamperproof, resealable bags
- Certified Non-Pyrogenic & DNase/RNase-free - Gamma sterile
- Maximum RCF speed of 25,000 x g, Temperature range: -196°C to +121°C



Flat-top O-ring caps
protect against
evaporation

Perfect
leakproof seal at
any temperature

Excellent Clarity

Skirted tubes
are self-standing

PROGENE
CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
71-050-CS	0.5 ml, Conical, Graduated	500/Bag	79.95
71-050-SKCS	0.5 ml, Skirted	500/Bag	79.95
71-150-CS	1.5 ml, Conical, Graduated	500/Bag	79.95
71-150-SKCS	1.5 ml, Skirted	500/Bag	79.95
71-200-CS	2.0 ml, Conical, Graduated	500/Bag	79.95
71-200-SKCS	2.0 ml, Skirted, Graduated	500/Bag	79.95

SIMPORT 13 x 100 MM PS TUBES, **STERILE**

- Certified Non-pyrogenic and sterile by gamma radiation
- 8 ml capacity - Tubes are made of optically clear polystyrene and can be centrifuged up to 3000 Xg

Simport
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-T417-4S	13 X 100 mm PS Screw-Cap Tubes, Sterile	125/Bag, 1000/Cs	290.95

TUBES AND VIALS

PROGENE® CRYOGENIC VIALS, **STERILE**

- Optimized for vapor-phase liquid nitrogen storage
- Temperature range -196°C to 121°C
- Vials feature a white writing surface
- Made from virgin polypropylene and caps
- Includes secure O-ring seals
- Packaged in heavy-duty, tamper-proof, resealable bags



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
Internal-Thread			
71-ITCV-2SK	2.0 ml Internal-Thread, Skirted	50/Bag, 1000/Cs	199.95
71-ITCV-4SK	4.0 ml Internal-Thread, Skirted	50/Bag, 1000/Cs	199.95
External-Thread			
71-ETCV-2SK	2.0 ml External-Thread, Skirted	50/Bag, 1000/Cs	199.95
71-ETCV-5SK	5.0 ml External-Thread, Skirted	50/Bag, 1000/Cs	229.95

PROGENE® IPA-FREE CONTROLLED-RATE FREEZING CONTAINERS

- Ideal for cryopreservation of most cells and cell lines
- Compatible with standard 1.0 or 2.0 ml cryogenic vials
- Alcohol-free - Ensures controlled-rate -1°C/minute cell freezing in a -80°C freezer
- Polyethylene foam container and thermo-conductive core improves the rate of heat removal and provide reproducible cell cryopreservation



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
170-7577-P	12-Vial Cooling Container, Purple	Each	109.95

PROGENE® 5 & 10 ML TRANSPORT TUBES, **STERILE**

- Sterile, RNase, DNase, Pyrogen free
- Meets 95kPa transportation guidelines
- Choose from 5 ml (16 x 60 mm) and 10 ml (16 x 80 mm) formats
- Polypropylene construction with white printed graduations and marking area

PROGENE

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C1811	5 ml Transport Tubes	100/Bag, 1000/Cs	207.95
48-C1812	10 ml Transport Tubes	100/Bag, 1000/Cs	247.95

**BUY 2,
GET 1
FREE**

equal or lower
price



SSIbio® CRYOFREEZE® VIALS

- Made in the USA from virgin, medical grade
- FDA approved polypropylene with HDPE caps
- High-contrast graduations with large writing area



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

SSIbio

CAT. #	DESCRIPTION	PACKAGING	PRICE
External Thread Vials with Lip Seal Cap			
SSI-6212-S0-45000	1.0 ml Vial - Stargrip Bottom	100/Pk	39.95
SSI-6321-S0-45000	1.8 ml Vial - Round Bottom	100/Pk	39.95
SSI-6222-S0-45000	1.8 ml vial - Stargrip Bottom	100/Pk	39.95
SSI-6242-S0-45000	4.5 ml vial - Stargrip Bottom	100/Pk	48.95
Internal Thread Vials with Silicone O-ring Seal Cap			
SSI-6112-S0-45000	1.0 ml Vial - Stargrip Bottom	100/Pk	39.95
SSI-6121-S0-45000	1.8 ml Vial - Round Bottom	100/Pk	39.95
SSI-6122-S0-45000	1.8 ml Vial - Stargrip Bottom	100/Pk	39.95
SSI-6131-S0-45000	3.6 ml Vial - Round Bottom	100/Pk	52.95
SSI-6132-S0-45000	3.6 ml Vial - Stargrip Bottom	100/Pk	52.95
SSI-6141-S0-45000	4.5 ml Vial - Round Bottom	100/Pk	51.95
SSI-6142-S0-45000	4.5 ml Vial - Stargrip Bottom	100/Pk	51.95

CryoFreeze® Grip Rack

SSI-5160-0GS	40-Place Workstation	Each	32.95
--------------	----------------------	------	--------------

TUBES AND VIALS

PROGENE® CENTRIFUGE TUBES, STERILE

- Choose from bulk, recyclable paper racks or reusable polypropylene racks
- Manufactured from medical-grade virgin polypropylene with flat polyethylene (HDPE) caps - Can be used down to -80°C
- Easy-to-read black graduations are accurate to $\pm 2\%$ - 1 ml increments (15's) or 2.5 ml increments (50's)
- Both the graduations and writing areas are chloroform-resistant



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	RCF	PACKAGING	PRICE
71-1500-B	15 ml, Bulk	12,000 Xg	25/Bag, 500/Cs	89.95
71-1500-R	15 ml, Paper Racks	12,000 Xg	25/Rack, 500/Cs	109.95
71-1500-PP	15 ml, Polypropylene Racks	12,000 Xg	50/Rack, 300/Cs	126.95
71-5000-B	50 ml, Bulk	12,000 Xg	25/Bag, 500/Cs	129.95
71-5000-R	50 ml, Paper Racks	12,000 Xg	25/Rack, 500/Cs	139.95
71-5000-PP	50 ml, Polypropylene Racks	12,000 Xg	25/Rack, 300/Cs	159.95

LABCON™ PerformR™ CENTRIFUGE TUBES, STERILE

- Earth Friendly® - Made with Solar Power™
- Storage Temperature Range -80°C to +40°C (-112°F to 104°F)
- Certified RNase, DNase, DNA, ATP, Heavy Metal and Pyrogen-Free
- Low Carbon footprint - Packaging made from 100% renewable materials
- Easy-to-Read black graduations with white solvent-resistant labeling area
- Racked tubes are packaged in a reusable, autoclavable and freezable polypropylene racks



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

labcon

CAT. #	DESCRIPTION	RCF	PACKAGING	PRICE
3136-345	15 ml, Bulk	12,500 Xg	50/Bag, 500/Cs	139.95
3146-335	15 ml, Racked	12,500 Xg	25/Rack, 300/Cs	134.95
3186-345	50 ml, Bulk	12,500 Xg	50/Bag, 500/Cs	199.95
3196-335	50 ml, Racked	12,500 Xg	25/Rack, 300/Cs	199.95

MACROTUBE® VIVID™, STERILE

- Centrifuge up to 25,000 x g
- Use from -86°C to 110°C
- Unsurpassed Clarity
- Printed graduations every 0.2 and 0.5 ml
- Comes in snap apart racks (2 x 25 tubes)

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C2545	5 ml Screw-Cap Tubes, Sterile	10 Foam Racks of 50/Cs	169.95

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

**BUY 2,
GET 1
FREE**



TUBES AND VIALS

ULTIDENTBRAND CULTURE TUBES WITH 2-POSITION CAPS, **STERILE**

- Packaged in convenient resealable zip lock bags
- Polystyrene: Rated at 1400 Xg and temperatures to 70°C
- Polypropylene: Rated at 3000 Xg and temperatures to 120°C



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
170-T4100-S	12x75, Polypropylene	25/Bag, 500/Cs	119.95
170-T4200-S	12x75, Polystyrene	25/Bag, 500/Cs	119.95
170-T4300-S	17x95, Polypropylene	25/Bag, 500/Cs	149.95
170-T4400-S	17x95, Polystyrene	25/Bag, 500/Cs	149.95

ULTIDENTBRAND FLANGE PLUG CAPS

- Flexible flanges to ensure a leak-proof seal



UltiDent
Scientific


CAT. #	DESCRIPTION	PACKAGING	PRICE
170-CTC12MM-_*	12 mm Caps	1000/Pk	23.95
170-CTC13MM-_*	13 mm Caps	1000/Pk	25.95
170-CTC16MM-_*	16 mm Caps	1000/Pk	35.95

Color Options: B (Blue), BK (Black), N (Natural), G (Green), M (Mauve), R (Red), Y (Yellow), W (White)

ULTIDENTBRAND DISPOSABLE BOROSILICATE CULTURE TUBES

- Chemically inert - High heat tolerance - Convenient dispenser packaging

UltiDent
Scientific



CAT. #	DESCRIPTION	VOLUME	PACKAGING	PRICE
170-14100003B	10 x 75 mm	3 ml	250/Pk, 1000/Cs	39.95
170-14100005B	12 x 75 mm	5 ml	250/Pk, 1000/Cs	44.95
170-14100007B	13 x 100 mm	7 ml	250/Pk, 1000/Cs	49.95
170-14100010B	16 x 100 mm	14 ml	250/Pk, 1000/Cs	69.95
170-14106125B	16 x 125 mm	16 ml	100/Pk, 600/Cs	59.95
170-14106150B	16 x 150 mm	24 ml	250/Pk, 1000/Cs	94.95
170-14108150B	18 x 150 mm	28 ml	100/Pk, 600/Cs	66.95
170-14120150B	20 x 150 mm	36 ml	100/Pk, 600/Cs	99.95

ULTIDENTBRAND 16 x 100 MM TRANSPORT / CULTURE TUBES, **STERILE**

**NEW
PRODUCT**

- Leak-proof screw caps, packaged sterile in bags of 100
- Choose from polypropylene with a skirted bottom (14 ml) or in crystal clear polystyrene with a round bottom (12 ml)
- Tubes have easy to read graduations and a white marking area
- Tubes measure 16 x 100 mm



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-T7950	14 ml PP w/Skirted Bottom	100/Bag, 1000/Cs	334.95
48-T7955	12 ml PS w/Round Bottom	100/Bag, 1000/Cs	334.95

ULTIDENTBRAND BOROSILICATE MEDIA BOTTLES

- Bottles have permanent white enamel graduations and drip-free pour rings
- Made from 3.3 borosilicate glass, bottles and screw caps can be autoclaved
- Use for storage and mixing of media, buffers, aqueous and nonaqueous solutions



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
170-14170100	100 ml Media Bottles	10/Cs	43.95
170-14170250	250 ml Media Bottles	10/Cs	59.95
170-14170500	500 ml Media Bottles	10/Cs	85.95
170-14171000	1000 ml Media Bottles	10/Cs	119.95

TUBES AND VIALS

ULTIDENTBRAND BOROSILICATE VOLUMETRIC LABWARE



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
170-11010025	25 ml Low-Form Beaker	Each	5.95
170-11010050	50 ml Low-Form Beaker	Each	5.95
170-11010100	100 ml Low-Form Beaker	Each	5.95
170-11010150	150 ml Low-Form Beaker	Each	5.95
170-11010250	250 ml Low-Form Beaker	Each	5.95
170-11010300	300 ml Low-Form Beaker	Each	5.95
170-11010400	400 ml Low-Form Beaker	Each	5.95
170-11010500	500 ml Low-Form Beaker	Each	6.95
170-11010600	600 ml Low-Form Beaker	Each	6.95
170-11011000	1000 ml Low-Form Beaker	Each	17.95
170-11012000	2000 ml Low-Form Beaker	Each	20.95
170-11013000	3000 ml Low-Form Beaker	Each	22.95
170-11014000	4000 ml Low-Form Beaker	Each	27.95
170-11210050	50 ml Erlenmeyer Flask	Each	6.95
170-11210100	100 ml Erlenmeyer Flask	Each	6.95
170-11210150	150 ml Erlenmeyer Flask	Each	6.95
170-11210200	200 ml Erlenmeyer Flask	Each	6.95
170-11210250	250 ml Erlenmeyer Flask	Each	8.95
170-11210300	300 ml Erlenmeyer Flask	Each	8.95
170-11210500	500 ml Erlenmeyer Flask	Each	8.95
170-11211000	1000 ml Erlenmeyer Flask	Each	11.95
170-11212000	2000 ml Erlenmeyer Flask	Each	16.95
170-11213000	3000 ml Erlenmeyer Flask	Each	21.95
170-11215000	5000 ml Erlenmeyer Flask	Each	29.95
170-16010005	5 ml Hex-Base Cylinder	Each	8.95
170-16010010	10 ml Hex-Base Cylinder	Each	9.95
170-16010025	25 ml Hex-Base Cylinder	Each	10.95
170-16010050	50 ml Hex-Base Cylinder	Each	11.95
170-16010100	100 ml Hex-Base Cylinder	Each	14.95
170-16010250	250 ml Hex-Base Cylinder	Each	16.95
170-16010500	500 ml Hex-Base Cylinder	Each	21.95
170-16011000	1000 ml Hex-Base Cylinder	Each	28.95
170-16012000	2000 ml Hex-Base Cylinder	Each	34.95

ULTIDENTBRAND POLYPROPYLENE LABWARE

• Translucent polypropylene is autoclavable at 121°C for 20 minutes



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
11-601821	25 ml Griffin Beakers	Each	7.95
11-601822	50 ml Griffin Beakers	Each	8.95
11-601823	100 ml Griffin Beakers	Each	8.95
11-601824	250 ml Griffin Beakers	Each	10.95
11-601825	500 ml Griffin Beakers	Each	12.95
11-601826	1000 ml Griffin Beakers	Each	20.95
11-601827	2000 ml Griffin Beakers	Each	29.95
11-600818	3000 ml Griffin Beakers	Each	49.95
11-600819	5000 ml Griffin Beakers	Each	65.95
11-602560	10 ml Cylinders	Each	11.95
11-602561	25 ml Cylinders	Each	13.95
11-602562	50 ml Cylinders	Each	15.95
11-602563	100 ml Cylinders	Each	21.95
11-602564	250 ml Cylinders	Each	21.95
11-602565	500 ml Cylinders	Each	26.95
11-602566	1000 ml Cylinders	Each	39.95
11-602567	2000 ml Cylinders	Each	66.95
200-TCB-50	50 ml Tri-Cornered Beakers	25/Pack	16.95
200-TCB-100	100 ml Tri-Cornered Beakers	50/Pack	33.95
200-TCB-250	250 ml Tri-Cornered Beakers	25/Pack	22.95
200-TCB-400	400 ml Tri-Cornered Beakers	25/Pack	23.95
200-TCB-800	800 ml Tri-Cornered Beakers	25/Pack	31.95
200-TCB-1000	1000 ml Tri-Cornered Beakers	50/Pack	62.95
170-F5569-PP	Ø 55 x 69 mm Short Stem Plastic Funnels	25/Pack	9.95
170-F6578-PP	Ø 65 x 78 mm Short Stem Plastic Funnels	25/Pack	9.95

ULTIDENTBRAND DRYING RACK

- Holds up to 77 lbs.
- Made of durable, easy to clean polystyrene
- Includes: 72 Pegs, 2 Drying Baskets, Drainage Tubing and Mounting Hardware



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120791	Drying Rack 17.4 x 26 x 5.5 in	Each	328.95

PCR AND SEALING FILM

PROGENE® INDIVIDUAL THIN-WALL PCR® TUBES

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Evaporation-resistant, easy-to-open-and-close lids prevent sample loss
- Packaged in resealable bags to minimize risk of contamination



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C200-F	0.2 ml - Flat-Cap, Frosted	1000/Pk	42.95
87-C200-CF	0.2 ml - Flat-Cap, Optical	1000/Pk	42.95
87-C200-D	0.2 ml - Dome-Cap	1000/Pk	49.95
87-C200-NC	0.2 ml - No Cap	1000/Pk	64.95
87-C500-F	0.5 ml - Flat-Cap, Frosted	1000/Pk	42.95

PROGENE® 8 & 12-STRIP THIN-WALL PCR® TUBES & CAPS

- Ultra-thin walls produce efficient thermal transfers during amplifications



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C100-8-TB	0.1 ml 8-Strip Tubes, Low-Profile	125/Pk	85.95
87-C200-8-TB	0.2 ml 8-Strip Tubes	125/Pk	65.95
87-C200-12-TB	0.2 ml 12-Strip Tubes	80/Pk	89.95
87-C8-DOME	8-Strip Domed Caps	125/Pk	19.95
87-C8-QPCR	8-Strip Optical qPCR Caps	125/Pk	20.95
87-C12-DOME	12-Strip Domed Caps	80/Pk	22.95
87-C12-QPCR	12-Strip Optical qPCR Caps	80/Pk	22.95

PROGENE® 8-STRIP TUBES WITH INDIVIDUAL-ATTACHED CAPS

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Angular geometry caps prevent lids and hinges from interfering with each other
- **Rigid:** Robust tube linkages make the strips rigid and easy to handle
- **Tear-Strip:** Easily customize strips into smaller sections



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C100-8-IND	0.1 ml Clear - Optical qPCR Caps, Tear Strip	120/Pk	84.95
87-C100-8-IND-W	0.1 ml White - Optical qPCR Caps, Rigid	120/Pk	88.95
87-C200-8-IND	0.2 ml Clear - Optical qPCR Caps, Tear Strip	120/Pk	79.95
87-C200-8RGD-IND	0.2 ml Clear - Optical qPCR Caps, Rigid	120/Pk	79.95
87-C200-8RGD-IND-W	0.2 ml White - Optical qPCR Caps, Rigid	120/Pk	85.95
87-C200-8D-IND	0.2 ml Clear - Domed Caps, Tear Strip	120/Pk	79.95



Tear Strip



Rigid

PCR AND SEALING FILM

Find the Perfect Plate for Your Thermal Cycler!

CONSULT OUR **PCR PLATE SELECTION GUIDE**
ON PAGE 18 FOR MORE INFORMATION



PROGENE® THIN-WALL PCR® PLATES

- Made from high-quality molecular grade polypropylene plastic
- Broad thermal cycler compatibility
- Ultra-thin wells ensure optimal heat transfer and high reaction efficiency during amplifications
- High-contrast, black alphanumeric indexing for easy well identification



87-C96-NS



87-C96-ABI-1



87-C96-ABI-2



87-C96-FAST



87-C96-RR



87-C96-SK

CERTIFIED
RNase, DNase & Pyrogen-Free

PROGENE



CAT. #	DESCRIPTION	PACKAGING	PRICE
96-Well PCR Plates			
87-C96-NS	0.2 ml - Non-Skirted	10/Pk, 50/Cs	129.95
87-C96-LPNS	0.1 ml - Non-Skirted, Low-Profile	20/Pk	63.95
87-C96-ABI-1	0.2 ml - Half-Skirt for ABI, Hybrid	10/Pk, 50/Cs	149.95
87-C96-ABI-2	0.2 ml - Half-Skirt for ABI, Universal	10/Pk, 50/Cs	129.95
87-C96-FAST	0.1 ml - Half-Skirt for ABI Fast PCR	10/Pk, 50/Cs	149.95
87-C96-RR	0.2 ml - Half-Skirt, Raised Rim	10/Pk, 50/Cs	149.95
87-C96-SK	0.1 ml - Full-Skirt, Low-Profile	10/Pk, 50/Cs	129.95
87-C96-CT	0.2 ml - Chimney-Top	10/Pk, 50/Cs	159.95

384-Well PCR Plates

87-C384-1N	Full Skirt, 1-Notch	10/Pk, 50/Cs	230.95
87-C384-2N	Full Skirt, 2-Notch	10/Pk, 50/Cs	230.95

SIMPORT SUPERFLEX™ “BREAK-AWAY” 96 WELL PCR PLATES

Compatible with all
instruments that fit non-
skirted, low profile plates

Can be easily separated
into 8 strips or higher



Pierceable Cap Strips

CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-T320-96N	0.2ml -Non-Skirted , Low-Profile	10/Bag, 50/Case	214.95
17-T321-96N	8 x 12 Flat Optical Pierceable Cap Strips	10/Bag, 50/Case	197.95

SIMPORT AMPLATE™ NON-SKIRTED PCR® PLATES

- Can be cut to the desired format
- Printed alphanumeric grid for easier identification



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-T323-24N	0.2 ml - 24-Well Pre-Cut Plates	10/Bag, 50/Cs	46.95
17-T323-48N	0.2 ml - 48-Well Pre-Cut Plates	10/Bag, 50/Cs	68.95
17-T323-96N	0.2 ml - 96-Well Plates	10/Bag, 50/Cs	144.95

PCR AND SEALING FILM

PROGENE® THIN-WALL LIGHTCYCLER® PLATE

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Our low-profile plates are fully LightCycler® compatible



PROGENE

CERTIFIED

RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C96-LC480-WP	96W White Plates - No Film	10/Pk	32.95

PROGENE® AUTORADIOGRAPHY FILM

- High contrast, high sensitivity, low background, light blue film
- Especially useful with chemiluminescence, works well with 14C, 32P, 125I, and 35S



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
39-20810	8" x 10" Autoradiography Film	100/Pk	169.95

PROGENE® ROTOR-GENE® PCR® STRIP TUBES

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Compatible with QIAGEN®/Corbett Rotor-Gene® Instruments
- Made in the USA



PROGENE

CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C100-4-RG	0.1 ml 4-Strip with Caps for Rotor-Gene®	250/Pk	130.95

SIMPORT ROTOR-GENE® PCR® STRIPS

- RotoCycler™ strips are perfectly designed to perform on Rotor-Gene™ instruments
- Frosted caps offer a convenient area for labelling



Simport

CERTIFIED

RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
17-T319-4N	0.1 ml 4-Strip with Caps for Rotor-Gene®	250/Pk	99.95

PROGENE® PCR SEALING MAT

- Pierceable and reusable silicone creates an evaporation-resistant seal



PROGENE

CERTIFIED

RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C96-MAT	PCR Sealing Mat	5/PK	60.95

PROGENE® ROTOR-GENE® COOLER/WORKSTATION

- Maintains a sample temperature of less than 4°C for almost 4 hours (with lid on rack) when at room temperature



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W96-RG	64-Place RotorGene® Cooler with Lid	2/Pk	109.95




PLATESEAL™ PCR/qPCR OPTICAL SEALING FILM

- Polyester-based sealing film
- Ideal for PCR, qPCR (real-time) and robotics
- Choose from standard or X-Large sized film
- Temperature Range: -40°C to +120°C



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PETST-100	Clear Film, Polyester	100/Pk	99.95
PS-PETSTXL-100	Clear Film, Polyester, X-Large	100/Pk	105.95


PLATESEAL™ PCR SEALING FOIL

- Easy to pierce with pipet tips or robotic probes for direct sample recovery
- Ideal for standard PCR and long-term storage
- Temperature Range: -80°C to +120°C



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-FOILP-100	Easy-to-Pierce Foil	100/Pk	114.95
PS-FOILP-50S	Easy-to-Pierce Foil, Sterile	50/Pk	85.95


PLATESEAL™ qPCR PRESSURE-ACTIVATED SEALING FILM

- Highly chemically resistant, non-tacky, pressure-activated silicone adhesive
- For qPCR (real-time), Long-Term Storage and Protein Crystallization - Temperature Range: -80°C to +100°C



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PPO-100	qPCR, Pressure-Activated Adhesive	100/Pk	278.95
PS-PPO-50S	qPCR, Pressure-Activated Adhesive, Sterile	50/Pk	171.95


PLATESEAL™ EXTRA-THICK SEALING FOIL

- Extra thick, chemically-resistant aluminum minimizes shredding upon removal
- Ideal for standard PCR and long-term storage
- Temperature Range: -80°C to +120°C



CERTIFIED

RNase, DNase & Pyrogen-Free




CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-FOIL-100	Chemically Resistant Foil	100/Pk	123.95
PS-FOIL-50S	Chemically Resistant Foil, Sterile	50/Pk	80.95

NOVASBIO® PCR/qPCR ADHESIVE SEALING FILM

- Certified RNase, DNase and Pyrogen-free
- Temperature Range: -80°C to +120°C
- Use for standard PCR and real-time PCR applications
- Optically clear polyester film with a strong, 50 µm thick acrylic adhesive

CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
65-SF-UC-550	qPCR Optical Adhesive Seal	100/Pk	159.95

PLATESEAL™ FILM ACCESSORIES



CAT. #	DESCRIPTION	PACKAGING	PRICE
PSA-RRLL-1	Film Roller	Each	53.95
PSA-PP-5	Sealing Applicator	5/Pk	19.95

PCR AND SEALING FILM

PROGENE[®] PCR & qPCR Plates

Broad compatibility with popular thermal cyclers

- Contact us for samples!

Ultra-thin, uniform wells

- Optimal heat transfer and higher reaction efficiency.

A variety of plate sealing options

- Film, foil, silicone mats, and strip caps prevent evaporation.

**We have a plate for
virtually every Thermal Cycler.**

**Don't see your Thermal Cycler listed?
Please contact us!**


Cat No.	Skirting	Volume	Plate Description	Color	Notch	Compatibility	Packaging	Pricing
87-C96-NS	Non-skirted	0.2 mL	Standard 96-Well Plate (Cutttable)	Clear	A12	qPCR & PCR	10/Pk, 50/Cs	\$129.95
87-C96-LPNS	Non-skirted	0.1 mL	Low-Profile 96-Well Plate	Clear	H12	qPCR & PCR	20/Pk	\$63.95
87-C96-ABI-1	Semi-skirted	0.2 mL	Standard ABI [®] Style 96-Well Plate (Flat Top)	Clear	A12	qPCR & PCR	10/Pk, 50/Cs	\$149.95
87-C96-RR	Semi-skirted	0.2 mL	ABI [®] Shape 96-Well Plate (Raised Top Ridge)	Clear	A12	qPCR & PCR	10/Pk, 50/Cs	\$149.95
87-C96-FAST	Semi-skirted	0.1 mL	ABI [®] Style FAST [®] , Low-Profile 96-Well Plate	Clear	A1	qPCR	10/Pk, 50/Cs	\$149.95
87-C96-ABI-2	Semi-skirted	0.2 mL	Standard-Shape 96-Well Plate	Clear	A12	PCR	10/Pk, 50/Cs	\$129.95
87-C96-LC480-WP	Semi-skirted	0.1 mL	LightCycler [®] , Low-Profile, 96-Well Plate	White	H12	qPCR	10/Pk	\$32.95
87-C96-SK	Full-skirted	0.1 mL	Low-Profile 96-Well Plate	Clear	A12	qPCR & PCR	10/Pk, 50/Cs	\$129.95
87-C384-2N	Full-skirted	30 µL	Two-Notch Type, 384-Well Plate	Clear	A24, P24	qPCR & PCR	10/Pk, 50/Cs	\$230.95
87-C384-1N	Full-skirted	30 µL	One-Notch Type, 384-Well Plate	Clear	A24	qPCR & PCR	10/Pk, 50/Cs	\$230.95



87-C96-NS



87-C96-LPNS



87-C96-ABI-1



87-C96-RR



87-C96-FAST

PCR AND SEALING FILM



High-contrast,
black alphanumeric
indexing.



Raised wells
ensure optimal
sealing to prevent
evaporation.



- Best fit
- Also fits

Applied Biosystems®

Roche®

BIO-RAD®

Eppendorf®

7300	7500	7500 Fast	7700	7900HT	7900HT Fast	GeneAmp® 9700	StepOnePlus™	Veriti®	ViiA™ 7	QuantStudio® Fast	QuantStudio® 0.2 mL	LightCycler®	PTC-200, 220, 225	iCycler® & MyiQ	MyCycler™	C1000™ / S1000™	IQ5®	T100	Mastercycler®	Mastercycler® Gradient	Mastercycler® ep / Pro	Mastercycler® ep / Pro384	Mastercycler® Nexus
------	------	-----------	------	--------	-------------	---------------	--------------	---------	---------	-------------------	---------------------	--------------	-------------------	-----------------	-----------	-----------------	------	------	---------------	------------------------	------------------------	---------------------------	---------------------

Cat No.

87-C96-NS	●	●		●	●	●	●						●	●	●	●	●	●	●	●	●	●	●
87-C96-LPNS													●			●			●	●		●	●
87-C96-ABI-1	●	●		●	●	●	●													●			
87-C96-RR	●	●		●	●	●	●	●	●											●			
87-C96-FAST			●		●	●	●	●	●														
87-C96-ABI-2				●		●	●					●	●	●		●	●		●	●		●	●
87-C96-LC480-WP										●													
87-C96-SK												●	●			●							
87-C384-2N										●											●		
87-C384-1N				●		●						●			●								



87-C96-ABI-2



87-C96-LC480-WP



87-C96-SK



87-C384-2N



87-C384-1N

CELL CULTURE AND MICROPLATES

PROGENE® SEROLOGICAL PIPETTES, **STERILE**

- Made of polystyrene plastic
- Certified non-pyrogenic and DNase-/RNase-Free
- Sterile, individually wrapped

Features a premium material barrier plug

Bi-directional (positive and negative) graduations provide choice of ascending and descending scales


Paper/plastic wrap with color-coding on paper allows for easy size identification

Printed indicator arrow for quickly locating negative graduations

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	VOLUME	GRADUATION	COLOR	PACKAGING	PRICE
85-7125	1 ml	1/100	Yellow	500/Cs	99.95
85-7126	2 ml	1/50	Green	500/Cs	99.95
85-7127	5 ml	1/10	Blue	200/Cs	42.95
85-7128	10 ml	1/10	Orange	200/Cs	42.95
85-7129	25 ml	2/10	Red	200/Cs	96.95
85-7130	50 ml	5/10	Purple	100/Cs	132.95
85-7131	100 ml	1/100	Pink	50/Cs	149.95



SURESTRAIN™ CELL STRAINERS

- No "air-lock" effect - Faster flow through
- Each case include a 15 ml volume reducing adapter



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C4040	40 µm Cell Strainer, Blue	1/Pk, 50/Cs	141.95
48-C4070	70 µm Cell Strainer, White	1/Pk, 50/Cs	141.95
48-C4100	100 µm Cell Strainer, Yellow	1/Pk, 50/Cs	141.95

PROGENE® CELL STRAINERS

- Ideal for sample preparation and pre filtration

PROGENE

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
229481	40 µm Strainer, Blue	1/Pk, 50/Cs	79.95
229483	70 µm Strainer, White	1/Pk, 50/Cs	79.95
229485	100 µm Strainer, Yellow	1/Pk, 50/Cs	79.95
229480	Pestle, Green	1/Pk, 50/Cs	63.95

**PROMO
PRICE**



CELL CULTURE AND MICROPLATES

PROGENE® CELL CULTURE FLASKS, STERILE

- Sterilized by gamma irradiation
- Manufactured in Class 100,000 clean-room and TC-Treated by vacuum gas plasma

PROGENE

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE



**PROMO
PRICE**

CAT. #	DESCRIPTION	PACKAGING	PRICE
Tissue Culture Treated			
229321	12.5 cm ² Vent Cap	10/Bag, 200/Cs	209.95
229331	25 cm ² Vent Cap	10/Bag, 200/Cs	129.95
229341	75 cm ² Vent Cap	5/Bag, 100/Cs	129.95
229351	182 cm ² Vent Cap	5/Bag, 40/Cs	129.95
229371	225 cm ² Vent Cap	5/Bag, 25/Cs	156.95
229361	300 cm ² Vent Cap	3/Bag, 18/Cs	191.95
Non-Treated Suspension Culture			
229500	25 mL Vent Cap	10/Bag, 200/Cs	209.95
229510	50 mL Vent Cap	10/Bag, 200/Cs	215.95
229520	250 mL Vent Cap	5/Bag, 100/Cs	215.95
229530	600 mL Vent Cap	5/Bag, 40/Cs	170.95
229540	850 mL Vent Cap	3/Bag, 18/Cs	182.95

PROGENE® CELL CULTURE PLATES, **STERILE**

- Sterilized by gamma irradiation
- Manufactured in Class 100,000 clean-room and TC-Treated by vacuum gas plasma
- No sidewall growth
- Raised rims on wells with matching rings on plate lids reduce evaporation
- Supplied in easy-peel packaging for sterile product integrity



**PROMO
PRICE**

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
Tissue Culture Treated			
229105	6-Well with Lid, Flat Bottom	1/Pk, 50/Cs	89.95
229111	12-Well with Lid, Flat Bottom	1/Pk, 50/Cs	89.95
229123	24-Well with Lid, Flat Bottom	1/Pk, 50/Cs	89.95
229147	48-Well with Lid, Flat Bottom	1/Pk, 50/Cs	89.95
229195	96-Well with Lid, Flat Bottom	1/Pk, 50/Cs	89.95
229190	96-Well with Lid, U-Bottom	1/Pk, 100/Cs	239.95
Non-Treated Suspension Culture			
229506	6-Well with Lid, Flat Bottom	1/Pk, 100/Cs	190.95
229512	12-Well with Lid, Flat Bottom	1/Pk, 100/Cs	194.95
229524	24-Well with Lid, Flat Bottom	1/Pk, 100/Cs	194.95
229548	48-Well with Lid, Flat Bottom	1/Pk, 100/Cs	207.95
229596	96-Well with Lid, Flat bottom	1/Pk, 100/Cs	202.95

PROGENE® CELL CULTURE DISHES, STERILE

- Sterilized by gamma irradiation
- Manufactured in Class 100,000 clean-room and TC-Treated by vacuum gas plasma
- Our dishes offer distortion-free optical clarity and optimal gas exchange




CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
229661	60 mm x 15 mm	20/Bag, 500/Cs	176.95
667660	60 mm x 15 mm with Grip Ring	10/Bag, 500/Cs	186.95
229621	100 mm x 20 mm	10/Bag, 300/Cs	176.95
229620	100 mm x 20 mm with Grip Ring	10/Bag, 300/Cs	222.95
667651	150 mm x 20 mm	5/Bag, 100/Cs	206.95
229652	150 mm x 25 mm	5/Bag, 60/Cs	154.95

PLATESEAL™ BREATHABLE FILM

- Designed to facilitate air and CO² exchange
- Ideal for cell, tissue and bacterial cultures
- Temperature Range: -20°C to +80°C



PlateSeal™

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-RA-100	Breathable Film	100/Pk	123.95
PS-RA-50S	Breathable Film, Sterile	50/Pk	94.95

CELL CULTURE AND MICROPLATES

PROGENE® VACUUM FILTRATION SYSTEM, **STERILE**

**PROMO
PRICE**

- Low binding PES membrane for high flow rate
- 0.10 µm membrane for mycoplasma removal, 0.22 µm membrane for filtration of media, 0.45 µm PES membrane for solution clarification
- Supplied complete with separate O-ring cap for receiver bottle
- Easy-to-read molded graduations marks on bottle top filter with bold easy-to-read printed membrane and pore size for easy identification
- Angled hose barb adapter fits three size hose diameters and prevents the filter from tipping



PROGENE

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
71-PES-150-22	150 ml, PES, 0.22 µm	1/Bag, 12/Cs	89.95
71-PES-250-22	250 ml, PES, 0.22 µm	1/Bag, 12/Cs	99.95
71-PES-500-22	500 ml, PES, 0.22 µm	1/Bag, 12/Cs	129.95
71-PES-1000-22	1000 ml, PES, 0.22 µm	1/Bag, 12/Cs	199.95
71-PES-150-45	150 ml, PES, 0.45 µm	1/Bag, 12/Cs	89.95
71-PES-250-45	250 ml, PES, 0.45 µm	1/Bag, 12/Cs	99.95
71-PES-500-45	500 ml, PES, 0.45 µm	1/Bag, 12/Cs	149.95
71-PES-1000-45	1000 ml, PES, 0.45 µm	1/Bag, 12/Cs	199.95

PROGENE® BARBED VENT FILTERS

- 15 mm to 25 mm stepped Barb
- Hydrophobic membrane acts as a barrier to aqueous solutions and aerosols
- Filter housing is made of polypropylene and the membrane material is PTFE



PROGENE

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
71-235050	50 mm 0.2 µm PTFE Barbed Vent Filter	20/Pk	84.95

PROGENE® SYRINGE FILTERS

- **PES (Polyethersulfone)** - Low affinity for proteins and extractables - Faster flow rate than PVDF
- **CA (Cellulose Acetate)** - Low protein binding, reduced extractables and low hold-up volume
- **MCE (Mixed Cellulose Ester)** - Filtration of aqueous solutions - Effectively binds trace proteins
- **PVDF (Polyvinylidene fluoride)** - Low-binding non-aggressive aqueous and mild organic solutions
- **NYLON** - Universal application for analytical procedures - Chemical filtration
- **PTFE (Polytetrafluoroethylene)** - Hydrophobic - Strong acid solvent filtration



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
Sterile Blister Packs, Individual-wrapped			
229746	PES, 0.22 µm, 13 mm	75/Pk	79.95
229747	PES, 0.22 µm, 30 mm	30/Pk	73.95
229748	PES, 0.45 µm, 13 mm	75/Pk	79.95
229749	PES, 0.45 µm, 30 mm	30/Pk	73.95
229762	CA, 0.22 µm, 30 mm	30/Pk	78.95
229763	CA, 0.45 µm, 30 mm	30/Pk	78.95
229750	MCE, 0.22 µm, 13 mm	75/Pk	79.95
229751	MCE, 0.22 µm, 30 mm	30/Pk	73.95
229752	MCE, 0.45 µm, 13 mm	75/Pk	79.95
229753	MCE, 0.45 µm, 30 mm	30/Pk	73.95
229756	PTFE, 0.22 µm, 13 mm	75/Pk	146.95
229757	PTFE, 0.22 µm, 30 mm	30/Pk	104.95
229758	PTFE, 0.45 µm, 13 mm	75/Pk	146.95
229759	PTFE, 0.45 µm, 30 mm	30/Pk	104.95
229740	PVDF, 0.10 µm, 13 mm	75/Pk	88.95
229741	PVDF, 0.10 µm, 30 mm	30/Pk	73.95
229742	PVDF, 0.22 µm, 13 mm	75/Pk	88.95
229743	PVDF, 0.22 µm, 30 mm	30/Pk	73.95
229744	PVDF, 0.45 µm, 13 mm	75/Pk	88.95
229745	PVDF, 0.45 µm, 30 mm	30/Pk	73.95
Non-Sterile Bulk Packs			
229771	PES, 0.45 µm, 30 mm	100/Pk	190.95
229774	PVDF, 0.45 µm, 30 mm	100/Pk	190.95
229775	Nylon, 0.22 µm, 30 mm	100/Pk	190.95
229776	Nylon, 0.45 µm, 30 mm	100/Pk	190.95
229777	PTFE, 0.22 µm, 13 mm	100/Pk	141.95
229778	PTFE, 0.22 µm, 30 mm	100/Pk	253.95
229779	PTFE, 0.45 µm, 13 mm	100/Pk	141.95
229780	PTFE, 0.45 µm, 30 mm	100/Pk	253.95

CELL CULTURE AND MICROPLATES

PROGENE® ROLLER BOTTLES

- Sterilized by gamma irradiation
- All models feature "easy-to-read" black printed graduations
- Choose from traditional or Expanded Surface Roller Bottles (ESRB) formats
- ESRB Ribbed surface increases surface area from 850cm² to 1900cm²



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

**BUY 3,
GET 1
FREE**

equal or lower price



CAT. # DESCRIPTION

Tissue Culture Treated

229382	490 cm ² , Non-Vented Cap	1/Bag, 24/Cs	170.95
229383	490 cm ² , Vented Cap	1/Bag, 24/Cs	176.95
229384	850 cm ² , Non-Vented Cap	1/Bag, 12/Cs	114.95
229385	850 cm ² , Vented Cap	1/Bag, 12/Cs	109.95
229386	1900 cm ² ESRB, Non-Vented Cap	1/Bag, 12/Cs	125.95
229387	1900 cm ² ESRB, Vented Cap	1/Bag, 12/Cs	130.95

Non-Treated Suspension Culture

229582	1000 ml, Non-Vented Cap	1/Bag, 24/Cs	154.95
229584	2000 ml, Non-Vented Cap	1/Bag, 12/Cs	96.95

PROGENE® ROBOTIC RESERVOIRS

- Designed to assist in repetitive or robotic liquid transfer
- Made of polypropylene and chemically inert for combinatorial chemistry applications
- Pyramid and Trough Bottom designs promote funneling of liquid to maximize reagent recovery



CERTIFIED
RNase, DNase & Pyrogen-Free



**BUY 3,
GET 1
FREE**

equal or lower price

CAT. #	DESCRIPTION	PACKAGING	PRICE
229561	8 Well, Trough-Bottom	5/Bag, 25/Cs	144.95
229562	12 Well, Trough-Bottom	5/Bag, 25/Cs	144.95
229563	96 Well Pyramid-Bottom	5/Bag, 25/Cs	127.95
229564	384 Well Pyramid-Bottom	5/Bag, 25/Cs	127.95

PROGENE® BIO-REACTION TUBES, **STERILE**

- Packaged sterile (gamma irradiated) in zip-closure bags
- Convenient and economical option for small scale, high density suspension cell culture
- 0.22 µm hydrophobic membrane is ideal for biological reactions needing free gas exchange



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE




CAT. #	DESCRIPTION	RCF	PACKAGING	PRICE
229471	15 ml, Conical-Bottom	12,000 Xg	10/Bag, 100/Cs	138.95
229475	50 ml, Conical-Bottom	12,000 Xg	10/Bag, 100/Cs	153.95
229479	50 ml, Self-Standing	6,000 Xg	10/Bag, 100/Cs	186.95

PROGENE® CELL SCRAPERS & LIFTERS



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE



CAT. #	DESCRIPTION	PACKAGING	PRICE
71-800-25	Cell Scraper, 23 cm with 20 mm Blade	1/Pk, 100/Cs	119.95
229305	Cell Lifter, Flat & Narrow Blade	1/Pk, 100/Cs	123.95
229306	Cell Lifter, Flat & J-Hook Blade	1/Pk, 100/Cs	123.95

PROGENE® MEDIA STORAGE BOTTLES, **STERILE**

- For reduced gas permeability, 229788 is manufactured from medical grade virgin PETG
- Ideal for media prep, storage and stem cell, mammalian and bacterial culture applications



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
229781	150 ml Polystyrene, Round	1/Bag, 24/Cs	104.95
229782	250 ml Polystyrene, Round	1/Bag, 24/Cs	98.95
229784	500 ml Polystyrene, Round	1/Bag, 24/Cs	124.95
229785	1000 ml Polystyrene, Round	1/Bag, 24/Cs	232.95
229788	500 ml PETG, Square	1/Bag, 24/Cs	182.95

CELL CULTURE AND MICROPLATES

BRANDTECH BRANDPLATES® 96-WELL MICROPLATES



- **pureGrade™: Non-treated, medium binding.** suitable for homogeneous assays, screening and storage. Choose **pureGrade™ S** for a sterile surface.
- **immunoGrade™: High-binding. Ideal for the immobilization of IgG.** High binding capacity for molecules with hydrophilic and hydrophobic regions. The best choice for standard ELISA applications.
- **hydroGrade™: Hydrophilic, ideal for immobilization of hydrophilic molecules.** Plates with a hydroGrade™ surface can be used in liquid phase assays with hydrophobic components.
- **lipoGrade™: Hydrophobic, ideal for immobilization of hydrophobic molecules.** Microplates with lipoGrade™ surfaces can be used in liquid phase assays with hydrophilic components.
- **cellGrade™: Ideal for the cultivation of adherent cell lines.** Chemical groups on the surface of the plate bind to serum compounds, stimulating the growth of immobilized cells.

CAT. #	DESCRIPTION	VOL./WELL	PACKAGING	PRICE
--------	-------------	-----------	-----------	-------

pureGrade™ - Non Sterile (20 stacks of 5 plates)

34-781600	Clear, U-bottom	330 µl	100/Case	214.95
34-781601	Clear, V-bottom	360 µl	100/Case	214.95
34-781602	Clear, F-bottom	350 µl	100/Case	214.95
34-781603	Clear, C-bottom	350 µl	100/Case	214.95
34-781605	Solid White, F-bottom	350 µl	100/Case	369.95
34-781608	Solid Black, F-bottom	350 µl	100/Case	369.95

pureGrade™ S - Sterile

34-781660	Clear, U-bottom	330 µl	1/Pack, 50/Case	232.95
34-781661	Clear, V-bottom	360 µl	1/Pack, 50/Case	232.95
34-781662	Clear, F-bottom	350 µl	1/Pack, 50/Case	232.95
34-781663	Clear, C-bottom	350 µl	1/Pack, 50/Case	232.95
34-781665	Solid White, F-bottom	350 µl	1/Pack, 50/Case	352.95
34-781668	Solid Black, F-bottom	350 µl	1/Pack, 50/Case	352.95

immunoGrade™ - Non Sterile (20 stacks of 5 plates)

34-781720	Clear, U-bottom	330 µl	100/Case	352.95
34-781721	Clear, V-bottom	360 µl	100/Case	352.95
34-781722	Clear, F-bottom	350 µl	100/Case	352.95
34-781723	Clear, C-bottom	350 µl	100/Case	352.95
34-781724	Solid White U-bottom	330 µl	100/Case	541.95

hydroGrade™ - Non Sterile (20 stacks of 5 plates)

34-781780	Clear, U-bottom	330 µl	100/Case	390.95
34-781782	Clear, F-bottom	350 µl	100/Case	390.95

lipoGrade™ - Non Sterile (20 stacks of 5 plates)

34-781840	Clear, U-bottom	330 µl	100/Case	232.95
34-781841	Clear, V-bottom	360 µl	100/Case	232.95
34-781842	Clear, F-bottom	350 µl	100/Case	232.95

cellGrade™ - Sterile

34-781960	Clear, U-bottom	330 µl	1/Pack, 50/Case	264.95
34-781961	Clear, V-bottom	360 µl	1/Pack, 50/Case	264.95
34-781962	Clear, F-bottom	350 µl	1/Pack, 50/Case	264.95
34-781965	Solid White, F-bottom	350 µl	1/Pack, 50/Case	264.95
34-781968	Solid Black, F-bottom	350 µl	1/Pack, 50/Case	264.95

PLATESEAL™ ELISA FILM


- Ideal for ELISA, incubation and general use
- Clear, polypropylene film for sealing both polystyrene and polypropylene plates
- Excellent seal that removes easily and cleanly
- Temperature Range: -40°C to +120°C

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PET-100	Incubation Film	100/Pk	64.95
PS-PET-100S	Incubation Film, Sterile	100/Pk	91.95

PlateSeal™

CERTIFIED

RNase, DNase & Pyrogen-Free



HISTOLOGY, PATHOLOGY AND CUVETTES

ULTIDENTBRAND MICROSCOPE SLIDES & COVER GLASS

- All slides feature 45° clipped corners
- Our cover glass is made of first quality "Hydrolytic Class" glass
- Positive Charged Slides: Electrostatically attracts frozen tissue sections and cytology preparations



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Microscope Slides 3" X 1"			
170-7107P	Plain Slides, Beveled Edges and Safety Corners	72/Pk	9.95
170-7107A-S	Frosted 1-End & 2-Sides, Ground Edges and Safety Corners	72/Pk	10.95
Positive Charged Slides			
170-8105-B+	Color-Frost + Charged Slides, Blue	72/Pk	27.95
170-8105-G+	Color-Frost + Charged Slides, Green	72/Pk	27.95
170-8105-P+	Color-Frost + Charged Slides, Pink	72/Pk	27.95
170-8105-W+	Color-Frost + Charged Slides, White	72/Pk	27.95
170-8105-Y+	Color-Frost + Charged Slides, Yellow	72/Pk	27.95
Cover Glass - #1			
170-C1818	18 x 18 mm	1-Oz/Pk	6.95
170-C2222	22 x 22 mm	1-Oz/Pk	6.95
170-C2240	22 x 40 mm	1-Oz/Pk	6.95
170-C2250	22 x 50 mm	1-Oz/Pk	6.95
170-C2424	24 x 24 mm	1-Oz/Pk	6.95
170-C2432	24 x 32 mm	1-Oz/Pk	6.95
170-C2440	24 x 40 mm	1-Oz/Pk	6.95
170-C2450	24 x 50 mm	1-Oz/Pk	6.95
170-C2460	24 x 60 mm	1-Oz/Pk	6.95

ULTIDENTBRAND SLIDE ACCESSORIES

- Choose from cardboard or polypropylene construction slide mailers
- Cardboard mailers feature grooves to facilitate slide removal
- Polypropylene mailers close securely with a snap



Slide Mailer Polypropylene




Slide Mailer Cardboard



Rotary Slide Dispenser



Slide Portfolio



Slide Holders Polyethylene



Staining Tray Set

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Slide Mailers			
170-60001001	1-Slide Mailer Cardboard	100/Box	19.95
170-60001002	2-Slide Mailer Cardboard	50/Box	21.95
170-60001003	3-Slide Mailer Cardboard	50/Box	22.95
170-60001004	4-Slide Mailer Cardboard	50/Box	26.95
170-31011001	1-Slide Mailer Polypropylene	100/Box	19.95
170-31011002	2-Slide Mailer Polypropylene	100/Box	22.95
170-31011003	3-Slide Mailer Polypropylene	100/Box	23.95
170-31011103	3-Slide Mailer with Screw-Cap, PP	100/Box	26.95
170-31011005	5-Slide Mailer with Snap-Cap, PP	100/Box	31.95
Slide Portfolios			
170-60001020	20-Slide Portfolio Cardboard	5/Pk	27.95
Rotary Slide Dispenser			
170-31011000	Rotary Slide Dispenser	Each	31.95
Staining Jars			
200-4500	5-Place Staining Jar with Lid ABS	Each	5.95
200-4501	10-Place Staining Jar with Lid ABS	Each	7.95
Slide Holders, Polyethylene			
33-15981-A	40-Slide Holder	Each	61.95
33-15981-B	78-Slide Holder	Each	75.95
33-15981-C	15-Slide Holder	Each	56.95
Slide Staining Tray Set With Lid			
33-15951A	8-Slide Staining Tray Set	4 Base Trays + 1 Lid	57.95

SIMPORT EASYDIP™ SLIDE STAINING SYSTEM

- Easy to clean and no metals to corrode - Ideal for special stains, frozen sections
- Molded from Acetal - Resistant to alcohol and xylene (but not phenol, iodine or ferric chloride)
- System made up of two components: a square staining jar and a 12-position vertical slide rack
- Wide stable base with recessed interior allowing for a smaller reagent volume of only 80 ml



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-M900-12AS	Staining Jars, Assorted Colors	5/Pk	94.95
17-M905-12DGY	12-Place Slide Staining Rack	6/Pk	80.95
17-M906	6-Place Aluminum Holder	Each	90.95

ULTIDENTBRAND SLIDE STORAGE BOXES

- Durable ABS construction
- Molded, rigid boxes securely hold slides in grooved slots
- Boxes feature padded interiors with numbered slots and slide inventory sheets
- 100-Place deluxe boxes feature a rust-resistant, nickel plated clasp and pin hinges




CAT. #	DESCRIPTION	PACKAGING	PRICE
20-99418-_*	100-Place Deluxe with Hinged lid, cork	Each	12.95

Color Options: B (Blue), - BK (Black), - G (Green), - R (Red), - W (White), - Y (Yellow)

SIMPORT HISTOSETTE® I CASSETTES

- Choose from cassettes with a writing angle of 30° or 45°
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage
- One-piece cassettes snap apart and can be opened or closed as often as necessary
- Cassettes with 45° angles make them suitable for use with cassette labeling instruments



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-M490_*	Tissue Cassettes, 30° Writing Area	500/Box, 1500/Cs	174.95
17-M498_*	Tissue Cassettes, 45° Writing Area	500/Box, 1500/Cs	174.95
17-M491_*	Biopsy Cassettes, 30° Writing Area	500/Box, 1500/Cs	174.95
17-M499_*	Biopsy Cassettes, 45° Writing Area	500/Box, 1500/Cs	174.95

Color Options: -2 (White), -3 (Pink), -4 (Green), -5 (Yellow), -6 (Blue), -7 (Peach), -8 (Tan), -9 (Grey), -10 (Lilac), -11 (Orange), -12 (Aqua)

SIMPORT CASSETTE STORAGE DRAWERS

- Store up to 165 embedding rings or 250 cassettes per drawer
- Choose from High Impact Polystyrene (HIPS) or heavy-duty waterproof cardboard construction
- Each drawer measures 40.5 cm x 23 cm x 5.1 cm



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-M495-6	HIPS Storage Drawers	6/Cs	196.95
17-M495-7	Cardboard Storage Drawers	12/Cs	179.95

SIMPORT FORMALIN-FILLED CONTAINERS

- Designed for collection, transport and storage of histology specimens
- Half filled with 10% Neutral Buffered Formalin - Penetrates quickly, but fixes slowly
- Leak-resistant jars are translucent and specimens can be viewed without having to open the lid



CAT. #	DESCRIPTION	PACKAGING	PRICE
17-M961-20FMA	20 ml Container with Formalin	24/Tray, 96/Cs	90.95
17-M961-40FMA	40 ml Container with Formalin	24/Tray, 96/Cs	50.95
17-M961-60FMA	60 ml Container with Formalin	24/Tray, 96/Cs	55.95
17-M961-90FMA	90 ml Container with Formalin	24/Tray, 96/Cs	70.95
17-M961-120FMA	120 ml Container with Formalin	24/Tray, 96/Cs	78.95

HISTOLOGY, PATHOLOGY AND CUVETTES

ULTIDENTBRAND SPECIMEN TRANSPORT BAGS

- Zipper closure allows repeated access to specimen
- Constructed of durable, leak-resistant, 2-mil-thick polyethylene
- Bags are printed with biohazard symbol and include a document pouch



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
11-4918	6" x 9" Bag, Clear - Coded Red	1000/Cs	96.95
11-4950A	6" x 9" Stat Bag, Red	1000/Cs	80.95
11-4950B	6" x 9" Bag, Clear - Coded Yellow	1000/Cs	80.95

TERUMO SYRINGES & NEEDLES, STERILE



TERUMO

CAT. #	DESCRIPTION	PACKAGING	PRICE
10-SS-01T	1cc Tuberculin Syringes - No Needles	100/Pk	20.95
10-SS-03L	3cc Luer-Lock Syringes - No Needles	100/Pk	16.95
10-SS-05L	5cc Luer-Lock Syringes - No Needles	100/Pk	30.95
10-SS-10L	10cc Luer-Lock Syringes - No Needles	100/Pk	35.95
10-SS-20L	20cc Luer-Lock Syringes - No Needles	50/Pk	35.95
10-SS-30L	30cc Luer-Lock Syringes - No Needles	25/Pk	35.95
10-SS-60L	60cc Luer-Lock Syringes - No Needles	25/Pk	45.95

Hypodermic Needles

10-NN-1925R	19 Gauge x 1	100/Pk	11.95
10-NN-2025R	20 Gauge x 1	100/Pk	11.95
10-NN-2038R	20 Gauge x 1 ½	100/Pk	11.95
10-NN-2125R	21 Gauge x 1	100/Pk	11.95
10-NN-2138R	21 Gauge x 1 ½	100/Pk	11.95
10-NN-2225R	22 Gauge x 1	100/Pk	11.95
10-NN-2238R	22 Gauge x 1 ½	100/Pk	11.95
10-NN-2338R	23 Gauge x 1 ½	100/Pk	11.95
10-NN-2538R	25 Gauge x 1 ½	100/Pk	11.95
10-NN-2713R	27 Gauge x ½	100/Pk	11.95

PROGENE® INOCULATING LOOPS & SPREADERS, STERILE

- Made from polypropylene
- Color coded for easy identification
- L-Shaped Cell Spreader is ideal for evenly spreading bacterial cultures

PROGENE

**BUY 3,
GET 1
FREE**

equal or lower
price



CAT. #	DESCRIPTION	PACKAGING	PRICE
229601	1 µl Inoculating Loop, White	20/Pouch, 2000/Cs	177.95
229602	1 µl Inoculating Loop, White	1/Pouch, 600/Cs	89.95
229610	10 µl Inoculating Loop, Blue	20/Pouch, 2000/Cs	177.95
229612	10 µl Inoculating Loop, Blue	1/Pouch, 600/Cs	89.95
229615	Inoculating Needle, Yellow	20/Pouch, 2000/Cs	177.95
229616	L-Shaped Spreader, Green	1/Pouch, 100/Cs	34.95
229617	L-Shaped Spreader, Green	10/Pouch, 500/Cs	74.95

ULTIDENTBRAND RIGID INOCULATION LOOPS, STERILE

- Gamma sterilized - Made from medical grade of polystyrene
- Strong, rigid feel to provide comfortable, stable and uniform streaking
- Loops have an integrated needle on the other end



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
11-2811	1 µl loop with needle, Blue	20/Pk, 500/Cs	71.95
11-2815	10 µl loop with needle, Yellow	20/Pk, 500/Cs	71.95

PARAGON SCALPELS, STERILE

- All blades are sterilized by irradiation
- Single-use, they are individually-wrapped in foil pouches
- Blades 10-15 are compatible with scalpel handle #3 and 20-24 with scalpel handle #4



CAT. #	DESCRIPTION	PACKAGING	PRICE
02-90010-*	Stainless Steel Blades	1/Pouch, 100/Box	39.95
02-90000-*	Disposable Scalpels with Plastic Handle	1/Pouch, 10/Box	18.95

Size options: -10, -11, -12, -15, -20, -21, -22, -23, -24

HISTOLOGY, PATHOLOGY AND CUVETTES

BRANDTECH® SPECTROPHOTOMETER CUVETTES



**NEW
PRODUCT**

- Choose between polystyrene or UV-grade methacrylate (PMMA) construction
- Cells feature 10 mm optical paths and feature 2 recessed windows to prevent scratching
- Cuvettes are automatically sorted by a robot so that each package only contains cuvettes from the same mold cavity, ensuring the lowest possible cuvette-to-cuvette variance
- Ultra-micro cuvettes are specifically engineered for the analysis of DNA, RNA, and proteins and are available individually wrapped and Certified DNA-, DNase-, RNase-free - They must be selected with a window height that matches the beam height of your spectrophotometer



CAT. #	DESCRIPTION	PACKAGING	PRICE
Ultra-micro UV, 220-900 nm			
34-759250	70-850 µl, 8.5 mm Window Height	100/Tray	124.95
34-759215	70-850 µl, 8.5 mm Window Height, Sterile	1/Pouch, 100/Pack	199.95
34-759230	70-850 µl, 15 mm Window Height	100/Tray	124.95
34-759235	70-850 µl, 15 mm Window Height, Sterile	1/Pouch, 100/Pack	199.95
Polystyrene, 340-900 nm			
34-759075D	1.5 ml Semi-micro	100/Tray, 500/Cs	141.95
34-759070D	4.5 ml Macro	100/Tray, 500/Cs	129.95
Methacrylate (PMMA), 300-900 nm			
34-759080D	4.5 ml Macro	100/Tray, 500/Cs	186.95
34-759170	4.5 ml Macro, UV-Range	100/Tray	134.95
34-759085D	3.0 ml Semi-micro	100/Tray, 500/Cs	161.95
34-759165	3.0 ml Semi-micro, UV-Range	100/Tray	124.95

ULTIDENTBRAND PEG TUBE RACK

- Use as a drying rack for glass plates or upside-down tubes
- Fiberglass reinforced polypropylene racks are ideal for bench work or tube storage



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-24312A	Peg Rack for 10-13 mm Tubes, White	2/Pk	34.95
33-24312B	Peg Rack for 10-13 mm Tubes, Blue	2/Pk	34.95

ULTIDENTBRAND CUVETTES

- Choose from optical polystyrene or UV grade PMMA - Engineered for high reproducibility - Packaged in dust-free PE foam trays of 100 cuvettes



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Polystyrene (340-800 nm)			
11-112117	1.5 ml Semi-Micro	100/Tray, 500/Cs	82.95
11-111117	4.5 ml Macro	100/Tray, 500/Cs	78.95
11-111137	4.5 ml Macro 4-Window	100/Tray, 500/Cs	79.95
Methacrylate (280-800 nm)			
11-112147	1.5 ml Semi-Micro	100/Tray, 500/Cs	121.95
11-111147	4.5 ml Macro	100/Tray, 500/Cs	126.95
11-111157	4.5 ml Macro 4-Window	100/Tray, 500/Cs	115.95

PROGENE® BACTERIAL PETRI DISHES, STERILE

**PROMO
PRICE**

- Sterilized by gamma irradiation and Certified Non-pyrogenic
- Precision molded from crystalline polystyrene with super-flat surface
- Venting ribs built into the underside of the lid prevent buildup of condensation
- Stackable



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
229638	35 x 10 mm Bacterial Petri Dishes	10/Bag, 960/Cs	256.95
229663	60 x 15 mm Bacterial Petri Dishes with Grip Ring	10/Bag, 500/Cs	201.95
85-7453	100 x 15 mm Bacterial Petri Dishes	10/Bag, 500/Cs	92.95

ULTIDENTBRAND GLASS DROPPERS, 9 CM

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
21-1440	Glass Droppers	144/Pk	19.95



GLOVES AND APPAREL

AURELIA® SONIC POWDER-FREE NITRILE GLOVES

aurelia 2.2 MIL

OPTIMAL STRETCHABILITY



BUY 6,
GET 4
FREE
equal or lower
price

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93775	Sonic® Stretch Nitrile, X-Small	Cobalt Blue	300/Box	24.95
262-93776	Sonic® Stretch Nitrile, Small	Cobalt Blue	300/Box	24.95
262-93777	Sonic® Stretch Nitrile, Medium	Cobalt Blue	300/Box	24.95
262-93778	Sonic® Stretch Nitrile, Large	Cobalt Blue	300/Box	24.95
262-93779	Sonic® Stretch Nitrile, X-Large	Cobalt Blue	270/Box	24.95

AURELIA® TRANSFORM™ 100 POWDER FREE NITRILE GLOVES

aurelia 3.2 MIL

OPTIMAL STRETCHABILITY



BUY 8,
GET 2
FREE
equal or lower
price

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-9889A5	Transform™ 100 Stretch Nitrile, X-Small	Ice Blue	100/Box	7.95
262-9889A6	Transform™ 100 Stretch Nitrile, Small	Ice Blue	100/Box	7.95
262-9889A7	Transform™ 100 Stretch Nitrile, Medium	Ice Blue	100/Box	7.95
262-9889A8	Transform™ 100 Stretch Nitrile, Large	Ice Blue	100/Box	7.95
262-9889A9	Transform™ 100 Stretch Nitrile, X-Large	Ice Blue	100/Box	7.95

AURELIA® AMAZING™ POWDER-FREE NITRILE GLOVES

aurelia 2.0 MIL

UNIQUE STRETCHY FORMULATION



CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-92885	Amazing™ Nitrile, X-Small	Purple	300/Box	24.95
262-92886	Amazing™ Nitrile, Small	Purple	300/Box	24.95
262-92887	Amazing™ Nitrile, Medium	Purple	300/Box	24.95
262-92888	Amazing™ Nitrile, Large	Purple	300/Box	24.95
262-92889	Amazing™ Nitrile, X-Large	Purple	270/Box	24.95

AURELIA® PERFORM™ POWDER-FREE NITRILE GLOVES

aurelia 2.5 MIL

OPTIMAL STRETCHABILITY



CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-92395	Perform™ Stretch Nitrile, X-Small	Teal	200/Box	17.95
262-92396	Perform™ Stretch Nitrile, Small	Teal	200/Box	17.95
262-92397	Perform™ Stretch Nitrile, Medium	Teal	200/Box	17.95
262-92398	Perform™ Stretch Nitrile, Large	Teal	200/Box	17.95
262-92399	Perform™ Stretch Nitrile, X-Large	Teal	200/Box	17.95

AURELIA® PROTÉGÉ® POWDER-FREE NITRILE GLOVES

aurelia 4 MIL

ULTRA SOFT AND STRETCHY



CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93995	Protégé™ Nitrile, X-Small	Blue	100/Box	7.95
262-93996	Protégé™ Nitrile, Small	Blue	100/Box	7.95
262-93997	Protégé™ Nitrile, Medium	Blue	100/Box	7.95
262-93998	Protégé™ Nitrile, Large	Blue	100/Box	7.95
262-93999	Protégé™ Nitrile, X-Large	Blue	100/Box	7.95

AURELIA® RENEW POWDER-FREE NITRILE GLOVES

Renew examination gloves biodegrade 10X
faster than standard nitrile gloves

aurelia 2.2 MIL

BIODEGRADABLE



CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93285	Renew® Biodegradable Nitrile, X-Small	Green Apple	300/Box	21.95
262-93286	Renew® Biodegradable Nitrile, Small	Green Apple	300/Box	21.95
262-93287	Renew® Biodegradable Nitrile, Medium	Green Apple	300/Box	21.95
262-93288	Renew® Biodegradable Nitrile, Large	Green Apple	300/Box	21.95
262-93289	Renew® Biodegradable Nitrile, X-Large	Green Apple	300/Box	21.95

GLOVES AND APPAREL

AURELIA® ROBUST® POWDER-FREE NITRILE GLOVES

aurelia 5 MIL



**BUY 8,
GET 2
FREE**
equal or lower
price

OPTIMAL STRETCHABILITY

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93895	Robust™ Thick Nitrile, X-Small	Blue	100/Box	8.95
262-93896	Robust™ Thick Nitrile, Small	Blue	100/Box	8.95
262-93897	Robust™ Thick Nitrile, Medium	Blue	100/Box	8.95
262-93898	Robust™ Thick Nitrile, Large	Blue	100/Box	8.95
262-93899	Robust™ Thick Nitrile, X-Large	Blue	100/Box	8.95

AURELIA® BOLD® POWDER-FREE NITRILE GLOVES

aurelia 5 MIL



**BUY 8,
GET 2
FREE**
equal or lower
price

CHEMOTHERAPY-TESTED

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-73995	Bold® Thick Nitrile, X-Small	Black	100/Box	9.95
262-73996	Bold® Thick Nitrile, Small	Black	100/Box	9.95
262-73997	Bold® Thick Nitrile, Medium	Black	100/Box	9.95
262-73998	Bold® Thick Nitrile, Large	Black	100/Box	9.95
262-73999	Bold® Thick Nitrile, X-Large	Black	100/Box	9.95
262-73990	Bold® Thick Nitrile, 2X-Large	Black	100/Box	9.95

AURELIA® ROBUST PLUS® EXTENDED-CUFF 12"

NITRILE POWDER-FREE GLOVES

aurelia 4.5 MIL



EXTENDED CUFF

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-63886	Robust™ Plus Nitrile, Small	Blue	100/Box	11.95
262-63887	Robust™ Plus Nitrile, Medium	Blue	100/Box	11.95
262-63888	Robust™ Plus Nitrile, Large	Blue	100/Box	11.95
262-63889	Robust™ Plus Nitrile, X-Large	Blue	100/Box	11.95

AURELIA® VIBRANT™ POWDER-FREE LATEX GLOVES

aurelia



UNIQUE STRETCHY FORMULATION

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-9822-XS	Vibrant™ Latex, X-Small	Beige	100/Box	9.95
262-9822-S	Vibrant™ Latex, Small	Beige	100/Box	9.95
262-9822-M	Vibrant™ Latex, Medium	Beige	100/Box	9.95
262-9822-L	Vibrant™ Latex, Large	Beige	100/Box	9.95
262-9822-XL	Vibrant™ Latex, X-Large	Beige	100/Box	9.95

PROGENE® SURGICAL MASK ASTM LEVEL 2

- Pleated 3 layer design
- Latex and fiberglass free
- Bacterial Filtration Efficiency (BFE) ≥ 99%
- Sub-micron particulate Efficiency @ 0.1 µm (PFE) ≥ 99. Resistance to penetration by synthetic blood @ 120 mmHg



PROGENE

**BUY 3,
GET 1
FREE**
equal or lower
price

CAT. #	DESCRIPTION	PACKAGING	PRICE
WLM2002	Surgical Masks Level 2, Blue	50/Box	5.95

AURELIA® SURGICAL MASK



aurelia

CAT. #	DESCRIPTION	PACKAGING	PRICE
262-2110	Aurelia Surgical Masks Level 1, Blue	50/Box	4.95
262-2120	Aurelia Surgical Masks Level 2, Blue	50/Box	5.95
262-2130	Aurelia Surgical Masks Level 3, Blue	50/Box	6.95

GLOVES AND APPAREL

ULTIDENTBRAND SMS COVERALLS

- Made of three-layer SMS (spunbond/meltblown/spun-bond) material - Offers some splash resistance while keeping users cool and comfortable



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
20-428-L	Coverall, White, Large	1/Bag, 25/Cs	99.95
20-428-XL	Coverall, White, X-Large	1/Bag, 25/Cs	99.95
20-428-XXL	Coverall, White, 2X-Large	1/Bag, 25/Cs	109.95

ULTIDENTBRAND MICROPOROUS LAB COATS

- Latex free, microporous fabric is made of two layers of spunbonded polypropylene with a thin microporous film layer
- Breathable, fluid & abrasion resistant
- 2-Pockets & collar



UltiDent
Scientific

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-6023	Microporous Lab Coat XL, 99 x 140cm	White	50/Cs	129.95
262-6024	Microporous Lab Coat XXL, 104 x 146cm	White	50/Cs	129.95

ULTIDENTBRAND DISPOSABLE LAB COATS

- Disposable latex free lab coats with knit collar and cuffs protect lab personnel from fluids and cross contamination
- Made of non-woven material: Spunbonded Meltblown: Spunbond (SMS) Polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
77-9027	SMS Lab Coat, Medium	Blue	50/Cs	109.95
77-9028	SMS Lab Coat, Large	Blue	50/Cs	109.95
77-9029	SMS Lab Coat, X-Large	Blue	50/Cs	109.95

ULTIDENTBRAND DISPOSABLE ISOLATION GOWNS

- Made of spunbond polypropylene and are latex-free
- Isolation gowns are made for a comfortable fit with elastic cuffs, waist and neckties



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
77-900-L	Gowns, Yellow, Large	10/Bag, 50/Cs	69.95

ULTIDENTBRAND LATEX-FREE SHOE COVERS

- Non-conductive and non-skid with sewn seams provide increased tear strength



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
32-2132	Shoe Covers Non-Skid, Blue, Reg.	300/Box	39.95
32-2134	Shoe Covers Non-Skid, Blue, XL	300/Box	44.95

ULTIDENTBRAND MICROPOROUS FLUID RESISTANT SHOE COVERS

- Latex free, microporous fabric is made of two layers of spunbonded polypropylene with a thin microporous film layer



UltiDent
Scientific

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
77-2618	Microporous Shoe Covers, 42x19 cm	White	300/Cs	129.95

GLOVES AND APPAREL

ULTIDENTBRAND PROTECTIVE SLEEVES

- Polyethylene coated with elastic cuffs
- Water resistant, but also allows moisture to evaporate



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
02-20300	16" Polyethylene, Blue	100/Bag, 2000/Cs	99.95

ULTIDENTBRAND MICROPOROUS FLUID RESISTANT SLEEVES

- Latex free, microporous fabric is made of two layers of spunbonded polypropylene with a thin microporous film layer
- 18" length with elastic at both ends



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
77-2407	Microporous Sleeves, White	500/cs	139.95

ULTIDENTBRAND BOUFFANT CAPS & BEARD COVER

- Full hair coverage
- Lightweight and breathable
- Made of light weight spunbond polypropylene
- Latex and lint free



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
20-4402-W	Bouffant Caps White, 21"	500/Cs	39.95
20-4403-W	Bouffant Caps White, 24"	500/Cs	40.95
20-P21-W	Pleated Bouffant Caps White, 21"	500/Cs	34.95
20-4400	Beard Covers, 18", White	1000/Cs	69.95

ULTIDENTBRAND FACE SHIELDS

- Wrap around design offers full facial coverage
- Top guard foam headband provides fluid protection & relieves pressure
- Excellent visibility thanks to double sided anti-fog
- Comfortable, light weight, anti-static treated
- Can be used with face mask and/or glasses
- FDA Certified, meets European(CE)
- Face shields measure 34 x 22 cm



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
262-1050	UltidentBrand Face Shield	40/Cs	49.95

PROSURFACE + DISINFECTANT SPRAY AND WIPES

- ProSurface+ eliminates 48 pathogens, including bacteria, viruses and fungi, without damaging surfaces. Powerful one-step cleaning and disinfection in a best-in-class contact time of just one minute.



Medicom

CAT. #	DESCRIPTION	PACKAGING	PRICE
02-40060C	ProSurface+ Wipes, 6" x 6.75"	160/Tub	12.95
02-40061C	ProSurface+ Wipes, 9" x 12"	65/Tub	14.95
02-40066C	ProSurface+ Spray	32oz/Bottle	14.95

ULTIDENTBRAND TACKY MATS

- Helps remove dirt and bacteriacarrying particles from shoes, casters and cart wheels before entry into the clean environment
- Anti-bacterial low density polyethylene with water based acrylic adhesive



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
TM-2436B	24" X 36" Blue	30/Mat & 4 Mats/Cs	89.95
TM-3645B	36" X 45" Blue	30/Mat & 4 Mats/Cs	156.95
TM-3672B	36" X 72" Blue	30/Mat & 4 Mats/Cs	199.95



Featured Products

Offers valid from April 1 to August 31, 2023

30+
YEARS
OF GLOVE EXPERTISE

Say

HELLO

to Medicom Gloves



AssureTouch Power™

Nitrile Medical Examination Gloves

12" Extended cuffs
for powerful protection



02-9999A	Extra Small	02-9999D	Large
02-9999B	Small	02-9999E	Extra Large
02-9999C	Medium		

4 mils¹, 6.2 g | 100 gloves/box

\$14.95 /box of 100

AssureTouch Pure™

Nitrile Medical Examination Gloves

Strong protection that's
safe for skin



02-9970A	Extra Small	02-9970D	Large
02-9970B	Small	02-9970E	Extra Large
02-9970C	Medium		

2.76 mils¹, 3.2 g | 200 gloves/box

\$17.95 /box of 200

AssureTouch Protect™

Nitrile Medical Examination Gloves

The ideal balance between
protection and tactile sensitivity



02-9990A	Extra Small	02-9990D	Large
02-9990B	Small	02-9990E	Extra Large
02-9990C	Medium		

3 mils¹, 3.8 g | 100 gloves/box

\$8.95 /box of 100

Ocean Pacific Ninja™

Nitrile Medical Examination Gloves

Exceptional protection
when you need it most



02-OPXS-NI	Extra Small	02-OP-L-NI	Large
02-OPS-NI	Small	02-OPXL-NI	Extra Large
02-OPM-NI	Medium		

4.8 mils¹, 5.3 g | 200 gloves/box

\$18.95 /box of 200

¹Mils are estimated based on a range of measurements. Medicom reserves the right to modify pricing and/or promotions at any time and without prior notice.

©2023 AMD Medicom Inc. All rights reserved.

33 - To order, Call 1-800-463-4163 or 514-335-3433 | sales@ultident.com | www.ultident.com

PROGENE® POLYSTYRENE RESERVOIRS



CERTIFIED
RNase, DNase & Pyrogen-Free

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
Non-Sterile			
825-1-10	10 ml	100/Cs	39.95
825-1-25	25 ml	100/Cs	40.95
825-1-50	50 ml	100/Cs	42.95
825-10-100	100 ml	100/Cs	42.95
Sterile, Individually-Wrapped			
825-1-10-1S	10 ml	1/Pk, 100/Cs	95.95
825-1-25-1S	25 ml	1/Pk, 100/Cs	95.95
825-1-50-1S	50 ml	1/Pk, 100/Cs	105.95
825-10-100-1S	100 ml	1/Pk, 100/Cs	110.95
Sterile, Multi-Pack			
825-1-10-5S	10 ml	5/Pk, 200/Cs	126.95
825-1-25-5S	25 ml	5/Pk, 200/Cs	126.95
825-1-50-5S	50 ml	5/Pk, 200/Cs	139.95
825-10-100-5S	100 ml	5/Pk, 200/Cs	139.95

ULTIDENTBRAND 50 ML PVC RESERVOIRS

- Clear reservoirs allow maximum visibility
- "V" shaped well design allows maximum reagent recovery
- Drip free pour-off spouts at every corner for spill free pouring



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-20521A	50 ml PVC Reservoirs	100/Pk	93.95

PROGENE® ANTI-STATIC POLYSTYRENE WEIGH BOATS

- Use for liquid or solid samples during weighing procedures
- Smooth surface and rounded corners provide for easy removal of powdered and granular samples



PROGENE

CAT. #	DESCRIPTION	CAPACITY	DEMENSIONS(MM)	PACKAGING	PRICE
825-3620	Square, Small	7 ml	41 x 41 x 8 H	500/Pk	15.95
825-3621	Square, Medium	100 ml	78 x 78 x 25 H	500/Pk	37.95
825-3622	Square, Large	330 ml	140 x 140 x 25 H	500/Pk	77.95
825-20	Hexagonal, Small	20 ml	OD: 45, Depth: 10	500/Pk	16.95
825-50	Hexagonal, Medium	50 ml	OD: 75, Depth: 20	500/Pk	35.95
825-200	Hexagonal, Large	200 ml	OD: 120, Depth: 22	500/Pk	58.95
825-350	Hexagonal, X-Large	350 ml	OD: 140, Depth: 30	500/Pk	70.95

ULTIDENTBRAND WEIGHBUCKET™

- The weigh boat for liquids
- Made of polyethylene - Safety transports liquids, solids and powders
- Dimensions excluding rim, spout and handle: 53 x 33mm (D x H)



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-B6602-W	40 ml WeightBucket™ for Liquids	500/Pk	59.95

BIOBUCKET™ 5L COLLAPSIBLE LAB BUCKET

- Collapsible bucket saves storage space
- Attached handle simplifies transportation
- Molded base grip allows for easy pouring and cleaning



CAT. #	DESCRIPTION	PACKAGING	PRICE
48-B8300	BioBucket™ 5L Lab Bucket	Each	65.95

LABORATORY ESSENTIALS AND RACKS

ULTIDENTBRAND 80-WELL MICROTUBE RACKS

- 80-place polypropylene fraction collector style
- Holds tubes up to 2.0 ml, 11 x 26 mm maximum



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
24-R80-AS	80-Place Workstations, Assorted	5/Pk	24.99

ULTIDENTBRAND 50-PLACE FREEZER BOXES

- With hinged lid, Polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-10601-A	50-Place, Assorted	5/Pk	38.95
33-10601-N	50-Place, Natural	5/Pk	38.95

PROGENE® 96-PLACE PCR RACKS

- Lid clicks into place

PROGENE



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W96-RK	96-Place PCR Racks, Assorted	5/Pk	45.99

ULTIDENTBRAND REVERSIBLE RACKS


- One side holds 0.5ml and the other 1.5/2.0 mL tubes
- Bright colors enable organization



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W32-FLP-A	32-Place w/Clear Lid, Assorted	5/Pk	59.99
87-W96-FLP-A	96-Place w/Clear Lid, Assorted	5/Pk	67.99

ULTIDENTBRAND 4-WAY TUBE RACKS



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-29022-A	4Way, Assorted	5/Pk	44.99
33-29022-N	4Way, Natural	5/Pk	44.99

ULTIDENTBRAND REVOLVER WORKSTATION

- Rotate and lock mechanism




UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
5405-59	Revolver Workstation	4/Pk	49.99

ULTIDENTBRAND EXPANDING MCT STORAGE BOX

- Fits standard 2" freezer racks - Stackable and autoclavable
- Holds up to 64 x 1.5/2.0 ml MCT's and cryogenic vials up to 2.0 ml
- Facilitate sample processing - Expands for convenience and contracts for storage



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120018	64-Place Expanding Storage Box	Each	38.95

ULTIDENTBRAND PP FLOATING RACKS

- Floating racks keep microtubes afloat in a water bath or can be pushed into crushed ice



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-2132A	16-Place Square - 10.2 x 10.2 cm	4/Pk	41.95
33-2134A	20-Place Round - 9.8 cm dia.	4/Pk	41.95

LABORATORY ESSENTIALS AND RACKS

ULTIDENTBRAND RESERVOIR CHILLERS

- Keeps reagents at less than 5°C (41°F) for approximately 3 hours
- Non-toxic blend of water and bacteriostatic agent chills and freezes quickly
- Made from sonic welded, engineered polypropylene that ensures a leak free design



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120730	Chiller for 50/55 ml Reservoirs, Purple	2/Pk	108.95
33-120734	Chiller for 100 ml Reservoirs, Blue	2/Pk	108.95

ULTIDENTBRAND ALUMINUM WEIGHING BOATS & DISHES

- Sturdy and lightweight with flat-bottom
- Crimped wall for rigidity



UltiDent
Scientific

CAT. #	DESCRIPTION	CAPACITY	DIMENSIONS	PACKAGING	PRICE
33-14521A	Weigh Boat, Small	20 ml	Dia. 43mm H. 12 mm	100/Pk	11.95
33-14521B	Weigh Boat, Medium	39 ml	Dia. 57 mm H. 14 mm	100/Pk	12.95
33-14521C	Weigh Boat, Large	50 ml	Dia. 70 mm H. 16 mm	100/Pk	17.95
33-14522	Weigh Dish, Large	30 ml	Dia. 101 mm H. 8 mm	50/Pk	21.95

ULTIDENTBRAND 4-WAY FLIPPER RACK WITH TRAY

- Each flipper can be independently oriented
- 4-Way rack can be used either in a tray or independently



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
5400-39	4-Way Flipper Rack with Tray	4/Pk	61.95

ULTIDENTBRAND 6-WAY CUBE RACKS

- Use for 4 x 50 ml, 10 x 15 ml, 12 x 12 mm and 16 x 1.5/2.0 ml tubes
- Made of polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-29050-A	6-Way Cube Rack	5/Pk	69.95

ULTIDENTBRAND 15 & 50 ML TRANSFORMER CUBE

- Includes grids for 16 x 15 ml and 9 x 50 ml
- Made of polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-23271A	Transformer Cube	5/Pk	84.95

ULTIDENTBRAND CLINICAL 4 WAY RACK

- Holds all brands of 13- and 16-mm tubes and 15 ml centrifuge tubes
- Customize the rack to fit the number of wells required, infinite units can be linked together
- Made of polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120683	Clinical 4 Way Rack, Assorted	5/Pk	122.95

ULTIDENTBRAND CARDBOARD FREEZER BOXES

- Made of durable, moisture resistant coated cardboard
- All boxes are supplied complete with installed grids



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
21-BOX-2-81	2" Box with 81-Place Grid	12/Pk	24.95
21-BOX-2-100	2" Box with 100-Place Grid	12/Pk	24.95
21-BOX-3-81	3" Box with 81-Place Grid	12/Pk	24.95
21-BOX-3-100	3" Box with 100-Place Grid	12/Pk	24.95
21-BOX-PCR	144-Place Box for 0.2 ml PCR® Tubes	6/Pk	24.95
21-BOX-15ML	36-Place Box for 15 mL Centrifuge Tubes	6/Pk	34.95
21-BOX-50ML	16-Place Box for 50 mL Centrifuge Tubes	6/Pk	34.95

ULTIDENTBRAND 81-PLACE FREEZER BOXES POLYPROPYLENE



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-23255-N	81-Place, Natural	5/Pk	43.99
33-23255-A	81-Place, Assorted	5/Pk	43.99

ULTIDENTBRAND 100-PLACE FREEZER BOXES WITH HINGED LID, POLYPROPYLENE



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-1064	100-Place, Natural	5/Pk	45.99
33-1065	100-Place, Assorted	5/Pk	45.99

ULTIDENTBRAND STAINLESS RACKS FOR UPRIGHT FREEZERS

- Models listed below are designed for standard 5¼" x 5¼" x 2" standard boxes
- Drawer racks let you reach boxes at the back of the freezer - No need to remove the entire rack
- Also available: Racks for 3" boxes, 96 and 384-well plates, vertical racks for chest freezers and more.



UltiDent
Scientific

CAT. #	DESCRIPTION	DIMENSION	SIZE	PRICE
Standard Racks for 2" boxes				
2525-UF-332	9-Place	415×168×140	3 x 3	139.95
2525-UF-342	12-Place	415×225×140	3 x 4	175.95
2525-UF-352	15-Place	415×280×140	3 x 5	211.95
2525-UF-432	12-Place	550×168×140	4 x 3	182.95
2525-UF-442	16-Place	550×225×140	4 x 4	231.95
2525-UF-452	20-Place	550×280×140	4 x 5	242.95
2525-UF-462	24-Place	550×321×140	4 x 6	344.95
2525-UF-542	20-Place	679×225×140	5 x 4	325.95
2525-UF-552	25-Place	679×280×140	5 x 5	400.95

Drawer Racks for 2" Boxes

2525-UFD-332	9-Place	419×181×140	3 x 3	186.95
2525-UFD-342	12-Place	419×240×140	3 x 4	247.95
2525-UFD-442	16-Place	559×240×140	4 x 4	324.95
2525-UFD-472	28-Place	559×419×140	4 x 7	571.95
2525-UFD-552	25-Place	680×300×140	5 x 5	571.95

Drawer Racks for 15 & 50 ml Tubes

2525-UFD-LT15-1	15 ml x 80 Tubes, 1 Drawer	559×126×140	---	235.95
2525-UFD-LT15-2	15 ml x 160 Tubes, 2 Drawers	559×252×140	---	365.95
2525-UFD-LT50-1	50 ml x 39 Tubes, 1 Drawer	559×126×140	---	235.95
2525-UFD-LT50-2	50 ml x 78 Tubes, 2 Drawers	559×252×140	---	357.95

Racks for 100-Place Hinged-Top Freezer Boxes

2525-UFHT-23	6-Place	318×176×149	2x3	101.95
2525-UFHT-24	8-Place	318×235×149	2x4	130.95
2525-UFHT-25	10-Place	318×293×149	2x5	142.95
2525-UFHT-33	9-Place	476×176×149	3x3	166.95
2525-UFHT-34	12-Place	476×235×149	3x4	194.95
2525-UFHT-35	15-Place	476×293×149	3x5	227.95
2525-UFHT-44	16-Place	634×235×149	4x4	255.95

Other formats and sizes available.
Please contact your Territory Manager for pricing.

LABORATORY ESSENTIALS AND RACKS

COOLCADDY™ PCR WORKSTATION

- Keep your reagents and samples cool without using crushed ice
- Store overnight in freezer and it's ready to hold temperature between 2° and 4°C for up to 5 hours
- **Capacity:** 9 x 1.5 ml; 8 x 2.0 ml Cryogenic vials; 5 x 0.5 ml 1 x 96W PCR Plate & 1 x 50 ml reservoir



CAT. #	DESCRIPTION	PACKAGING	PRICE
48-R4015	CoolCaddy™ PCR workstation	Each	179.95

ULTIDENTBRAND 50-PLACE GRIP RACK

- Accepts all tubes from 12 to 18 mm in diameter
- Made of polypropylene, they are resistant to most laboratory chemicals and solvents
- Innovative design "grips" various diameters of culture and centrifuge tubes or secure handling



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
200-4100	Multi-Purpose Grip Rack	1/Pk	17.95

ULTIDENTBRAND 15 & 50 ML TUBE RACK

- Fits in standard freezer rack-144 x 144 x 127mm
- Use from -90° C min to +120° C max (autoclavable)
- Made of polypropylene



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W15ML	15 ml 36-Place, Green	2/Pk	40.95
87-W50ML	50 ml 16-Place, Blue	2/Pk	40.95

ULTIDENTBRAND UNIVERSAL TUBE RACK

- Made of chemically-resistant POM
- Racks can be linked together for additional capacity
- Each slot adapts to accommodate 12-30 mm diameter tube



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120184	Multi-Tube Rack	2/Pk	16.95

SSIbio® CRYOFREEZE® BOXES

- Freezer boxes are made of tough and durable polycarbonate
- Built-in drains and vents prevent pressure build-up



SSIbio

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
SSI-5700-19	25-Place for 1 to 2 ml Vials	Assorted	5/Pk	42.95
SSI-5720-19	81-Place for 1 to 2 ml Vials	Assorted	5/Pk	64.95
SSI-5730-19	81-Place for 3.6 to 4.5 ml Vials	Assorted	5/Pk	86.95
SSI-5750-19	100-Place for 1 to 2 ml Vials	Assorted	5/Pk	69.95

ULTIDENTBRAND SNAP-TOGETHER 15 & 50 ML RACK

- Designed for stable use in water baths
- Holds 30 x 15 and 20 x 50 ml tubes simultaneously
- Shipped flat, easy-to-assemble ABS racks snap together securely



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
200-4200	15 and 50 ml Tube Rack Orange	Each	32.95

ULTIDENTBRAND WIRELESS TUBE RACKS

- Racks sink in a water bath - Excellent temperature and chemical resistance
- Fully Autoclavable - Molded from reinforced nylon for strength and durability



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Half Size			
11-456445_	36 Place for 13 mm Tubes	Each	14.95
11-456545_	36 Place for 16 mm Tubes	Each	14.95
11-456645_	20 Place for 15 ml Centrifuge Tubes	Each	14.95
11-456845_	9 Place for 50 ml Centrifuge Tubes	Each	24.95
Full Size			
11-456440_	72 Place for 13 mm Tubes	Each	19.95
11-456540_	72 Place for 16 mm Tubes	Each	21.95
11-456640_	40 Place for 15 ml Centrifuge Tubes	Each	34.95
11-456840_	24 Place for 50 ml Centrifuge Tubes	Each	35.95

*Available in 7 colors: B (Blue), G (Green), M (Magenta), O (Orange), R (Red), Y (Yellow), W (White)

ULTIDENTBRAND ONE RACK™ TUBE RACKS

- Durable racks are molded in a single, continuous piece - no assembly required
- Choose from Delrin® fiberglass-reinforced polyoxymethylene (POM) or polypropylene construction
- POM is an engineered polymer developed by DuPont with exceptional chemical and heat resistance



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Half Rack Models - Polyoxymethylene (POM)			
33-27501_	36-Place for 13 mm Tubes 104 x 104 x 59 mm	8/Cs	156.95
33-27502_	36-Place for 16 mm Tubes 127 x 127 x 70 mm	8/Cs	163.95
33-27503_	20-Place for 20 mm Tubes 100 x 127 x 83 mm	8/Cs	178.95
33-27504_	16-Place for 25 mm Tubes 120 x 122 x 92 mm	8/Cs	178.95
33-27505_	9-Place for 30 mm Tubes 110 x 110 x 85 mm	8/Cs	178.95
Full Rack Models - Polyoxymethylene (POM)			
33-27511_	72-Place for 13 mm Tubes 104 x 202 x 59 mm	8/Cs	231.95
33-27512_	72-Place for 16 mm Tubes 127 x 250 x 70 mm	8/Cs	231.95
33-27513_	40-Place for 20 mm Tubes 100 x 252 x 83 mm	8/Cs	231.95
33-27514_	40-Place for 25 mm Tubes 120 x 300 x 92 mm	8/Cs	251.95
33-27515_	24-Place for 30 mm Tubes 110 x 282 x 85 mm	8/Cs	251.95
Full Rack Models - Polypropylene (PP)			
33-27551_	72-Place for 13 mm Tubes 104 x 202 x 59 mm	8/Cs	150.95
33-27552_	72-Place for 16 mm Tubes 127 x 250 x 70 mm	8/Cs	150.95

*Available in 4 colors: A (White), B (Blue), G (Green), D (Magenta)

ULTIDENTBRAND SILICONE LAB MATS

- Protect countertops from hot items up to 200°C (392°F)
- Chemical resistance: chemically inert and does not react with most chemicals
- Mats measure 35 cm x 60 cm, 2 mm thick - Autoclavable at 121° for 15 minutes



UltiDent
Scientific

CAT. #	DESCRIPTION	PRICE
33-120506	Side One: Solid Yellow - Side Two: Blue design on grey background	76.95
33-120507	Side One: Solid Purple - Side Two: Grey design on purple background	76.95
33-120599	Side One: Solid Blue - Side Two: Solid white with blue rim	76.95

ULTIDENTBRAND SILICONE LAB WORKSTATIONS

- Lip around the edge can hold approx. 250 ml of liquid
- Small footprint 30.5 x 37.9 x 0.2 cm - Rolls up compactly for storage
- Reversible, Side One: combination gridded and smooth surface, Side Two: smooth surface only



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120747	Silicone Lab Mat 30.5 cm x 37.9 cm; Blue	Each	63.95
33-120748	Silicone Lab Mat 30.5 cm x 37.9 cm; Purple	Each	63.95

ULTIDENTBRAND LABORATORY TIMERS

- **33-24490D:** Four independent channels with visual and audible Alarms - Spring Clip and Magnetic Back
- **33-120365:** Count up or down to 24 hours 59 minutes - 4 timer channels - Stands up or attaches to a magnetic surface - Whiteboard and marking pen included
- **33-24490M:** Three Channels to count down - Start independently or simultaneously
- **33-24670:** Four Channel Timer/Clock - Toggle between channels and clock with top buttons - Spring clip with magnetic back - Timing Capacity: 99 Hrs.




UltiDent
Scientific

CAT. #	DESCRIPTION	PRICE
33-24490D	4-Channel Timer/Stopwatch	35.95
33-24490M	Triple Timer	46.95
33-120365	4-Channel Timer With Whiteboard & Pen	41.95
33-24670	4-Channel Timer/Clock	35.95
33-6594	Hand Tally Counter	18.95

ULTIDENTBRAND MEGA RACKS

- High-capacity storage for your 13-16 mm blood tubes
- Choose from single (120 Place) or double (240 Place) formats
- Fully assembled, polypropylene tube racks are fully autoclavable
- Alphanumeric identification for quick location and retrieval of tubes




UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-24323-W	120-Place, White	Each	57.95
33-120114-W	240-Place, White	Each	112.95

PROGENE® BENCHTOP COOLERS

- Maintains a sample temperature of less than 4°C for almost 4 hours (with lid on rack) when at room temperature.
- Color change plastic enables you to see how cold your samples are.
- Color change racks have a dramatic change to a contrasting color at 7°C.



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W24-MCT	24-Place MCT Cooler with Lid	Each	85.95
87-W96-RG	64-Place RotorGene® Cooler with Lid	2/Pk	105.95
87-W96-PCR	96-Place PCR Cooler with Lid	2/Pk	100.95

ULTIDENTBRAND MINI COOLERS

- Polycarbonate construction with rubber non-slip feet
- 0.5 ml tube adapters supplied as standard with the
- 0.5 - 2.0 mL tube size versions
- Ideal for keeping PCR® reagents, enzymes, and samples cool while working at the bench



UltiDent
Scientific

CAT. #	DESCRIPTION	TEMPERATURE	PRICE
33-120070	12 x 0.5 - 2 ml Tubes	0° C x 4 h	213.95
33-120071	32 x 0.5 - 2 ml Tubes	0° C x 5 h	359.95
33-120073	12 x 0.5 - 2 ml Tubes	-20° C x 2.5 h	213.95
33-120074	32 x 0.5 - 2 ml Tubes	-20° C x 3.5 h	359.95
33-120076	96 x 0.2 ml PCR® Tubes	-20° C x 3.5 h	392.95

SSIbio® ISOFREEZE® MCT FLIPPER® RACK

- Maintains a temperature between -10° C and -20° C for up to 3 hours
- One side holds 20 x 1.5/2.0 ml tubes, and the other side 20 x 0.5 ml tubes



SSIbio

CAT. #	DESCRIPTION	PACKAGING	PRICE
5610-43	20-Place Freezer Rack	Each	111.95

LABORATORY ESSENTIALS AND RACKS

KIMTECH SCIENCE® WIPERS



KIMTECH

CAT. #	DESCRIPTION	PACKAGING	PRICE
14-05511	Precision Wipers - Small	280/Box	4.95
14-05514	Precision Wipers - Large	140/Box	7.95
30-34120	Delicate Task Wipers - Small	280/Box	6.95
30-34256	Delicate Task Wipers - Large	140/Box	12.95

ULTIDENTBRAND DELICATE TASK WIPERS

**NEW
PRODUCT**

- Available by the box, and as a case of 60 boxes
- Manufactured from low-lint, 100% virgin cellulose fiber
- Easily wipe up liquid and dust - Ideal for most delicate tasks



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
170-5180	UltidentBrand Wipers, 10.5 x 19.5 cm	280/Box	2.95

PARAFILM® M LABORATORY SEALING FILM & DISPENSERS



CAT. #	DESCRIPTION	PACKAGING	PRICE
36-PM996	Roll - 4" x 125'	Each	39.95
36-PM999	Roll - 4" x 250'	Each	69.95
33-120667	Parafilm® Top Box Dispenser, Blue	Each	43.95
33-120668	Parafilm® Top Box Dispenser, Purple	Each	43.95
33-234525-B	Parafilm® M ABS Dispenser	Each	72.95

ULTIDENTBRAND ALUMINUM FOIL ROLLS & PLASTIC WRAP



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
99-45100	Aluminum Foil - 45cm x 152m (18" x 500')	Each	39.95
99-45200	Plastic Wrap - 27.9cm x 762m (11" x 2500')	Each	49.95
33-23534A	Mini Aluminum Foil Roll, 10.2 x 1270 cm (4" x 500')	Each	25.95
33-23536A	Mini Aluminum Foil Roll, 15.2 x 1270 cm (6" x 500')	Each	31.95

ULTIDENTBRAND BENCH PROTECTORS

- Protects surfaces and keeps spills from spreading
- Polypropylene backing is waterproof and resistant to most chemicals



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
06-2050	20" x 50'	One Roll/Dispenser Box	79.95
06-20300	20" x 300'	Two Rolls/Cs	309.95

ULTIDENTBRAND BENCH PROTECTOR PADS

- Spunbond nonwoven topsheet with waterproof polypropylene backing



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
10-7105	17" x 24" Disposable Pads	300/CS	94.95

ULTIDENTBRAND AUTOCLAVABLE BIOHAZARD BAGS

**NEW
PRODUCT**



- Made from puncture resistant, high-molecular weight, high-density 4 mil polyethylene.
- Bright red color for quick identification
- Clearly labeled with bilingual biohazard warnings
- Wire stand holds 8 x 12 inch bags



**UltiDent
Scientific**

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-10320	8 x 12", Red	100/Cs	19.95
33-10321	14 x 19", Red	200/Cs	49.95
33-10322	19 x 23", Red	200/Cs	64.95
33-10323	25 x 33", Red	150/Cs	74.95
33-10044	Wire Stand	Each	31.95

ULTIDENTBRAND LABTAPE™

- Use with pencil, ball point pen or lab markers
- Fully autoclavable at +121°C - without melting or degrading



**UltiDent
Scientific**

CAT. #	DESCRIPTION	PACKAGING	PRICE
½" Wide x 60 Yards (13 mm x 55 m)			
53-PAT-13_*	Choose Color	Each	16.95
53-PAT-13AS	Rainbow Pack	10/Pk	149.95
¾" Wide x 60 Yards (18 mm x 55 m)			
53-PAT-18_*	Choose Color	Each	18.95
53-PAT-18AS	Rainbow Pack	10/Pk	159.95
1" Wide x 60 Yards (24 mm x 55 m)			
53-PAT-24_*	Choose Color	Each	19.95
53-PAT-24AS	Rainbow Pack	10/Pk	169.95



BENCHBIN™ BENCHTOP BIOHAZARD BIN

- Accepts 8 x 12 in. bags
- Easy push open/close spring lid
- Safely and conveniently dispose of tips and tubes



CAT. #	DESCRIPTION	PACKAGING	PRICE
48-A8401B	BenchBin™	Each	48.95

PROGENE® 15 & 50 ML POLYPROPYLENE RACKS

- Autoclavable, reusable and recyclable
- Widely spaced design allows tubes to be retrieved more easily



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
229419	25-Place rack for 15 ml tubes	Each	15.95
229429	25-Place rack for 50 ml tubes	Each	15.95

ULTIDENTBRAND TAPE DISPENSERS

- Dispensers: Heavy duty metal frame



**UltiDent
Scientific**

CAT. #	DESCRIPTION	PACKAGING	PRICE
53-TDSP-1	1 Inch Wide	Each	64.95
53-TDSP-2	2 Inches Wide	Each	69.95
53-TDSP-3	3 Inches Wide	Each	74.95
53-TDSP-6	5 Inches Wide	Each	105.95

ULTIDENTBRAND SEROLOGICAL PIPET RACK

- Four angled shelves protects and organizes glass and plastic pipettes
- Comes with two removable attachments that hold pipette bulbs, fillers, and pumps
- Magnetized rack mounts securely to most metal surfaces



**UltiDent
Scientific**

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-20612	Pipet Rack, Blue	Each	99.95
33-20613	Pipet Rack, Red	Each	99.95
33-20612M	Pipet Rack, Blue - Magnetized	Each	107.95

LABORATORY ESSENTIALS AND RACKS

ULTIDENTBRAND LAB MARKERS

- Alcohol-resistant and water-resistant

UltiDent
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
53-SM-1B	Science-Marker™, Black	Each	4.95
53-SM-1G	Science-Marker™, Green	Each	4.95

ULTIDENTBRAND LABORATORY NOTEBOOKS

- Archival-quality pages are numbered, acid-free - Pages include spaces for title, date, project number, etc.

UltiDent
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
33-8610-A	100 Page Notebook - Grid Style Black	Each	37.95
33-8610-B	200 Page Notebook - Grid Style Blue	Each	46.95
33-8610-D	200 Page Notebook - Lined Style Burgundy	Each	45.95

BLUELINE® LAB NOTEBOOK

- 10 1/2 x 8 inches - 200 pages with blue lines and red margin
- Made with 50-Percent of post-consumer recycled paper

UltiDent
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
9055-A91	200 Page Notebook - Lined Blue	Each	13.95

ULTIDENTBRAND 4-PLACE PIPETTOR WORKSTATION

- Store pipettes and tips together in one convenient stand
- Includes dual bins include hinged lids

UltiDent
Scientific



CAT. #	DESCRIPTION	PACKAGING	PRICE
33-234631	Acrylic Pipettor Workstation	Each	183.95

ULTIDENTBRAND 500 ML WASH BOTTLES



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Plain			
33-120245	White Cap	6/Pk	41.95
33-120246	Red Cap	6/Pk	41.95
33-120247	Yellow Cap	6/Pk	41.95
33-120248	Green Cap	6/Pk	41.95
33-120249	Blue Cap	6/Pk	41.95
Labeled			
33-120250	Red Cap - Acetone	6/Pk	66.95
33-120251	Clear Cap - Distilled Water	6/Pk	66.95
33-120252	Yellow Cap - Isopropanol	6/Pk	66.95
33-120253	White Cap - Ethanol	6/Pk	66.95
33-120254	Green Cap - Methanol	6/Pk	66.95
33-120255	Blue Cap - Bleach	6/Pk	66.95
33-120256	Label Your Own	6/Pk	66.95
33-120257	Assorted Pack (1 of each)	6/Pk	66.95

ULTIDENTBRAND STYROFOAM FREEZER BOXES

- Designed for the long-term storage of 1.5 or 2.0 ml microtubes
- 50-Place racks measure 210 x 110 x 71 mm -
- 100-Place racks measure 336 x 96 x 73 mm



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-4283	50-Place with Lid White	10/Pk	69.95
33-4282	100-Place with Lid White	10/Pk	80.95

BUY ANY 3 **VISTALAB™** PRODUCTS, GET 1 FREE (EQUAL OR LOWER PRICED)

VISTALAB™ OVATION® MECHANICAL (M) SINGLE-CHANNEL PIPETTES

- Mechanical 2-stroke low-force plunger operation with built-in microtube cap opener
- Mechanical precise 'click stop' volume adjustment
- Autoclavable for easy sterilization - Shelf hanger included with each pipette



CAT. #	DESCRIPTION	COLOUR	PACKAGING	PRICE
Mechanical (M) Single-Channel				
VL-1070-0002	0.2-2.0 µl	Red	Each	443.95
VL-1070-0010	1.0-10 µl	Orange	Each	443.95
VL-1070-0020	2.0-20 µl	Yellow	Each	443.95
VL-1070-0100	10-100 µl	Green	Each	443.95
VL-1070-0200	20-200 µl	Blue	Each	443.95
VL-1070-1000	100-1000 µl	Purple	Each	443.95
Mechanical (M) Single-Channel Starter Sets				
VL-1070-0998	10 µl, 100 µl and 1000 µl		Set of 3	1,165.95
VL-1070-0999	2 µl, 20 µl, 200 µl and 1000 µl		Set of 4	1,550.95

VISTALAB™ OVATION® ELECTRONIC PIPETTES



Quick Set Electronic Pipettes

- Ergonomically correct pipetting - Spring-loaded tip ejection for low thumb force
- Superior accuracy and precision - Easy in-lab calibration adjustment with no tools
- Unmatched patented ergonomic design - Tips load onto pipettors with an audible "click"
- QS pipettes allow you to program up to 5 frequently used aspirating and dispensing volumes at the press of a single button - Other volumes can also be selected throughout the pipettor's volume range

ESC, E8 & E12 Electronic Pipettes

- Individually serialized - Calibration certificate included
- Ergonomically correct pipetting - Spring-loaded tip ejection for low thumb force
- Superior accuracy and precision - Easy in-lab calibration adjustment with no tools
- Unmatched patented ergonomic design - Tips load onto pipettors with an audible "click"
- ESC, E8 and E12 pipettes allow you to program up to 3 frequently used aspirating and dispensing volumes at the press of a single button - Other volumes can also be selected throughout the pipettor's volume range
- 5 & 10 ml macro models offer up to 2 presets and can accept VistaLab macro tips and serological pipets



CAT. #	DESCRIPTION	COLOUR	PACKAGING	PRICE
Quick Set Electronic (QS)				
VL-1057-0002	0.2-2 µl	Red	Each	730.95
VL-1057-0010	1-10 µl	Orange	Each	730.95
VL-1057-0020	2-20 µl	Yellow	Each	730.95
VL-1057-0100	10-100 µl	Green	Each	730.95
VL-1057-0200	20-200 µl	Blue	Each	730.95
VL-1057-1000	100-1000 µl	Purple	Each	730.95
Quick Set Electronic (QS) Starter Set				
VL-1057-0999	20 µl, 200 µl and 1000 µl		Set of 3	2,045.95
Single-Channel Electronic (ESC)				
VL-1065-0020	0.5-20 µl	Yellow	Each	1,024.95
VL-1065-0125	2-125 µl	Green	Each	1,024.95
VL-1065-0250	5-250 µl	Blue	Each	1,024.95
VL-1065-1250	25-1250 µl	Purple	Each	1,024.95
VL-1065-05ML	100-5000 µl	Grey	Each	1,135.95
VL-1065-10ML	0.2-10 ml	Red	Each	1,162.95
Multi-8 Electronic (E8)				
VL-1160-0020	0.5-20 µl	Yellow	Each	2,150.95
VL-1160-0125	2.0-125 µl	Green	Each	2,150.95
VL-1160-0250	5.0-250 µl	Blue	Each	2,150.95
VL-1160-1250	25-1250 µl	Purple	Each	2,150.95
Multi-12 Electronic (E12)				
VL-1060-0020	0.5-20 µl	Yellow	Each	2,336.95
VL-1060-0125	2.0-125 µl	Green	Each	2,336.95
VL-1060-0250	5.0-250 µl	Blue	Each	2,336.95
VL-1060-0850	25-1250 µl	Purple	Each	2,336.95

Find the right tip
for your Ovation®
Pipette

PAGE 45

Pipette volume					
Ovation ESC, E8 & E12	0.5 - 20 µL	2 - 125 µL 5 - 250 µL	25 - 850 µL 25 - 1250 µL	n/a	n/a
Ovation QS & M	0.2 - 2 µL 1 - 10 µL 2 - 20 µL	10 - 100 µL 20 - 200 µL	100 - 1000 µL	n/a	n/a
Use this Tip	25 µL	250 µL	1250 µL	5 mL Tips	10 mL Tips

BUY ANY 3 **VISTALAB™** PRODUCTS, GET 1 FREE (EQUAL OR LOWER PRICED)

VISTALAB™ PIPETTE TIPS

- RNase/DNase free, Non-pyrogenic and trace metal certified
- **VistaBulk™**: Choose from Econo-Pak™ Bag or Protectainer™ packaging
- **VistaStak™**: 16x12 multi-channel layout - 192 tips per layer
- **Stacked Rack**: Racks of 200 tips per layer - Non-sterile options only
- Use genuine VistaLab™ tips with Ovation pipettes for perfect fit and reliable performance, every time



VistaLab™
TECHNOLOGIES
A subsidiary of CELLTREAT™ Scientific Products

CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
VistaBulk™			
VL-4058-1000	25 µl, Non-sterile	1000/Bag	65.95
VL-4058-2000	250 µl, Non-sterile	1000/Bag	63.95
VL-4058-3000	1250 µl, Non-sterile	1000/Bag	81.95
VL-4058-5000	5 ml, Non-sterile	250/Bag	71.95
VL-4058-6000	10 ml, Non-sterile	250/Bag	82.95
Protectainer™			
VL-4025	250 µl, Non-sterile	1000/Container	82.95
VL-4026	1250 µl, Non-sterile	1000/Container	86.95
VistaStak™			
VL-4060-9024	25 µl, Non-sterile	192 Tips x 3 Layers = 576	60.95
VL-4060-9025	250 µl, Non-sterile	192 Tips x 5 Layers = 960	93.95
VL-4060-9026	1250 µl, Non-sterile	192 Tips x 3 Layers = 576	69.95
VistaRak™			
VL-4060-1002	25 µl, Non-sterile	192 Tips x 5 Layers = 960	89.95
VL-4060-2004	250 µl, Non-sterile	192 Tips x 5 Layers = 960	89.95
VL-4060-3004	1250 µl, Non-sterile	192 Tips x 5 Layers = 960	107.95
VistaRak™ Filter Tips Sterile			
VL-4070-1002LR	25 µl, Low Retention, Non-sterile	192 Tips x 5 Layers = 960	106.95
VL-4070-2004LR	250 µl, Low Retention, Non-sterile	192 Tips x 5 Layers = 960	106.95
VL-4070-3004LR	1250 µl, Low Retention, Non-sterile	192 Tips x 5 Layers = 960	191.95
VL-4060-1032	25 µl, Sterile	192 Tips x 5 Layers = 960	118.95
VL-4060-2132	250 µl, Sterile	192 Tips x 5 Layers = 960	120.95
VL-4060-3132	1250 µl, Sterile	192 Tips x 5 Layers = 960	138.95
VL-4070-1032LR	25 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	142.95
VL-4070-2132LR	250 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	142.95
VL-4070-3132LR	1250 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	209.95
VistaRak™ Filter Tips Sterile			
VL-4070-1333LR	2 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	278.95
VL-4070-1332LR	25 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	217.95
VL-4070-2332LR	250 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	213.95
VL-4070-3332LR	1250 µl, Low Retention, Sterile	192 Tips x 5 Layers = 960	260.95

VISTALAB™ ALI-Q™ 2 ALIQUOTING PIPETTE CONTROLLERS

- Fast and easy programmable aliquoting - No need to watch the lines
- Ideal for cell culture, dispensing media, filling plates, tubes, cryovials, etc.
- 2-in-1: Conventional pipetting plus repeat pipetting all in a single pipette controller
- Use any size and any brand of serological pipets - Aspirate and dispense as you're used to
- When you first receive your ali-Q™ 2, press the purple button to wake it out of shipping mode
- Sensitive touch-control buttons - Pressure sensitive buttons that smoothly control pipetting speed every time

- **ali-Q™ 2**: Fixed high speed equal to speed 7 on VS model
- **ali-Q™ 2 LS (Low speed)**: Fixed low speed equal to speed 4 on VS model
- **ali-Q™ 2 VS (Variable speed)**: Custom adjustable: 7 speeds; 1 is slowest, 7 is fastest

	ali-Q™ 2	ali-Q™ 2 LS	ali-Q™ 2 VS
Aliquoting Volume Range	0.5 - 5.0 mL	0.3 - 3.0 mL	0.5 - 5.0 mL
Aliquot Speed Control Dial	No	No	Yes
On-board Calibration Software	Yes	Yes	Yes
Removable Kickstand	Yes	Yes	Yes

VistaLab™
TECHNOLOGIES
A subsidiary of CELLTREAT™ Scientific Products

CAT. #	DESCRIPTION	RANGE	PACKAGING	PRICE
VL-2100-0005	ali-Q™ 2 - High Speed	0.5 - 5.0 ml	Each	885.95
VL-2100-0003	ali-Q™ 2 LS - Low Speed	0.3 - 3.0 ml	Each	885.95
VL-2100-1005	ali-Q™ 2 VS - Variable Speed	0.5 - 5.0 ml	Each	1,039.95
VL-9020-4003	Replacement PTFE 0.22 µm Filter, Non-Sterile		5/Pk	67.95
VL-9020-4004	Replacement PTFE 0.22 µm Filter, Sterile		5/Pk	89.95
VL-9020-4005	Replacement PTFE 0.45 µm Filter, Non-Sterile		5/Pk	67.95
VL-9020-4006	Replacement PTFE 0.45 µm Filter, Sterile		5/Pk	89.95



TRANSFERPETTE® S SINGLE-CHANNEL PIPETTES

**BUY 2,
GET 1
FREE**

valid to June
30th 2023

BRANDTECH.

CAT. #	VOLUME	PRICE
34-705868	0.1-1 µl	510.95
34-705869	0.1-2,5 µl	510.95
34-705870	0.5-10 µl	510.95
34-705872	2-20 µl	510.95
34-705873	5-50 µl	510.95
34-705874	10-100 µl	510.95
34-705878	20-200 µl	510.95
34-705880	100-1000 µl	510.95
34-705882	500-5000 µl	510.95
34-705884	1000-10000 µl	510.95



TRANSFERPETTE® S MULTI-CHANNEL PIPETTES

**BUY 2,
GET 1
FREE**

equal or lower
price

BRANDTECH.

CAT. #	VOLUME	PRICE
34-705900	Multi-8, 0.5-10 µl	1,292.95
34-705906	Multi-8, 5-50 µl	1,292.95
34-705908	Multi-8, 10-100 µl	1,292.95
34-705910	Multi-8, 20-200 µl	1,292.95
34-705912	Multi-8, 30-300 µl	1,292.95
34-705920	Multi-12, 0.5-10 µl	1,598.95
34-705926	Multi-12, 5-50 µl	1,598.95
34-705928	Multi-12, 10-100 µl	1,598.95
34-705930	Multi-12, 20-200 µl	1,598.95
34-705932	Multi-12, 30-300 µl	1,598.95



PROMO

BUY 3 **ACCU-JET® S** PIPETTE CONTROLLERS,
GET YOUR CHOICE OF AN - **ACCU-JET® S**, OR
- **TRANSFERPETTE® S** SINGLE-CHANNEL
FOR FREE.

ACCU-JET® S PIPETTE CONTROLLER

- **Simple:** Enhanced comfort and precise control for delicate cell culture work
- **Easy:** Adjustable power and precision at your fingertips
- **Efficient:** Keep pipettes clean & sterile with integrated rest position



BRANDTECH.

CAT. #	COLOR	PRICE
34-26350	Anthracite	749.95
34-26351	Berry	749.95
34-26352	Petrol	749.95
34-26353	Amethyst	749.95

HALO™ SINGLE CHANNEL PIPETTOR

**BUY 3,
GET 1
FREE**

equal or lower
price

Easy in-lab
calibration - Key
supplied with
every pipette

Features soft
ejection and
universal tip
compatibility



The four digit
micrometer and
click-based volume
adjustment wheel
system prevents
volume drift

Ultra-lightweight
series of pipettors -
200 µl model weighs
only 68 gm

CAT. #	VOLUME	ACCURACY	REPRODUCIBILITY	PRICE
48-H6800-2	0.2-2 µl	+12.0-2.5%	<6.00-2.00%	279.95
48-H6800-10	1-10 µl	+2.5-1.0%	<1.50-0.80%	279.95
48-H6800-20	2-20 µl	+3.0-0.9%	<2.00-0.40%	279.95
48-H6800-50	5-50 µl	+2.0-0.6%	<2.00-0.30%	279.95
48-H6800-100	10-100 µl	+3.0-0.8%	<1.50-0.15%	279.95
48-H6800-200	20-200 µl	+3.0-0.6%	<1.00-0.15%	279.95
48-H6800-1M	100-1000 µl	+3.0-0.6%	<0.60-0.20%	279.95
48-H6800-1MXL	200-1250 µl	+3.0-0.6%	<0.60-0.20%	279.95
48-H6800-5M	0.5-5 ml	+2.4-0.6%	<0.60-0.15%	279.95
48-H6800-10M	1-10 ml	+3.0-0.6%	<0.60-0.20%	279.95
48-H6800-CR	Carousel Stand for up to 6 Pipettors			109.95

FREE HANDYSTEP® S MANUAL WITH PURCHASE OF 500 PD-TIP II

HANDYSTEP® REPEATING PIPETTES & SYRINGE TIPS

- **HandyStep® S manual:** Can pipette up to 90 pre-set volumes ranging from 2 µl to 50 ml
- **HandyStep® touch models:** Both models feature touchscreen operation, automatic tip size detection, automatic tip ejection, multi-language software and the ability to program any volume from 1 µl to 50 ml and include a universal USB-C AC adapter and shelf/rack mount
- **HandyStep® touch:** Functions include multi-dispensing, auto-dispensing, pipetting and favorites
- **HandyStep® touch S:** Multi-dispensing, auto-dispensing pipetting, sequential dispensing, multi-aspiration, titration and favorites



CAT. #	DESCRIPTION	PACKAGING	PRICE
Repeating Pipettes			
34-2705110	HandyStep® S Manual	Each	678.95
34-705200	HandyStep® touch, White	Each	1,769.95
34-705210	HandyStep® touch S, Black	Each	2,196.95
34-705230	HandyStep® touch Stand	Each	250.95

CAT. #	DESCRIPTION	PACKAGING	PRICE
PD-Tip™ II Non-sterile			
34-705700	0.1 ml	100/Pk	197.95
34-705702	0.5 ml	100/Pk	197.95
34-705704	1.0 ml	100/Pk	197.95
34-705706	1.25 ml	100/Pk	197.95
34-705708	2.5 ml	100/Pk	197.95
34-705710	5 ml	100/Pk	197.95
34-705712	10 ml	100/Pk	197.95
34-705714	12.5 ml	100/Pk	197.95
34-705716	25 ml (adapter included)	50/Pk	131.95
34-705718	50 ml (adapter included)	25/Pk	66.95



CAT. #	DESCRIPTION	PACKAGING	PRICE
PD-Tip™ II Sterile			
34-705730	0.1ml	1/Pk, 100/Pk	269.95
34-705732	0.5ml	1/Pk, 100/Pk	269.95
34-705734	1.0ml	1/Pk, 100/Pk	269.95
34-705736	1.25ml	1/Pk, 100/Pk	269.95
34-705738	2.5ml	1/Pk, 100/Pk	269.95
34-705740	5ml	1/Pk, 100/Pk	269.95
34-705742	10ml	1/Pk, 100/Pk	269.95
34-705744	12.5ml	1/Pk, 100/Pk	269.95
34-705746	25ml	1/Pk, 25/Pk	82.95
34-705748	50ml	1/Pk, 25/Pk	82.95

ULTIDENTBRAND REPEATING SYRINGE TIPS

- Precision manufacturing under strict quality controls
- 25- & 50-ml sizes include an adapter - Sterile formats of all sizes are available
- **Compatible with:** Biohit® Mechanical Stepper, BrandTech™ HandyStep™ S, Capp, eLine®, Eppendorf® Repeater 4780, Multipette® plus and Multipette stream®/Xstream®, Gilson® Distriman®, Nichiryo NichiMate and Thermo® Finnpipe Stepper

CERTIFIED
RNase, DNase & Pyrogen-Free

UltiDent
Scientific


CAT. #	DESCRIPTION	PACKAGING	PRICE
11-3920	0.05 ml Tips	100/Pk	169.95
11-3921	0.1 ml Tips	100/Pk	169.95
11-3922	0.2 ml Tips	100/Pk	169.95
11-3923	0.5 ml Tips	100/Pk	169.95
11-3924	1.0 ml Tips	100/Pk	169.95
11-3925	1.25 ml Tips	100/Pk	169.95
11-3926	2.5 ml Tips	100/Pk	169.95
11-3927	5 ml Tips	100/Pk	169.95
11-3928	10 ml Tips	100/Pk	169.95
11-3929	12.5 ml Tips	100/Pk	169.95
11-3930	25 ml Tips	25/Pk	64.95
11-3931	50 ml Tips	25/Pk	74.95

Sterile, Individually wrapped formats available

Visit us at ULTIDENT.COM

ULTIDENTBRAND PIPETTE CONTROLLER

- Long-life Li-ion battery with quick charge - Up to 20 hours continuous use
- UV resistant housing - Bright LEDs indicate three operation modes: high, low, and gravity
- **Includes:** Charging stand and wall-mountable/magnetic holder, Nozzle adapter for use with small pipettes, Li-Ion rechargeable battery with AC adapter, Three replacement 0.45 µm filters and Five nosecones



UltiDent
Scientific

CAT. #	DESCRIPTION	PRICE
33-3000	Pipette Controller	299.95

MC-24™ TOUCH HIGH SPEED MICROCENTRIFUGE

- Unique air-flow system, keeps samples cool
- Full Color Touch Screen with mySpin™ 1 button program recall
- **Speed:** 200 - 13,500 rpm (100 - 16,800 x g)
- **Capacity:** 24 x 1.5/2.0 ml tubes, 2 strips of 0.2 ml tubes (16 x 0.2ml)
- **Timer:** 20 sec. to 99 minutes, with momentary or continuous function



Benchmark

CAT. #	DESCRIPTION	PRICE
48-C2417	MC-24™ Touch High Speed Microcentrifuge with Combi-Rotor	3,384.95

MYFUGE™ 12 MINI-CENTRIFUGE

- No need to change rotors - Conveniently handles 12 microtubes and 4 PCR strips simultaneously
- **Speed:** Fixed at 2,000 x g - Simply close the lid and the unit quickly ramps up to 5500 rpm
- **Capacity:** 12 x 1.5/2.0ml tubes, 32 x 0.2ml tubes and/or 4 x PCR Strips
- Starts and stops with closing/opening of the lid



Benchmark

CAT. #	DESCRIPTION	PRICE
48-C1012	MyFuge™12 Mini-Centrifuge with Combi-Rotor	587.95

MYFUGE 5™ MICROCENTRIFUGE

- Includes adapters for 5 ml blood tubes
- Rotor holds 4 x 5 ml and 4 x 1.5/2.0 ml tubes
- Quickly accelerates to the fixed speed of 5,500 rpm



Benchmark

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C1005	MyFuge 5™ Microcentrifuge	Each	699.95

PLATEFUGE™ MICROFUGE

- Quickly spins down droplets in PCR plates
- Unique swing-out rotor design prevents spillage
- Accepts all popular PCR plates and most microtiter plates



Benchmark

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C2000	PlateFuge™ Microcentrifuge	Each	1,147.95

GUSTO® HIGH-SPEED MINI CENTRIFUGE

- ideal for protocols calling for fast spins with a higher RCF requirement. Ideal for PCR®, microfilter cell separation, clinical & diagnostic applications, and HPLC protocols.



CAT. #	DESCRIPTION	PRICE
33-10050	Gusto® High Speed Mini Centrifuge	1,472.95

SPROUT® PLUS MINI-CENTRIFUGE

- **Package Includes:** 6 x 1.5/2.0 ml rotor, 16 x 0.2 ml PCR rotor, 0.2 ml adapters, 0.5 ml adapters and storage case



CAT. #	DESCRIPTION	PRICE
33-120610	Sprout®, 6,000 rpm - Max RCF: 2,000 Xg, Green	435.95
33-120611	Sprout®, 6,000 rpm - Max RCF: 2,000 Xg, Blue	435.95

MAGFUGE® CENTRIFUGE & MAGNETIC STIRRER

- High-speed Centrifuge and Magnetic Stirrer in ONE
- Digital display, one button programming and tool free rotor exchange
- **MagFuge® Centrifuge:** RCF: 16-9,783 x g (12,500 RPM), Capacity: 12 x 1.5/2.0 ml tubes or 6 x 5 ml tubes
- **MagFuge® Magnetic Stirrer:** Speed Range: 50-2500 rpm, Capacity: 3L maximum
- **Package Includes:** 5 ml Rotor, 1.5/2.0 ml Rotor, Stirrer Rotor; 12 x 0.2 ml adapters, 12 x 0.5 ml adapters, 2 Silicone mats for stirrer lid & 2 Rare-Earth PTFE coated stir bars



CAT. #	DESCRIPTION	PRICE
33-120581	MagFuge® Centrifuge & Magnetic Stirrer	1,482.95

PROGENE® MINI CENTRIFUGE

- Ideal for quick spin downs - Simply close the lid and the rotor instantly reaches 6,000 rpm
- **Includes:** Large capacity 8-place microtube rotor, 16-place x 0.2 ml PCR tube and strip



CAT. #	DESCRIPTION	PRICE
48-1008-C	Mini Centrifuge with Clear Lid	299.95
48-1008-B	Mini Centrifuge with Blue Lid	299.95
48-1008-G	Mini Centrifuge with Green Lid	299.95
48-1008-P	Mini Centrifuge with Purple Lid	299.95
48-1008-R	Mini Centrifuge with Red Lid	299.95

MINI MAGNETIC STIRRER



- Capacity: Up to 3L
- Speed: 200 to 2200 rpm (Infinitely adjustable)
- Quiet operation with low vibration - Small but powerful



CAT. #	DESCRIPTION	PRICE
33-120155	Mini Magnetic Stirrer Grey/Purple	166.95



MINIMAG™ MAGNETIC STIRRER

- The 3.4 inch square work surface accepts a wide variety of beakers, flasks, bottles and other laboratory vessels
- Square footprint of under 5 inches
- Adjustable speed control, 100 to 1500rpm
- Effectively mixes volumes up to 1.5L
- The stirring speed is continuously adjustable from 100 to 1500rpm.
- Incubator and cold-room safe for use from +4°C to +45°C.



Benchmark

CAT. #	DESCRIPTION	PRICE
48-S1005	Mini Mag™ Magnetic Stirrer	159.95

MAGNETIC INDUCTION STIRRER

- This stirrer has no moving parts - 100% wear-free and maintenance free
- Magnetic coil technology works on the inductive principle with AC as its driving force
- **Multiple stirring modes:** Clockwise, counterclockwise and oscillate - Adjustable time: 30 seconds to 60 minutes or continuous - Speed: 10 to 2,000 rpm - Capacity: 3L total liquid volume, diameter of beaker not to exceed 165 mm
- **Package includes:** 2 Stir bars 35 mm x 12 mm, silicone anti-skid plate cover and low voltage, double insulated power adapter



Heathrow Scientific
design | function | innovation

CAT. #	DESCRIPTION	PRICE
33-120584	Grey with Purple Accents	538.95
33-120585	Grey with Blue Accents	538.95

ANALOG HOTPLATE STIRRERS

- Durable, chemical resistant white ceramic work surface, 7.1 x 7.1 in.
- Space-efficient design (8 X 9 in. footprint)
- Advanced microprocessor controls with convenient turn knobs allow quick, precise adjustment and maintenance of speed and temperature.
- Stir speed from 60 to 1500 rpm
- Temperature control up to 380°C (with safety indicator LEDs)
- Three models: heat-stir, heat only or stir only
- The hotplate/magnetic stirrer model (48-H4000-HS) comes complete with a support rod for mounting thermometers and temperature probes.

Benchmark

CAT. #	DESCRIPTION	PRICE
48-H4000-H	Hotplate, 7.5"x7.5"	440.95
48-H4000-S	Magnetic Stirrer, 7.5"x7.5"	440.95
48-H4000-HS	Hotplate Magnetic Stirrer, 7.5"x7.5"	665.95



Hotplate



Magnetic Stirrer



Hotplate Magnetic Stirrer

DIGITAL HOTPLATE STIRRERS

- Digital speed and temperature control
- Safety LED indicates hot surface (when top plate exceeds 50°C)
- Control actual sample temperature with optional probe attachment

Benchmark

CAT. #	DESCRIPTION	PRICE
48-H3770-H	Digital Hotplate, 7 x 7 inch	432.95
48-H3710-H	Digital Hotplate, 10 x 10 inch	752.95
48-H3770-S	Magnetic Stirrer, 7 x 7 inch	432.95
48-H3710-S	Magnetic Stirrer, 10 x 10 inch	752.95
48-H3770-HS	Digital Hotplate Magnetic Stirrer, 7 x 7 inch	657.95
48-H3710-HS	Digital Hotplate Magnetic Stirrer, 10 x 10 inch	1047.95



Digital Hotplate



Digital Magnetic Stirrer



Digital Hotplate Magnetic Stirrer

- **7 X 7" Digital Hotplate Stirrers:** feature 7x 7-inch white ceramic top plat with 5 L Maximum stirring capacity
- **10 X 10" Digital Hotplate Stirrers:** feature 10x 10-inch white ceramic top plat with 15 L Maximum stirring capacity

DISPENSETTE® S BOTTLETOP DISPENSERS

- Fast priming - Ideal for acids, bases and saline solutions
- Recirculation valve system enhances safety
- These digital Dispensette® S models feature Easy Calibration™
- Models for organic solvents and trace analysis also available



BRANDTECH.

CAT. #	DESCRIPTION	INCREMENTS (ml)	PRICE
34-4600311	0.1-1	0.005	980.95
34-4600321	0.2-2	0.01	980.95
34-4600331	0.5-5	0.02	980.95
34-4600341	1.10	0.05	980.95
34-4600351	2.5-25	0.1	1,306.95
34-4600361	5-50	0.2	1,322.95

**BUY 3,
GET 1
FREE**

equal or lower
price

BENCHMIXER™ & MORTEXER™ VORTEX MIXERS

- **BenchMixer™**: Powerful motor for INSTANT vortexing (115V) with continuous or "touch" operation - Variable speed from 200 to 3200 rpm - Balanced with Q-Drive technology, providing a smooth vortexing action, while minimizing noise and excessive vibration.
- **BenchMixer V2™**: Includes all the features of the original BenchMixer™ plus reversible flip top rubber head that can be attached to the mixer in two orientations.
- **Mortexer™**: Includes all the features of the original BenchMixer™ plus the Multi-Head™ which features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes.

Benchmark

CAT. #	DESCRIPTION	PRICE
48-BV1000	BenchMixer™ Vortex Mixer	369.95
48-BV1003	BenchMixer V2™ Vortex Mixer with flip top cup head	416.95
48-BV1005	Mortexer™ Vortex mixer with Multi-Head™	423.95



BenchMixer™



Mortexer™



BenchMixer V2™

VORTEXER – VARIABLE SPEED VORTEXING UNIT

- Provides holding slots for six different sizes of tubes and a standard assay plate, while the center pad area can be used for quick hand-held mixing of a vast range of tubes and vessels.
- Variable speed rotary dial provides pre-set speed control (1000, 2000, 3000 rpm) during vortex mixing

Heathrow Scientific
design | function | innovation



CAT. #	DESCRIPTION	PRICE
33-120209	Vortexer® Variable Speed	394.95

MINI VORTEX

- Touch start operation enables mixing with one hand
- Non-slip feet ensures the unit stays in place while mixing
- Speed: 4,500 rpm
- Orbit: 6 mm
- Capacity: Up to 50 ml tube



CAT. #	DESCRIPTION	PRICE
33-120567	Mini Vortex Gray/Purple	172.95
33-120598	Mini Vortex Gray/Blue	172.95

ORBI-SHAKER™ & ORBI-SHAKER™ JR

- Digital speed control up to 300 rpm
- Perfect horizontally circular 19mm orbit for aeration & mixing
- Instantly exchange flask clamps with MAGIC Clamp™ accessories
- **Orbi-Shaker™**: Touch screen display with 13 x 12 inch non-slip rubber coated platform - High capacity, up to 9 x 500ml, 5 x 1L or 2 x 2L
- **Orbi-Shaker™ Jr**: 9.5 x 8 inch non-slip rubber coated platform - Extremely compact, perfect for low to medium capacity applications

Benchmark

CAT. #	DESCRIPTION	PRICE
48-BT3001	Orbi-Shaker™ with 13 x 12 inch rubber mat platform	2,625.95
48-BT302	Orbi-Shaker™ Jr with 9.5 x 8 inch rubber mat platform	1,631.95



Orbi-Shaker™



Orbi-Shaker™ Jr

MULTI-THERM™ SHAKER

- Choose from two models: Heat/Cool or Heat only
- Exchangeable blocks (sold separately) accommodate tubes from 0.2 to 50 ml
- Simultaneous reading of actual and selected speed, time, and temperature settings
- Temperature controlled vortex (0 to 100°C), useful for a variety of molecular biology applications
- **Temp. Range**: 48-H5000-H: Ambient +5 to 100°C / 48-H5000-HC: Ambient -20 to 100°C
- **Speed**: 200-1500 rpm (3 mm orbit)

Benchmark



CAT. #	DESCRIPTION	PRICE
48-H5000-H	MultiTherm™ Shaker with Heating Only	3,391.95
48-H5000-HC	MultiTherm™ Shaker with Heating and Cooling	4,455.95
48-H5000-02	Block, 96 x 0.2 ml, or 1 X 96W PCR Plate	578.95
48-H5000-05	Block, 54 x 0.5 ml	486.95
48-H5000-CMB	Block, Combination 15 x 0.5 and 20 x 1.5 ml	486.95
48-H5000-15	Block, 35 x 1.5 ml	486.95
48-H5000-20	Block, 35 x 2.0 ml	486.95
48-H5000-12	Block, 24 x 12 mm	685.95
48-H5000-150	Block, 12 x 15 ml	685.95
48-H5000-500	Block, 6 x 50 ml	685.95
48-H5000-5MT	Block for 15 x 5 ml centrifuge tubes	685.95
48-H5000-MP	Block for One Micro Plate	685.95
48-H5000-DWMP	Block for One Deep Well Micro Plate	685.95

DIGITAL ORBITAL SHAKER

- Maximum load of 4.5 kg with an overload protection cut out
- **Timer range**: 0 to 49 hours, 1 second increments or continuous
- **Platform dimensions**: 27.9 W x 27.9 D cm
- Very low start speed for sensitive sample agitation with high top end speed
- **Speed range**: 20 to 300 rpm 1 rpm increments - Orbit diameter: 19.0 mm (0.75")

Heathrow Scientific
design | function | innovation



CAT. #	DESCRIPTION	PRICE
33-120460	Digital Orbital Shaker	1,444.95

BENCHROCKER™ 2D VARIABLE SPEED ROCKER

- Large 14 x 12 in. work surface
- Cold room and incubator safe (+4°C to +50°C)
- **Speed**: Variable 2 to 30 rpm
- **Tilt Angle**: Variable 0° to +/- 30°
- **Load Capacity**: 8 lbs. (3.6 kg)

Benchmark



CAT. #	DESCRIPTION	PRICE
48-BR2000	BenchRock™ 2D Variable Speed Rocker	1,145.95
48-BR2000-STACK	Optional Stacking Platform with Flat mat	176.95

BENCHWAVER™ 3D ROCKER

- Combines the basic rocking motion of traditional platform rockers with the circular motion of an orbital shaker.
- Instantly adjustable tilt from 0 to 10°
- Mixing speed can be adjusted from 5 to 105 rpm - Load Capacity: 20 lbs. / 9 kg
- Large, user-friendly LCD control panel
- Large 13.5 x 13 inch platform with a non-slip, rubber mat
- Select from 1 minute up to 99 hours of continuous operation
- Can be used in refrigerated or incubated environments from 2 to 60 degrees Celsius



Benchmark

CAT. #	DESCRIPTION	PRICE
48-B3D5000	BenchWaver™ 3d Rocker, variable speed	1,921.95

BENCHROCKER™ 3D NUTATING SHAKER WITH VARIABLE SPEED

- Includes 12 x 12 in. platform with 2 mats. (flat and dimpled-autoclavable)
- Adjustable speed & tilt
- Cold room and incubator safe (+4°C to +50°C)
- Maintenance free, brushless motor



Benchmark

CAT. #	DESCRIPTION	PRICE
48-B3D2300	BenchRocker™ 3D Nutating Shaker, variable speed	1,399.95

BIOMIXER™ & MINIMIXER™ 3D NUTATING SHAKER, FIXED SPEED

- Includes dimpled mat; ideal for foam-free mixing in blood & centrifuge tubes
- Tri-directional motion for thorough, uniform mixing
- Optional stacking platforms available
- Cold room and incubator safe (+4°C to +50°C)
- Stacking platforms can be purchased as accessories, effectively doubling the usable surface area of each model.
- Maintenance free, brushless motors capable of loads up to 2kg (4 lbs)

Benchmark

CAT. #	DESCRIPTION	PRICE
48-B3D1020	MiniMixer™ 3D Nutating Shaker, fixed speed	597.95
48-B3D1320	BioMixer™ 3D Nutating Shaker, fixed speed	743.95

BioMixer™ 3D
features a 12"x12" platform,
dimpled mat, 115V



MiniMixer™ 3D
features a 10"x7.5" platform,
dimpled mat, 115V

BLOTBOY™ & MINI BLOTBOY™ FIXED LOW SPEED 3D ROCKERS

- Perfect 3-dimensional motion for all blot related applications.
- Tilt angle and speed have been optimally set for gentle, but thorough mixing in gel trays, boxes and other flat vessels.
- Cold room and incubator safe
- Supplied with autoclavable flat mat - Optional stacking platforms available
- Designed to conserve bench space
- Maintenance free brushless motors

Benchmark

CAT. #	DESCRIPTION	PRICE
48-B3D1008	Mini BlotBoy™ 3D Rocker, fixed speed	614.95
48-B3D1308	BlotBoy™ 3D Rocker, fixed speed	743.95

BlotBoy™
features a 12"x12" platform
with flat mat



Mini BlotBoy™
features a 10"x7.5" platform
with flat mat

BENCHBLOTTER™ FIXED SPEED PLATFORM ROCKER

- Fixed speed platform rocker with 12x12" flat mat platform, simply press the power switch to start rocking
- Optimized for blot washing and gel staining with gentle mixing, 12 rpm +/- 8° tilt
- Cold room and incubator safe
- Optional stacking platforms available



Benchmark

CAT. #	DESCRIPTION	PRICE
48-BR1000	BenchBlotter™ Platform Rocker, fixed speed	709.95

ORBI-BLOTTER™ VARIABLE SPEED

- Gentle orbital motion for gel applications (Variable speed control from 3 to 70 rpm)
- Horizontally circular 19mm orbit for aeration & mixing
- Large 14 x 12 inch work surface
- 9 lbs (4 kg) load capacity
- Incubator and cold room safe



Benchmark

CAT. #	DESCRIPTION	PRICE
48-BT30	Orbi-Blotter™, variable speed	1,299.95

ROTO-THERM™ PLUS INCUBATED ROTATOR

- Includes tube holders for tubes from 0.5 ml to 50 ml as well as a 96 x 0.2 ml tube/strip/plate holder
- Suitable for resuspension, protein binding, nucleic acid purification, histochemistry, hybridization etc.
- Motion can be adjusted from vertical (end over end), to horizontal (rolling or hybridization style) mixing
- **Modes:** Advanced shaking/mixing modes: "Mix" function: Rotating in combination w/ oscillations (increased mixing), "Pause" function: Rotating with intermittent "stops" during rotation and "Rock" function: Back and forth, gentle rocking motion (see-saw)
- **Temperature:** Ambient +5 to 60°C
- **Speed:** Variable 5 to 70 rpm



Benchmark

CAT. #	DESCRIPTION	PRICE
48-H2024	Roto-Therm™ Plus Incubated Rotator	1,249.95

ROTOMINI™ & ROTOMINI™ PLUS TUBE ROTATORS

- Suitable for resuspension, protein binding, nucleic acid purification, histochemistry, hybridization etc.
- Includes tube holders for common tube sizes from 0.5 ml to 50 ml
- Motion can be adjusted from vertical (end over end), to horizontal (rolling or hybridization style) mixing
- **RotoMini™:** Fixed Speed 24 rpm. Max. Capacity: 44 x 1.5/2.0ml - 20 x 15ml - 6 x 50ml
- **RotoMini™ Plus:** Includes all features of the RotoMini™ and adds: Variable Speed 5 to 70 rpm; Advanced shaking/mixing modes: Oscillating, Pause and Rock, as well as a 96 x 0.2 ml tube/strip/plate holder

Benchmark

CAT. #	DESCRIPTION	PRICE
48-R2020	RotoMini™ Tube Rotator	578.95
48-R2024	RotoMini™ Plus Tube Rotator	681.95



BEADBATH DUO™

- Designed for precision temperature control when used with water, or with thermal beads
- Temp. Range: Water Mode: Ambient +5 to 80°C / Bead Mode: Ambient +5 to 70°C
- "MODE" button can be used to toggle between "Bead" or "Water" modes
- This combination bath that is suitable for accurate, and uniform temperature control in either condition - Get the best of both worlds, allowing conversion from one bath type to the other, with the simple press of a button



CAT. #	DESCRIPTION	PRICE
48-B2402	BeadBath Duo™ 2L with Lift-Off Lid	1,159.95
48-B2410	BeadBath Duo™ 10L with Hinged Lid	1,973.95



MYBATH™ DIGITAL WATER BATH

- Superior accuracy and performance in a compact footprint for your benchtop
- Seamless, stainless steel chamber for easy cleaning
- See-through cover that can be removed, or flipped back on a rear hinge (2L model lid is not hinged)
- Quik-CAL™ user calibration function
- Large, electronic 3-digit display adjusts temperature in 0.1°C increments (2-digit display / 1°C increments for 2L model)
- The 4, 8 and 12L models allow the user to digitally select and monitor settings
- 8 and 12L models feature a drain for hassle-free water changes
- 4L Bath: Includes three separate tube racks for 40 x 1.5/2.0, 41 x 15 ml and 15 x 50 ml tubes

CAT. #	DESCRIPTION	PRICE
48-B2000-2	MyBath™ 2L Digital Water Bath	957.95
48-B2000-4	MyBath™ 4L Digital Water Bath	1,467.95
48-B2000-8	MyBath™ 8L Digital Water Bath	1,615.95
48-B2000-12	MyBath™ 12L Digital Water Bath	1,926.95



ASPIRE™ OIL-FREE LABORATORY ASPIRATOR

- Engineered to remove media from tissue culture vessels or plates, and supernatants
- Includes single-channel handheld controller that utilizes universal 200 or 300 µl pipette tips
- Controller and adapters are made of POM (Polyoxymethylene), and are fully autoclavable
- Polycarbonate collection bottle incorporates a reliable Overflow Protection Valve that closes off the tube to the vacuum in the event of a full bottle
- Package Includes: Aspire™ Laboratory Aspirator with base and internal vacuum pump, 2L polycarbonate bottle with lid, silicon tubing, handheld vacuum controller and single channel adapters
- Specifications: Max. Aspiration: 14 ml / sec. - Max. Flow: 15L per min. - Vacuum Range: -30 to -750 mbar - Noise Level 59.6 db



CAT. #	DESCRIPTION	PRICE
48-V0020	Aspire Laboratory Aspirator	1,518.95
48-V1002-1SL	Optional single channel adapter, XL (80 mm) stainless steel	69.95
48-V1002-8	Optional 8-Channel Adapter for use with 200 µl Pipette Tips	338.95
48-V1002-S8	Optional 8-Channel Adapter, Fine Tip Stainless Steel	205.95





For over 30 years, Ultident Scientific has been a trusted leader in the distribution of lab supplies and equipment across Canada. We pride ourselves in delivering exceptional value and customer service, along with an ever-growing portfolio of industry-leading brands.

Contact Us

Ultident.com
Sales@ultident.com
1-800-463-4163

Ultident Scientific
381, rue Locke
St. Laurent, QC H4T 1X7



*Prices are subject to change without notice