

	1	2	3	4	5	6	7	8	9	10
PCR	●	●	●	●		●		●	●	●
qPCR	●	●	●	●	●		●	●	●	●
APPLIED BIOSYSTEMS® (ABI)										
3100, 3130XL, 3500, 3500XL, 3730, 3730XL										●
7300, 7500	●		●	●						
7500 Fast					●					
7700	●		●	●		●				
7900 HT (with 0.1/FAST block)					●					
7900 HT (with 0.2 mL block)			●	●		●				
7900 HT (with 384-Well block)										●
GeneAmp 2700, 2720			●	●		●				
GeneAmp 9600, 9700			●	●		●				●
GeneAmp 9800					●					●
MiniAmp, MiniAmp Plus			●	●		●				
ProFlex			●	●		●				●
QuantStudio 3, 5, 6, 7, 12 K Real-Time PCR System (with 0.1 mL/FAST block)					●					
QuantStudio 3, 5, 6, 7, 12 K Real-Time PCR System (with 0.2 mL block)			●	●		●				
QuantStudio 3, 5, 6, 7, 12 K Real-Time PCR System (with 384-Well block)										●
SimpliAmp			●	●		●				
Veriti 96-Well (with 0.1 mL/FAST block)					●					
Veriti 96-Well (with 0.2 mL block)			●	●		●				
Veriti 384-Well										●
ViiA7, 7900 HT FAST					●					●
AGILENT®										
AriaMx qPCR System		●			●			●		
Mx3000P, Mx3005P, Mx4000	●					●				
Robocycler Gradient						●		●		
Surecyclcr 8800	●									●
ANALYTIK JENA®										
Trobot 96, SpeedCycler2 SP, SPR	●					●		●		●
Alpha SC	●					●		●		●
Flexcycler2	●					●		●		●
T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series	●					●		●		●
qTOWER3 (G, touch), Toptical	●				●	●		●		●
BIOER® / BULLDOG®										
GeneExplorer	●	●				●				
GeneTouch, GeneTouch Plus	●				●	●		●		●
LifeECO	●					●		●		
LifeTouch	●				●	●		●		
LifePro	●				●	●		●		
BIOMETRA®										
Biometra T Advanced	●	●			●	●		●		●
Biometra T One	●	●			●	●		●		
Biometra T Professional series	●	●			●	●		●		
Biometra Uno	●	●			●	●		●		
BIO-RAD® / MJ RESEARCH®										
C1000, C1000 Touch, S1000	●	●			●	●		●		●
CFX Connect, CFX96, CFX96 Touch, Deep Well		●			●			●		●
CFX384, CFX384 Touch										●
Chromo4		●			●			●		
DNA Engine Opticon and Opticon 2		●			●			●		
DNA Engine, DNA Engine Tetrad, Tetrad 2					●	●		●		●
DNA Engine Dyad, Dyad Disciple					●	●		●		●
iCycler	●					●				
iQ, iQ-4, iQ-5, myiQ PCR System	●					●				
MJ Research PTC200, 220, 221, 225, 240					●	●		●		●
MyCycler	●									
PTC 200					●	●		●		●
T100	●					●				
CORBETT® / QIAGEN®										
Corbett Research (Qiagen) Palm Cycler	●					●				
Corbett Research (Qiagen) Palm Cycler 384										●
EPPENDORF®										
Mastercycler ep realplex	●	●			●	●		●		●
MasterCycler ep, ep gradient, Pro, Pro S	●	●			●	●		●		●
MasterCycler ep 384 / Pro 384										●
Nexus, Nexus gradient, Nexus SX1, Nexus GSX1	●	●			●	●		●		
GE / MEGABACE®										
GE Healthcare Amersham MegaBACE 4000										●
GE Healthcare Amersham MegaBACE 500, Mega-BACE1000 mark 2								●		
ROCHE®										
LightCycler 480							●		●	
TECHNE® / COLE-PARMER®										
Flexigene	●	●				●		●		●
Techne Genius, Genius Quad	●	●				●		●		●
Techne PCRmax Alpha cycler 1, 2, 4	●	●			●	●		●		●
Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite	●	●			●	●		●		●
Techne Quantica		●						●		●
Techne TC412, TC512	●	●				●		●		●
Techne TC-PLUS	●	●			●	●		●		●
Touchgene, Touchgene Gradient	●	●				●		●		●
OTHER BRANDS / MODELS										
Apollo ATC-201										
BIONEER AllInOneCycler	●					●				●
DLab TC1000-G	●					●				
G-Storm GS1, GS2, GS4	●	●			●	●		●		●
HiMedia® Prima-96	●					●				
Labnet MultiGene, OptiMax	●		●	●		●				
Idaho Technology LightScanner								●		●
peqLab peqSTAR	●	●			●	●		●		●
Sensoquest Labcycler	●	●			●	●		●		●
Takara Dice touch, Gradient	●					●				
ThermoFisher Hybaid MBS, Px2/PxE	●	●			●	●		●		
VWR Collection UNO96	●	●			●	●		●		●
VWR peqlab peqSTAR	●	●			●	●		●		



Standard-Shape (Cutttable),
96-Well, Non-Skirted,
Volume: 0.2 mL, A12 Cut Corner

87-C96-NS

1



Low Profile
96-Well, Non-Skirted,
Volume: 0.1 mL, H12 Cut Corner

87-C96-LPNS

2



Standard ABI® Style (Flat Top)
96-Well, Semi-Skirted, Volume: 0.2 mL,
A12 Cut Corner, ABI/Hybrid

87-C96-ABI-1

3



ABI® Shape (Raised Top Ridge)
96-Well, Semi-Skirted,
Volume: 0.2 mL, A12 Cut Corner

87-C96-RR

4



ABI® Style FAST®, Low-Profile
96-Well, Semi-Skirted,
Volume: 0.1 mL, A1 Cut Corner

87-C96-FAST

5



Standard-Shape
96-Well, Semi-Skirted, Volume: 0.2 mL,
A12 Cut Corner, ABI/Universal

87-C96-ABI-2

6



LightCycler®, Low-Profile
96-Well, White, Semi-Skirted,
Volume: 0.1 mL, H12 Cut Corner

87-C96-LC480-WP

7



Low-Profile
96-Well, Full-Skirted,
Volume: 0.1 mL, A12 Cut Corner

87-C96-SK

8



Two-Notch Type
384-Well, Full-Skirted,
Volume: 30 µL, A24 Cut Corner

87-C384-2N

9



One-Notch Type
384-Well, Full-Skirted,
Volume: 30 µL, A24, P24 Cut Corner

87-C384-1N

10



***Disclaimer:** This chart serves as a general guide and is subject to change. There may be additional compatibilities not listed in this chart. For the most accurate and up-to-date information, please contact your representative or e-mail us at sales@ultident.com to request a sample.