

CHEERS TO 2022

HAPPY NEW YEAR!

Valid to March 31st, 2022

ULTIDENTBRAND UNIVERSAL-FIT LOW-RETENTION FILTER TIPS

- Manufactured from super clear virgin polypropylene
- Certified Non-Pyrogenic & DNase/RNase-free - Gamma sterile
- Free of Nucleic Acids, Pyrogens/Endotoxins, PCR Inhibitors and trace Metals
- Tips are super hydrophobic, significantly reducing sample loss and increasing precision
- Provides high reproducibility

**PROMO
PRICE**

IN STOCK



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
50-6010	10 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6010-XL	10 µl Reach, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6011	20 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6012	100 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6013	200 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6014	300 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	64.95
50-6015	1000 µl, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	79.95
50-6015-XL	1250 µl Reach, Graduated Filter Tips, Sterile	10 Hinged Racks of 96	89.95



AURELIA® TRANSFORM™ 200 POWDER-FREE NITRILE GLOVES

aurelia

- Optimal stretchability, just like latex!
- Long-lasting comfort, durability and decreased hand fatigue
- Beaded cuff snugs wrist, eliminates roll-up, minimizes entry of liquids
- Textured fingertips for precision gripping
- Latex-free, hypoallergenic
- Non-sterile, ambidextrous
- Chemotherapy-tested

OPTIMAL STRETCHABILITY



3.2 MIL

200/BOX
LIMITED QUANTITIES
AVAILABLE

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-98895	Transform™ 200 Stretch Nitrile, X-Small	Ice Blue	200/Box	25.95
262-98896	Transform™ 200 Stretch Nitrile, Small	Ice Blue	200/Box	25.95
262-98897	Transform™ 200 Stretch Nitrile, Medium	Ice Blue	200/Box	25.95
262-98898	Transform™ 200 Stretch Nitrile, Large	Ice Blue	200/Box	25.95
262-98899	Transform™ 200 Stretch Nitrile, X-Large	Ice Blue	200/Box	25.95

Prices Subject to change without notice.



ULTIDENT.COM



Sales@ultident.com



1-800-463-4163 / 514-335-3433

PIPETTE TIPS

NoStick® UNIVERSAL-FIT FILTER BARRIER TIPS

- Clean Liquid Release
- Vertex NoStick® tips offer users the highest degree of fluid control
- Vertex tips prove that less is more when it comes to liquid handling!



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

CAT. #	DESCRIPTION	PACKAGING	PRICE
4117NSFS	10 µl, Graduated	10 Hinged Racks of 96	79.95
4137NSF	10 µl Reach, Graduated	10 Hinged Racks of 96	79.95
4237NAFS	20 µl, Graduated	10 Hinged Racks of 96	79.95
4237NBFS	100 µl, Graduated	10 Hinged Racks of 96	79.95
4237NSFS	200 µl, Graduated	10 Hinged Racks of 96	79.95
4927NSFS	300 µl, Graduated	10 Hinged Racks of 96	79.95
4337NSFS	1000 µl, Graduated	10 Hinged Racks of 100	93.95
4347NSFS	1250 µl Reach, Graduated	10 Hinged Racks of 96	95.95

PROGENE® UNIVERSAL-FIT TIPS

- Tips are fully autoclavable - Choose Research-Grade or Low-Retention formats
- All products made with colorants are guaranteed to be heavy metals-free
- Manufactured from virgin, medical-grade FDA-approved, bovine-free plastic resins in medically-clean facility
- All molds are diamond-polished to a glass-like finish ensuring optimum performance
- Bulk Formats: Packaged in resealable Zip-Lock® bags
- **Low-Retention Formats:** Minimal sample retention and maximum sample recovery



CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
Research-Grade			
87-A10	10 µl, Graduated	1000/Pk	15.95
87-A10XL	10 µl Reach, Graduated	1000/Pk	15.95
87-A200C	200 µl, Graduated, Clear	1000/Pk	13.95
87-A200Y	200 µl, Graduated, Yellow	1000/Pk	13.95
87-A1000B	1000 µl, Graduated, Blue	1000/Pk	19.95
Low-Retention			
87-A10-LR	10 µl, Graduated	1000/Pk	22.95
87-A10XL-LR	10 µl Reach, Graduated	1000/Pk	22.95
87-A200-LR	200 µl, Graduated	1000/Pk	17.95
87-A300-LR	300 µl Reach, Graduated	1000/Pk	26.95
87-A1000C-LR	1000 µl, Graduated	1000/Pk	26.95
87-A1250	1250 µl, Reach, Graduated	1000/Pk	27.95



MICROTUBES & CRYOGENIC VIALS

PROGENE® PREMIUM BOIL-PROOF MICROTUBES

- 0.6, 1.5, 1.7 & 2.0 mL tubes are rated to 20,000 Xg
- 5 mL tubes are rated to 25,000 Xg
- Temperature range from -80°C to +121°C - Boiling or freezing, tubes will not pop
- Caps "click" when closing and are still easy to open
- Crystal-clear polymer facilitates visualization of samples.
- Graduations at every 100 µL for easy volume estimation.
- Manufactured from virgin, medical-grade polypropylene



CERTIFIED
RNase, DNase & Pyrogen-Free

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-B600-C	0.6 ml Boil-Proof	1000/Pk	27.95
87-B600-C-LR	0.6 ml Boil-Proof, Low-Retention	500/Pk	28.95
87-B600-CS	0.6 ml Boil-Proof, Sterile	500/Pk	26.95
87-B600-AS	0.6 ml Boil-Proof, Asst'd Colors	1000/Pk	34.95
87-B150-C	1.5 ml Boil-Proof	500/Pk	14.95
87-B150-C-LR	1.5 ml Boil-Proof, Low-Retention	250/Pk	15.95
87-B150-CS	1.5 ml Boil-Proof, Sterile	250/Pk	15.95
87-B150-AS	1.5 ml Boil-Proof, Asst'd Colors	500/Pk	19.95
87-B150-NC	1.5 ml Microtubes, No-Cap	500/Pk	19.95
87-B175-C	1.7 ml Boil-Proof, Large Cap	500/Pk	19.95
87-B200-C	2.0 ml Boil-Proof	500/Pk	18.95
87-B200-C-LR	2.0 ml Boil-Proof, Low-Retention	250/Pk	19.95
87-B200-CS	2.0 ml Boil-Proof, Sterile	250/Pk	15.95
87-B200-AS	2.0 ml Boil-Proof, Asst'd Colors	500/Pk	22.95
87-B200-NC	2.0 ml Microtubes, No-Cap	500/Pk	20.95
48-B500-C	5 ml Snap-Cap	200/Pk	51.95
48-B500-CS	5 ml Snap-Cap, Sterile	50/Bag, 200/Pk	55.95
48-B500-AMS	5 ml Snap-Cap, Amber, Sterile	50/Bag, 200/Pk	59.95
48-C2530	5 ml Screw-Cap	500/Cs	134.95
48-C2540	5 ml Screw-Cap, Sterile	50/Rack, 500/Cs	169.95

PROGENE® CRYOGENIC VIALS

**NEW
PRODUCT**

- Optimized for vapor-phase liquid nitrogen storage
- Temperature range -196°C to 121°C
- Vials feature a white writing surface
- Made from virgin polypropylene and caps
- Includes secure O-ring seals
- Packaged in heavy-duty, tamper-proof, resealable bags
- Certified Non-Pyrogenic & ATP,
- DNase/RNase-free - Gamma sterile



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
Internal-Thread			
71-ITCV-2SK	2.0 ml Internal-Thread, Skirted	50/Bag, 1000/Cs	379.95
71-ITCV-4SK	4.0 ml Internal-Thread, Skirted	50/Bag, 1000/Cs	379.95
External-Thread			
71-ETCV-2SK	2.0 ml External-Thread, Skirted	50/Bag, 1000/Cs	379.95
71-ETCV-5SK	5.0 ml External-Thread, Skirted	50/Bag, 1000/Cs	379.95

PROGENE® 5 & 10 ML TRANSPORT TUBES

- Sterile, RNase, DNase, Pyrogen free
- Meets 95kPa transportation guidelines
- Choose from 5 ml (16 x 60 mm) and 10 ml (16 x 80 mm) formats
- Polypropylene construction with white printed graduations and marking area



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
48-C1811	5 ml Transport Tubes	100/Bag, 1000/Cs	215.95
48-C1812	10 ml Transport Tubes	100/Bag, 1000/Cs	251.95



ULTIDENT.COM



Sales@ultident.com



1-800-463-4163 / 514-335-3433

3

PROGENE® CENTRIFUGE TUBES

- Choose from bulk, recyclable paper racks or reusable polypropylene racks
- Graduations molded into the conical bottom for small volume quantification
- Manufactured from medical-grade virgin polypropylene with flat polyethylene (HDPE) caps
- Leak-proof tubes feature a large unerasable frosting white printed writing area
- Easy-to-read black graduations are accurate to $\pm 2\%$ - 1 ml increments (15's) or 2.5 ml increments (50's)
- Tubes also feature molded graduations at the conical bottom of each tube
- Both the graduations and writing areas are chloroform-resistant
- Packaged in secure, tamper-proof, resealable bags - Can be used down to -80°C



CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE

CAT. #	DESCRIPTION	RFC	PACKAGING	PRICE
Conical Tubes				
71-1500-B	15 ml, Bulk	12,000 Xg	25/Bag, 500/Cs	94.95
71-1500-R	15 ml, Paper Racks	12,000 Xg	25/Rack, 500/Cs	109.95
71-1500-PP	15 ml, Polypropylene Racks	12,000 Xg	50/Rack, 300/Cs	129.95
71-5000-B	50 ml, Bulk	12,000 Xg	25/Bag, 500/Cs	129.95
71-5000-R	50 ml, Paper Racks	12,000 Xg	25/Rack, 500/Cs	149.95
71-5000-PP	50 ml, Polypropylene Racks	12,000 Xg	25/Rack, 300/Cs	159.95



PROGENE® SEROLOGICAL PIPETTES

- Paper/Plastic wrap with color-coding on paper allows for easy size identification
- Bi-directional graduations provide choice of ascending and descending scales
- Printed indicator arrow for quickly locating negative graduations when using a pipette controller
- Features a premium material barrier plug
- Made of polystyrene plastic
- Sterile, individually wrapped

CERTIFIED
RNase, DNase & Pyrogen-Free
STERILE

PROGENE



CAT. #	VOLUME	PACKAGING	PRICE
85-7125	1 ml	500/Cs	99.95
85-7126	2 ml	500/Cs	99.95
85-7127	5 ml	200/Cs	42.95
85-7128	10 ml	200/Cs	42.95
85-7129	25 ml	200/Cs	96.95
85-7130	50 ml	100/Cs	132.95
85-7131	100 ml	50/Cs	149.95



PROGENE® CELL CULTURE PLATES

- **1-Well Plate:** Same footprint as standard multiple well plates - Use for high volume and robotic screening of DNA libraries and manual plating - Also suitable for dot-blotting bacterial clones on 11.4 x 7.3cm (4.5 x 2.875 in.) nylon membranes
- **4-Well Plate:** Smaller footprint with well diameter similar to that of a 24 well plate with expanded well spacing



CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

Tissue Culture Treated

229105	6-Well with Lid, Flat Bottom	1/Pk, 50/Cs	129.95
229111	12-Well with Lid, Flat Bottom	1/Pk, 50/Cs	133.95
229123	24-Well with Lid, Flat Bottom	1/Pk, 50/Cs	133.95
229147	48-Well with Lid, Flat Bottom	1/Pk, 50/Cs	137.95
229195	96-Well with Lid, Flat Bottom	1/Pk, 50/Cs	149.95

Non-Treated Suspension Culture

229506	6-Well with Lid, Flat Bottom	1/Pk, 100/Cs	224.95
229512	12-Well with Lid, Flat Bottom	1/Pk, 100/Cs	230.95
229524	24-Well with Lid, Flat Bottom	1/Pk, 100/Cs	229.95
229548	48-Well with Lid, Flat Bottom	1/Pk, 100/Cs	249.95
229596	96-Well with Lid, Flat bottom	1/Pk, 100/Cs	254.95

PROGENE® CELL SCRAPERS & LIFTERS

PROGENE

CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE



CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

71-800-25	Cell Scraper	1/Pk, 100/Cs	164.95
229305	Cell Lifter, Flat & Narrow Blade	1/Pk, 100/Cs	139.95
229306	Cell Lifter, Flat & J-Hook Blade	1/Pk, 100/Cs	139.95

PROGENE® CELL CULTURE FLASKS



- Sterilized by gamma irradiation
- Manufactured in Class 100,000 clean-room and TC-Treated by vacuum gas plasma

CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

Tissue Culture Treated

229331	25 cm ² Vent Cap	10/Bag, 200/Cs	263.95
229341	75 cm ² Vent Cap	5/Bag, 100/Cs	227.95
229351	182 cm ² Vent Cap	5/Bag, 40/Cs	205.95
229361	300 cm ² Vent Cap	3/Bag, 18/Cs	217.95

Non-Treated Suspension Culture

229500	25 mL Vent Cap	10/Bag, 200/Cs	231.95
229530	600 mL Vent Cap	5/Bag, 40/Cs	189.95
229540	850 mL Vent Cap	3/Bag, 18/Cs	201.95

PROGENE® CELL STRAINERS



CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

229481	40 µm Cell Strainer, Blue	1/Pk, 50/Cs	133.95
229483	70 µm Cell Strainer, White	1/Pk, 50/Cs	133.95
229485	100 µm Cell Strainer, Yellow	1/Pk, 50/Cs	133.95
229480	Pestle, Green	1/Pk, 50/Cs	73.95

PROGENE® CELL CULTURE DISHES

- Sterilized by gamma irradiation
- Manufactured in Class 100,000 clean-room and TC-Treated by vacuum gas plasma



CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE

PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
--------	-------------	-----------	-------

229661	60mm x 15mm	20/Bag, 500/Cs	194.95
667660	60mm x 15mm with Grip Ring	10/Bag, 500/Cs	269.95
229621	100mm x 20mm	10/Bag, 300/Cs	194.95
229620	100mm x 20mm with Grip Ring	10/Bag, 300/Cs	249.95
667651	150mm x 20 mm	5/Bag, 100/Cs	206.95
229652	150mm x 25 mm	5/Bag, 60/Cs	171.95



PROGENE® VACUUM FILTRATION

- Low binding PES membrane for high flow rate
- 0.10 µm membrane for mycoplasma removal, 0.22 µm membrane for filtration of media, 0.45 µm PES membrane for solution clarification
- Supplied complete with separate O-ring cap for receiver bottle
- Easy-to-read molded graduations marks on bottle top filter with bold easy-to-read printed membrane and pore size for easy identification
- Angled hose barb adapter fits three size hose diameters and prevents the filter from tipping



CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE



CAT. #	DESCRIPTION	PACKAGING	PRICE
71-PES-150-22	150 ml, PES, 0.22 µm	1/Bag, 12/Cs	130.95
71-PES-250-22	250 ml, PES, 0.22 µm	1/Bag, 12/Cs	152.95
71-PES-500-22	500 ml, PES, 0.22 µm	1/Bag, 12/Cs	184.95
71-PES-1000-22	1000 ml, PES, 0.22 µm	1/Bag, 12/Cs	279.95
71-PES-150-45	150 ml, PES, 0.45 µm	1/Bag, 12/Cs	131.95
71-PES-250-45	250 ml, PES, 0.45 µm	1/Bag, 12/Cs	154.95
71-PES-500-45	500 ml, PES, 0.45 µm	1/Bag, 12/Cs	189.95
71-PES-1000-45	1000 ml, PES, 0.45 µm	1/Bag, 12/Cs	284.95

PROGENE® SYRINGE FILTERS

- **PES (Polyethersulfone)** - Low affinity for proteins and extractables - Faster flow rate than PVDF
- **CA (Cellulose Acetate)** - Low protein binding, reduced extractables and low hold-up volume
- **MCE (Mixed Cellulose Ester)** - Filtration of aqueous solutions - Effectively binds trace proteins
- **PVDF (Polyvinylidene fluoride)** - Low-binding non-aggressive aqueous and mild organic solutions
- **NYLON** - Universal application for analytical procedures - Chemical filtration
- **PTFE (Polytetrafluoroethylene)** - Hydrophobic - Strong acid solvent filtration



CERTIFIED

RNase, DNase & Pyrogen-Free

STERILE



CAT. #	DESCRIPTION	PACKAGING	PRICE
Sterile Blister Packs, Individual-wrapped			
229746	PES, 0.22 µm, 13mm	75/Pk	79.95
229747	PES, 0.22 µm, 30mm	30/Pk	74.95
229748	PES, 0.45 µm, 13mm	75/Pk	79.95
229749	PES, 0.45 µm, 30mm	30/Pk	74.95
229762	CA, 0.22 µm, 30mm	30/Pk	79.95
229763	CA, 0.45 µm, 30mm	30/Pk	79.95
229750	MCE, 0.22 µm, 13mm	75/Pk	79.95
229751	MCE, 0.22 µm, 30mm	30/Pk	74.95
229752	MCE, 0.45 µm, 13mm	75/Pk	79.95
229753	MCE, 0.45 µm, 30mm	30/Pk	64.95
229756	PTFE, 0.22 µm, 13mm	75/Pk	149.95
229757	PTFE, 0.22 µm, 30mm	30/Pk	106.95
229758	PTFE, 0.45 µm, 13mm	75/Pk	149.95
229759	PTFE, 0.45 µm, 30mm	30/Pk	106.95
229740	PVDF, 0.10 µm, 13mm	75/Pk	89.95
229741	PVDF, 0.10 µm, 30mm	30/Pk	74.95
229742	PVDF, 0.22 µm, 13mm	75/Pk	89.95
229743	PVDF, 0.22 µm, 30mm	30/Pk	74.95
229744	PVDF, 0.45 µm, 13mm	75/Pk	89.95
229745	PVDF, 0.45 µm, 30mm	30/Pk	74.95
Non-Sterile Bulk Packs			
229771	PES, 0.45 µm, 30mm	100/Pk	195.95
229774	PVDF, 0.45 µm, 30mm	100/Pk	195.95
229775	Nylon, 0.22 µm, 30mm	100/Pk	195.95
229776	Nylon, 0.45 µm, 30mm	100/Pk	195.95
229777	PTFE, 0.22 µm, 13mm	100/Pk	144.95
229778	PTFE, 0.22 µm, 30mm	100/Pk	259.95
229779	PTFE, 0.45 µm, 13mm	100/Pk	144.95
229780	PTFE, 0.45 µm, 30mm	100/Pk	249.95



PROGENE® THIN-WALL INDIVIDUAL PCR® TUBES

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Evaporation-resistant, easy-to-open-and-close lids prevent sample loss
- Packaged in resealable bags to minimize risk of contamination



CERTIFIED

RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C200-F	0.2 ml - Flat-Cap, Frosted	1000/Pk	49.95
87-C200-CF	0.2 ml - Flat-Cap, Optical	1000/Pk	49.95
87-C200-D	0.2 ml - Dome-Cap	1000/Pk	54.95
87-C200-NC	0.2 ml - No Cap	1000/Pk	62.95
87-C500-F	0.5 ml - Flat-Cap, Frosted	1000/Pk	42.95

PROGENE® 8 & 12-STRIP THIN-WALL TUBES & CAPS

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Manufactured with virgin, medical grade polypropylene
- Standard and Low-Profile wells.
- Optically-clear caps for (real time) qPCR.



Tear Strip



Rigid

CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C100-8-TB	0.1 ml 8-Strip Tubes, Low-Profile	125/Pk	80.95
87-C200-8-TB	0.2 ml 8-Strip Tubes	125/Pk	64.95
87-C200-12-TB	0.2 ml 12-Strip Tubes	80/Pk	94.95
87-C8-DOME	8-Strip Domed Caps	125/Pk	19.95
87-C8-QPCR	8-Strip Optical qPCR Caps	125/Pk	21.95
87-C12-DOME	12-Strip Domed Caps	80/Pk	22.95
87-C12-QPCR	12-Strip Optical qPCR Caps	80/Pk	22.95
87-C200-8-IND	0.2 ml Clear - Optical qPCR Caps, Tear Strip	120/Pk	79.95
87-C200-8RGD-IND	0.2 ml Clear - Optical qPCR Caps, Rigid	120/Pk	79.95
87-C200-8D-IND	0.2 ml Clear - Domed Caps, Tear Strip	120/Pk	79.95

PROGENE® THIN-WALL PCR® PLATES

- Ultra-thin walls produce efficient thermal transfers during amplifications - Broad thermal cycler compatibility - High-contrast, black alphanumeric indexing for easy well identification



CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
96-Well PCR Plates			
87-C96-NS	0.2 ml - Non-Skirted	10/Pk, 50/Cs	124.95
87-C96-LPNS	0.1 ml - Non-Skirted, Low-Profile	20/Pk	63.95
87-C96-ABI-1	0.2 ml - Half-Skirt for ABI, Hybrid	10/Pk, 50/Cs	149.95
87-C96-ABI-2	0.2 ml - Half-Skirt for ABI, Universal	10/Pk, 50/Cs	124.95
87-C96-FAST	0.1 ml - Half-Skirt for ABI Fast PCR	10/Pk, 50/Cs	149.95
87-C96-RR	0.2 ml - Half-Skirt, Raised Rim	10/Pk, 50/Cs	159.95
87-C96-SK	0.1 ml - Full-Skirt, Low-Profile	10/Pk, 50/Cs	124.95
87-C96-CT	0.2 ml - Chimney-Top	10/Pk, 50/Cs	164.95
384-Well PCR Plates			
87-C384-1N	Full Skirt, 1-Notch	10/Pk, 50/Cs	254.95
87-C384-2N	Full Skirt, 2-Notch	10/Pk, 50/Cs	254.95
Lightcycler® Plates			
87-C96-LC480-WP	96W White Plates - No Film	10/Pk	35.95

PROGENE® ROTOR-GENE® PCR® STRIP TUBES

- Ultra-thin walls produce efficient thermal transfers during amplifications
- Compatible with QIAGEN®/Corbett Rotor-Gene® Instruments

CERTIFIED

RNase, DNase & Pyrogen-Free



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-C100-4-RG	0.1 ml 4-Strip with Caps for Rotor-Gene®	250/Pk	124.95



PROGENE® AUTORADIOGRAPHY FILM

- High contrast, high sensitivity, low background, light blue film
- Especially useful with chemiluminescence, works well with 14C, 32P, 125I, and 35S


CAT. #	DESCRIPTION	PACKAGING	PRICE
39-20810	8" x 10" Autoradiography Film	100/Pk	149.95

PROGENE



PCR/qPCR OPTICAL SEALING FILM

- Polyester-based sealing film
- Ideal for PCR, qPCR (real-time) and robotics
- Choose from standard or X-Large sized film
- Temperature Range: -40°C to +120°C




CERTIFIED
RNase, DNase & Pyrogen-Free

Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PETST-100	Clear Film, Polyester	100/Pk	79.95
PS-PETSTXL-100	Clear Film, Polyester, X-Large	100/Pk	89.95

PCR SEALING FOIL

- Easy to pierce with pipet tips or robotic probes for direct sample recovery
- Ideal for standard PCR and long-term storage
- Temperature Range: -80°C to +120°C




CERTIFIED
RNase, DNase & Pyrogen-Free

Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-FOILP-100	Easy-to-Pierce Foil	100/Pk	99.95
PS-FOILP-50S	Easy-to-Pierce Foil, Sterile	50/Pk	73.95

qPCR PRESSURE-ACTIVATED SEALING FILM

- Highly chemically resistant, non-tacky, pressure-activated silicone adhesive
- For qPCR (real-time), Long-Term Storage and Protein Crystallization - Temperature Range: -80°C to +100°C




CERTIFIED
RNase, DNase & Pyrogen-Free

Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PPO-100	qPCR, Pressure-Activated Adhesive	100/Pk	245.95
PS-PPO-50S	qPCR, Pressure-Activated Adhesive, Sterile	50/Pk	165.95

EXTRA-THICK SEALING FOIL

- Extra thick, chemically-resistant aluminum minimizes shredding upon removal
- Ideal for standard PCR and long-term storage
- Temperature Range: -80°C to +120°C




CERTIFIED
RNase, DNase & Pyrogen-Free

Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-FOIL-100	Chemically Resistant Foil	100/Pk	109.95
PS-FOIL-50S	Chemically Resistant Foil, Sterile	50/Pk	71.95

SEALSTRIPS™

- Ideal for PCR, qPCR and storage
- Seals 1 or 2 x 8-well rows at a time
- Economical and more user-friendly than PCR strip caps



CERTIFIED
RNase, DNase & Pyrogen-Free

Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PS-PPSS-200	Routine PCR Strips - Seals 2 Rows of 8-Wells	200/Pk	49.95
PS-PETOSS-400	Optical qPCR Strips - Seals 1 Row of 8-Wells	400/Pk	64.95
PS-FOILSS-300	Foil Strips - Seals 1 Row of 8-Wells	300/Pk	39.95

FILM ACCESSORIES



Plateal™

CAT. #	DESCRIPTION	PACKAGING	PRICE
PSA-RRL-1	Film Roller	Each	46.95
PSA-PP-5	Sealing Applicator	5/Pk	17.95



PROGENE® SURGICAL MASK ASTM LEVEL 2

- Pleated 3 layer design
- Latex and fiberglass free
- Bacterial Filtration Efficiency (BFE) ≥ 99%
- Sub-micron particulate Efficiency @ 0.1 µm (PFE) ≥ 99%
- Resistance to penetration by synthetic blood @ 120 mmHg
- ASTM Level 2



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
WLM2002	Surgical Masks Level 2, Blue	50/Box	5.95

AURELIA® SURGICAL MASK

**PROMO
PRICE**



aurelia

CAT. #	DESCRIPTION	PACKAGING	PRICE
262-2110	Aurelia Surgical Masks Level 1, Blue	50/Box	6.95
262-2120	Aurelia Surgical Masks Level 2, Blue	50/Box	7.95
262-2130	Aurelia Surgical Masks Level 3, Blue	50/Box	8.95

PROGENE® BOUFFANT CAPS & BEARD COVER

- Full hair coverage
- Lightweight and breathable
- Made of light weight spunbond polypropylene
- Latex and lint free



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
20-4402-W	Bouffant Caps White, 21"	500/Cs	39.95
20-4403-W	Bouffant Caps White, 24"	500/Cs	40.95
20-P21-W	Pleated Bouffant Caps White, 21"	500/Cs	34.95
20-4400	Beard Covers, 18", White	1000/Cs	69.95

ULTIDENTBRAND FACE SHIELD

- The face shield is made from PET for splash protection in medical environment
- Adjustable headband for customized fit - Comfortable and breathable foam band
- Each Sold individually
- FDA Certified, OSHA compliant, meets Conformité Européenne (CE)



UltiDent
Scientific

CAT. #	DESCRIPTION	PRICE
SH-FSV	Face Shield, Anti-Fog, 34 x 22 cm	2.95

ULTIDENTBRAND DISPOSABLE ISOLATION GOWNS

- Made of spunbond polypropylene and are latex-free
- Isolation gowns are made for a comfortable fit with elastic cuffs, waist and neckties



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
77-900-L	Gowns, Yellow, Large	10/Bag, 50/Cs	69.95

ULTIDENTBRAND LATEX-FREE SHOE COVERS

- Non-conductive and non-skid with sewn seams provide increased tear strength



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
32-2132	Shoe Covers Non-Skid, Blue, Reg.	300/Bx	39.95
32-2134	Shoe Covers Non-Skid, Blue, XL	300/Bx	44.95



AURELIA® TRANSFORM™ 100 POWDER-FREE NITRILE GLOVES

aurelia



OPTIMAL STRETCHABILITY

3.2 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-9889A5	Transform™ 100 Stretch Nitrile, X-Small	Ice Blue	100/Box	13.95
262-9889A6	Transform™ 100 Stretch Nitrile, Small	Ice Blue	100/Box	13.95
262-9889A7	Transform™ 100 Stretch Nitrile, Medium	Ice Blue	100/Box	13.95
262-9889A8	Transform™ 100 Stretch Nitrile, Large	Ice Blue	100/Box	13.95
262-9889A9	Transform™ 100 Stretch Nitrile, X-Large	Ice Blue	100/Box	13.95

AURELIA® SONIC POWDER-FREE NITRILE GLOVES

aurelia



OPTIMAL STRETCHABILITY

2.2 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93775	Sonic® Stretch Nitrile, X-Small	Cobalt Blue	300/Box	49.95
262-93776	Sonic® Stretch Nitrile, Small	Cobalt Blue	300/Box	49.95
262-93777	Sonic® Stretch Nitrile, Medium	Cobalt Blue	300/Box	49.95
262-93778	Sonic® Stretch Nitrile, Large	Cobalt Blue	300/Box	49.95
262-93779	Sonic® Stretch Nitrile, X-Large	Cobalt Blue	300/Box	49.95

AURELIA® BOLD® POWDER-FREE NITRILE GLOVES

aurelia



CHEMOTHERAPY-TESTED

5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-73996	Bold® Thick Nitrile, Small	Black	100/Box	19.95
262-73997	Bold® Thick Nitrile, Medium	Black	100/Box	19.95
262-73998	Bold® Thick Nitrile, Large	Black	100/Box	19.95
262-73999	Bold® Thick Nitrile, X-Large	Black	100/Box	19.95
262-73990	Bold® Thick Nitrile, 2X-Large	Black	100/Box	19.95

AURELIA® ROBUST PLUS™ EXTENDED-CUFF 12" NITRILE POWDER-FREE GLOVES

aurelia



EXTENDED CUFF

4.5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-63886	Robust™ Plus Nitrile, Small	Blue	100/Box	19.95
262-63887	Robust™ Plus Nitrile, Medium	Blue	100/Box	19.95
262-63888	Robust™ Plus Nitrile, Large	Blue	100/Box	19.95
262-63889	Robust™ Plus Nitrile, X-Large	Blue	100/Box	19.95

AURELIA® AMAZING™ POWDER-FREE NITRILE GLOVES

aurelia



UNIQUE STRETCHY FORMULATION

2.0 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-92885	Amazing™ Nitrile, X-Small	Purple	300/Box	49.95
262-92886	Amazing™ Nitrile, Small	Purple	300/Box	49.95
262-92887	Amazing™ Nitrile, Medium	Purple	300/Box	49.95
262-92888	Amazing™ Nitrile, Large	Purple	300/Box	49.95
262-92889	Amazing™ Nitrile, X-Large	Purple	300/Box	49.95

AURELIA® BLUSH™ POWDER-FREE NITRILE GLOVES

aurelia



ULTRA SOFT AND STRETCHY

2.5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-78885	Blush™ Nitrile, X-Small	Pink	200/Box	19.95
262-78886	Blush™ Nitrile, Small	Pink	200/Box	19.95
262-78887	Blush™ Nitrile, Medium	Pink	200/Box	19.95
262-78888	Blush™ Nitrile, Large	Pink	200/Box	19.95



AURELIA® PERFORM™ POWDER-FREE NITRILE GLOVES

aurelia



OPTIMAL STRETCHABILITY

2.5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-92395	Perform™ Stretch Nitrile, X-Small	Teal	200/Box	19.95
262-92396	Perform™ Stretch Nitrile, Small	Teal	200/Box	19.95
262-92397	Perform™ Stretch Nitrile, Medium	Teal	200/Box	19.95
262-92398	Perform™ Stretch Nitrile, Large	Teal	200/Box	19.95
262-92399	Perform™ Stretch Nitrile, X-Large	Teal	200/Box	19.95

AURELIA® TRANSFORM™ 200 POWDER-FREE NITRILE GLOVES

aurelia



OPTIMAL STRETCHABILITY

3.2 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-98895	Transform™ 200 Stretch Nitrile, X-Small	Ice Blue	200/Box	25.95
262-98896	Transform™ 200 Stretch Nitrile, Small	Ice Blue	200/Box	25.95
262-98897	Transform™ 200 Stretch Nitrile, Medium	Ice Blue	200/Box	25.95
262-98898	Transform™ 200 Stretch Nitrile, Large	Ice Blue	200/Box	25.95
262-98899	Transform™ 200 Stretch Nitrile, X-Large	Ice Blue	200/Box	25.95

AURELIA® ROBUST™ POWDER-FREE NITRILE GLOVES

aurelia



OPTIMAL STRETCHABILITY

5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-93895	Robust™ Thick Nitrile, X-Small	Blue	100/Box	16.95
262-93896	Robust™ Thick Nitrile, Small	Blue	100/Box	16.95
262-93897	Robust™ Thick Nitrile, Medium	Blue	100/Box	16.95
262-93898	Robust™ Thick Nitrile, Large	Blue	100/Box	16.95
262-93899	Robust™ Thick Nitrile, X-Large	Blue	100/Box	16.95

AURELIA® VIBRANT™ POWDER-FREE NITRILE GLOVES

aurelia



UNIQUE STRETCHY FORMULATION

4.5 MIL

CAT. #	DESCRIPTION	COLOR	PACKAGING	PRICE
262-9822-XS	Vibrant™ Latex, X-Small	Beige	100/Box	13.95
262-9822-S	Vibrant™ Latex, Small	Beige	100/Box	13.95
262-9822-M	Vibrant™ Latex, Medium	Beige	100/Box	13.95
262-9822-L	Vibrant™ Latex, Large	Beige	100/Box	13.95
262-9822-XL	Vibrant™ Latex, X-Large	Beige	100/Box	13.95

ULTIDENTBRAND COVERALLS

- Made of three-layer SMS (spunbond/meltblown/spunbond) material - Offers some splash resistance while keeping users cool and comfortable



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
20-428-L	Coverall, White, Large	1/Bag, 25/Cs	99.95
20-428-XL	Coverall, White, X-Large	1/Bag, 25/Cs	99.95
20-428-XXL	Coverall, White, 2X-Large	1/Bag, 25/Cs	109.95

PROSURFACE + DISINFECTANT SPRAY AND WIPES

- ProSurface+ disinfectants, with proprietary TotalClean technology, provide powerful one-step cleaning and disinfection in a best-in-class contact time of just one minute. 72.5% Ethyl Alcohol



Medicom

CAT. #	DESCRIPTION	PACKAGING	PRICE
02-40060C	ProSurface+ Wipes, 6x6.75	Each	22.95
02-40061C	ProSurface+ Wipes, 9x12	Each	29.95
02-40066C	ProSurface+ Spray, 32oz	Each	28.95



ULTIDENTBRAND MICROSCOPE SLIDES & COVER GLASS

- All slides feature 45° clipped corners
- Our cover glass is made of first quality "Hydrolytic Class" glass
- Positive Charged Slides: Electrostatically attracts frozen tissue sections and cytology preparations



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Microscope Slides 3" X 1"			
170-7107P	Plain Slides	72/Pk	9.95
170-7107A-S	Frosted 1-End & 2-Sides, 45° Clipped Corners	72/Pk	10.95
Positive Charged Slides			
170-8105-W+	Color-Frost + Charged Slides, White	72/Pk	19.95
Cover Glass - #1			
170-C1818	18 x 18 mm	1-Oz/Pk	6.95
170-C2222	22 x 22 mm	1-Oz/Pk	6.95
170-C2240	22 x 40 mm	1-Oz/Pk	6.95
170-C2250	22 x 50 mm	1-Oz/Pk	6.95
170-C2424	24 x 24 mm	1-Oz/Pk	6.95
170-C2432	24 x 32 mm	1-Oz/Pk	6.95
170-C2440	24 x 40 mm	1-Oz/Pk	6.95
170-C2450	24 x 50 mm	1-Oz/Pk	6.95
170-C2460	24 x 60 mm	1-Oz/Pk	6.95

ULTIDENTBRAND CUVETTES

- Choose from optical polystyrene or UV grade PMMA
- Engineered for high reproducibility
- Packaged in dust-free PE foam trays of 100 cuvettes



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Polystyrene (340-800 nm)			
11-112117	1.5 ml Semi-Micro	100/Tray, 500/Cs	71.95
11-111117	4.5 ml Macro	100/Tray, 500/Cs	71.95
11-111137	4.5 ml Macro 4-Window	100/Tray, 500/Cs	71.95
Methacrylate (280-800 nm)			
11-112147	1.5 ml Semi-Micro	100/Tray, 500/Cs	109.95
11-111147	4.5 ml Macro	100/Tray, 500/Cs	109.95
11-111157	4.5 ml Macro 4-Window	100/Tray, 500/Cs	109.95

ULTIDENTBRAND DISPOSABLE BOROSILICATE CULTURE TUBES

- Chemically inert
- High heat tolerance
- Convenient dispenser packaging
- Tubes feature well-rounded bottoms and fire-polished tops free of rough edges



UltiDent
Scientific

CAT. #	DESCRIPTION	VOLUME	PACKAGING	PRICE
170-14100003B	10 x 75 mm	3 ml	250/Pk, 1000/Cs	39.95
170-14100005B	12 x 75 mm	5 ml	250/Pk, 1000/Cs	44.95
170-14100007B	13 x 100 mm	7 ml	250/Pk, 1000/Cs	49.95
170-14100010B	16 x 100 mm	14 ml	250/Pk, 1000/Cs	69.95
170-14106125B	16 x 125 mm	16 ml	100/Pk, 600/Cs	54.95
170-14106150B	16 x 150 mm	24 ml	250/Pk, 1000/Cs	104.95
170-14108150B	18 x 150 mm	28 ml	100/Pk, 600/Cs	89.95
170-14120150B	20 x 150 mm	36 ml	100/Pk, 600/Cs	99.95

PROGENE® POLYSTYRENE RESERVOIRS

- Disposable polystyrene reservoirs have indented "spouts" in each corner
- Bottoms are a sharp tapered "V" shape to minimize waste and reagent consumption



PROGENE
CERTIFIED
RNase, DNase & Pyrogen-Free

CAT. #	DESCRIPTION	PACKAGING	PRICE
Non-Sterile			
825-1-10	10 ml	100/Cs	69.95
825-1-25	25 ml	100/Cs	74.95
825-1-50	50 ml	100/Cs	79.95
825-10-100	100 ml	100/Cs	79.95
Sterile, Individually-Wrapped			
825-1-10-1S	10 ml	1/Pk, 100/Cs	119.95
825-1-25-1S	25 ml	1/Pk, 100/Cs	119.95
825-1-50-1S	50 ml	1/Pk, 100/Cs	123.95
825-10-100-1S	100 ml	1/Pk, 100/Cs	131.95
Sterile, Multi-Pack			
825-1-10-5S	10 ml	5/Pk, 200/Cs	149.95
825-1-25-5S	25 ml	5/Pk, 200/Cs	149.95
825-1-50-5S	50 ml	5/Pk, 200/Cs	159.95
825-10-100-5S	100 ml	5/Pk, 200/Cs	159.95



ULTIDENTBRAND TRANSFER PIPETTES

- Molded from LDPE
- Inert to biological fluids and most acids - Can be sealed and frozen down to -196 °C



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Non-Sterile			
170-CTB700-1	5.0 ml Graduated	500/Pk	19.95
170-CTB700-212	3.0 ml Graduated, Small Bulb	500/Pk	19.95
170-CTB700-274	5.0 ml Graduated, Blood Bank	500/Pk	19.95
170-CTB700-3	7.0 ml Graduated	500/Pk	19.95
Sterile			
170-CTB700-1S	5.0 ml Graduated	20/Pk, 400/Cs	44.95
170-CTB700-1SPP	5.0 ml Graduated	1/Pk, 400/Cs	46.95
170-CTB700-3S	7.0 ml Graduated	20/Pk, 400/Cs	44.95
170-CTB700-3SPP	7.0 ml Graduated	1/Pk, 400/Cs	46.95
170-CTB700-233S	5.0 ml Graduated, Fine Tip	20/Pk, 400/Cs	49.95

ULTIDENTBRAND CARDBOARD FREEZER BOXES

- Made of durable, moisture resistant coated cardboard
- All boxes are supplied complete with installed grids



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
21-BOX-2-81	2" Box with 81-Place Grid	12/Pk	29.95
21-BOX-2-100	2" Box with 100-Place Grid	12/Pk	29.95
21-BOX-3-81	3" Box with 81-Place Grid	12/Pk	29.95
21-BOX-3-100	3" Box with 100-Place Grid	12/Pk	29.95
21-BOX-PCR	144-Place Box for 0.2 ml PCR® Tubes	6/Pk	29.95
21-BOX-15ML	36-Place Box for 15 mL Centrifuge Tubes	6/Pk	39.95
21-BOX-50ML	16-Place Box for 50 mL Centrifuge Tubes	6/Pk	39.95

ULTIDENTBRAND STAINLESS RACKS FOR UPRIGHT FREEZERS

- Models listed below are designed for standard 5¼" x 5¼" x 2" standard boxes
- Drawer racks let you reach boxes at the back of the freezer - No need to remove the entire rack

UltiDent
Scientific



CAT. #	DESCRIPTION	L x H x W (MM)	L x H	PRICE
Standard Racks for 2" boxes				
2525-UF-332	9-Place	415x168x140	3 x 3	139.95
2525-UF-342	12-Place	415x225x140	3 x 4	175.95
2525-UF-352	15-Place	415x280x140	3 x 5	211.95
2525-UF-432	12-Place	550x168x140	4 x 3	182.95
2525-UF-442	16-Place	550x225x140	4 x 4	231.95
2525-UF-452	20-Place	550x280x140	4 x 5	287.95
2525-UF-462	24-Place	550x321x140	4 x 6	344.95
2525-UF-542	20-Place	679x225x140	5 x 4	325.95
2525-UF-552	25-Place	679x280x140	5 x 5	400.95
Drawer Racks for 2" Boxes				
2525-UFD-332	9-Place	419x181x140	3 x 3	186.95
2525-UFD-342	12-Place	419x240x140	3 x 4	247.95
2525-UFD-442	16-Place	559x240x140	4 x 4	324.95
2525-UFD-472	28-Place	559x419x140	4 x 7	571.95
2525-UFD-552	25-Place	680x300x140	5 x 5	571.95
Drawer Racks for 15 & 50 ml Tubes				
2525-UFD-LT15-1	15 ml x 80 Tubes, 1 Drawer	559x126x140	---	235.95
2525-UFD-LT15-2	15 ml x 160 Tubes, 2 Drawers	559x252x140	---	365.95
2525-UFD-LT50-1	50 ml x 39 Tubes, 1 Drawer	559x126x140	---	235.95
2525-UFD-LT50-2	50 ml x 78 Tubes, 2 Drawers	559x252x140	---	357.95
Racks for 100-Place Hinged-Top Freezer Boxes				
2525-UFHT-23	6-Place	318x176x149	2x3	101.95
2525-UFHT-24	8-Place	318x235x149	2x4	130.95
2525-UFHT-25	10-Place	318x293x149	2x5	142.95
2525-UFHT-33	9-Place	476x176x149	3x3	166.95
2525-UFHT-34	12-Place	476x235x149	3x4	194.95
2525-UFHT-35	15-Place	476x293x149	3x5	227.95
2525-UFHT-44	16-Place	634x235x149	4x4	255.95



ULTIDENTBRAND WIRELESS TUBE RACKS

- Racks sink in a water bath - Excellent temperature and chemical resistance
- Fully Autoclavable - Molded from reinforced nylon for strength and durability



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
Half Size			
11-456445	36 Place for 13 mm Tubes	Each	20.95
11-456545	36 Place for 16 mm Tubes	Each	20.95
11-456645	20 Place for 15 ml Centrifuge Tubes	Each	20.95
11-456845	9 Place for 50 ml Centrifuge Tubes	Each	20.95
Full Size			
11-456440	72 Place for 13 mm Tubes	Each	29.95
11-456540	72 Place for 16 mm Tubes	Each	29.95
11-456640	40 Place for 15 ml Centrifuge Tubes	Each	29.95
11-456840	24 Place for 50 ml Centrifuge Tubes	Each	29.95

Available in 7 colors: B (Blue), G (Green), M (Magenta),
O (Orange), Y (Yellow), W (White)

PROGENE® 96-PLACE PCR RACKS

- Lid clicks into place


PROGENE



CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W96-RK	96-Place PCR Racks, Assorted	5/Pk	44.95

ULTIDENTBRAND 80-WELL MICROTUBE RACKS

- Fraction collector plate style.
- Versatile & efficient storage.
- Lots of color options (contact us for more)
- Fits Tubes up to 2.0 ml.



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
24-R80-AS	80-Place Workstations, Assorted	5/Pk	29.95

ULTIDENTBRAND LAB DRYING RACK & LAB MATS



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-120791	Drying Rack 17.4 x 26 x 5.5 in	Each	304.95
33-120747	Silicone Lab Mat 30.5 cm x 37.9 cm; Blue	Each	79.95
33-120748	Silicone Lab Mat 30.5 cm x 37.9 cm; Purple	Each	79.95

PROGENE® BENCHTOP COOLERS



PROGENE

CAT. #	DESCRIPTION	PACKAGING	PRICE
87-W24-MCT	24-Place MCT Cooler with Lid	Each	89.95
87-W96-RG	64-Place RotorGene® Cooler with Lid	2/Pk	109.95
87-W96-PCR	96-Place PCR Cooler with Lid	2/Pk	99.95

ULTIDENTBRAND FREEZER BOXES



UltiDent
Scientific

CAT. #	DESCRIPTION	PACKAGING	PRICE
33-10601-A	50-Place, Hinged lid, Assorted	5/Pk	39.95
33-10601-N	50-Place, Hinged lid, Natural	5/Pk	39.95
33-23255-N	81-Place, Natural	5/Pk	39.95
33-23255-A	81-Place, Assorted	5/Pk	39.95
33-1064	100-Place, Hinged lid, Natural	5/Pk	39.95
33-1065	100-Place, Hinged lid, Assorted	5/Pk	39.95



TRANSFERPETTE® S PIPETTES

**BUY 3, GET
1 FREE**

CAT. #	VOLUME	PRICE
Single-Channel Pipettes		
34-705868	0,1-1 µl	441.95
34-705869	0,1-2,5 µl	441.95
34-705870	0,5-10 µl	441.95
34-705872	2-20 µl	441.95
34-705873	5-50 µl	441.95
34-705874	10-100 µl	441.95
34-705878	20-200 µl	441.95
34-705880	100-1000 µl	441.95
34-705882	500-5000 µl	441.95

Multi-8 Pipettes

34-705900	Multi-8, 0,5-10 µl	1,129.95
34-705906	Multi-8, 5-50 µl	1,129.95
34-705908	Multi-8, 10-100 µl	1,129.95
34-705910	Multi-8, 20-200 µl	1,129.95
34-705912	Multi-8, 30-300 µl	1,129.95

Multi-12 Pipettes

34-705920	Multi-12, 0,5-10 µl	1,396.95
34-705926	Multi-12, 5-50 µl	1,396.95
34-705928	Multi-12, 10-100 µl	1,396.95
34-705930	Multi-12, 20-200 µl	1,396.95
34-705932	Multi-12, 30-300 µl	1,396.95

**BUY 2, GET
1 FREE**

ACCU-JET® PRO PIPETTE CONTROLLER

- Available in four attractive colors
- Variable motor speed - Fills a 25 mL pipette at max. speed in only 3 seconds

**BUY 3, GET
1 FREE**

CAT. #	COLOR	PRICE
34-26350	Anthracyte	579.95
34-26351	Berry	579.95
34-26352	Petrol	579.95
34-26353	Amethyst	579.95

**Buy 500 PD-TIPs™ syringe tips,
get a HandyStep S® Manual pipette for FREE!**

HANDYSTEP® S REPEATING PIPETTE

- HandyStep® S manual: Can pipette up to 90 pre-set volumes ranging from 2 µl to 50 ml

CAT. #	DESCRIPTION	PACKAGING	PRICE
34-2705110	HandyStep® S Manual	Each	609.95
PD-Tip™ II Non-sterile			
34-705700	0.1 ml	100/Pk	175.95
34-705702	0.5 ml	100/Pk	175.95
34-705704	1.0 ml	100/Pk	175.95
34-705706	1.25 ml	100/Pk	175.95
34-705708	2.5 ml	100/Pk	175.95
34-705710	5 ml	100/Pk	175.95
34-705712	10 ml	100/Pk	175.95
34-705714	12.5 ml	100/Pk	175.95
34-705716	25 ml (adapter included)	50/Pk	119.95
34-705718	50 ml (adapter included)	25/Pk	59.95

PD-Tip™ II Sterile

34-705730	0.1ml	1/Pk, 100/Pk	239.95
34-705732	0.5ml	1/Pk, 100/Pk	239.95
34-705734	1.0ml	1/Pk, 100/Pk	239.95
34-705736	1.25ml	1/Pk, 100/Pk	239.95
34-705738	205ml	1/Pk, 100/Pk	239.95
34-705740	5ml	1/Pk, 100/Pk	239.95
34-705742	10ml	1/Pk, 100/Pk	239.95
34-705744	12.5ml	1/Pk, 100/Pk	239.95

DISPENSETTE® S BOTTLETOP DISPENSERS

**BUY 3, GET
1 FREE**

- Fast priming - Ideal for acids, bases and saline solutions
- Recirculation valve system enhances safety
- Models for organic solvents and trace analysis also available

CAT. #	DESCRIPTION	INCREMENTS (ml)	PRICE
Dispensette® S, Digital Easy Calibration™			
34-4600311	0.1-1	0.005	849.95
34-4600321	0.2-2	0.01	849.95
34-4600331	0.5-5	0.02	849.95
34-4600341	1.10	0.05	849.95
34-4600351	2.5-25	0.1	1,129.95
34-4600361	5-50	0.2	1,145.95



ULTIDENTBRAND PIPETTE CONTROLLER

**PROMO
PRICE**

- Long-life li-ion battery with quick charge - Up to 20 hours continuous use
- UV resistant housing - Bright LEDs indicate three operation modes: high, low, and gravity
- **Includes:** Charging stand and wall-mountable/magnetic holder, Nozzle adapter for use with small pipettes, Li-Ion rechargeable battery with AC adapter, Three replacement 0.45 µm filters and Five nosecones

UltiDent
Scientific



CAT. #	DESCRIPTION	PRICE
33-3000	Pipette Controller	299.95

PROGENE® MINI CENTRIFUGE

- Ideal for quick spin downs - Simply close the lid and the rotor instantly reaches 6,000 rpm
- **Includes:** Large capacity 8-place microtube rotor, 16-place x 0.2 ml PCR tube and strip rotor and adapters



PROGENE

CAT. #	DESCRIPTION	PRICE
48-1008-C	Mini Centrifuge with Clear Lid	299.95
48-1008-B	Mini Centrifuge with Blue Lid	299.95
48-1008-G	Mini Centrifuge with Green Lid	299.95
48-1008-P	Mini Centrifuge with Purple Lid	299.95
48-1008-R	Mini Centrifuge with Red Lid	299.95

MAGNETIC INDUCTION STIRRER

- This stirrer has no moving parts - 100% wear-free and maintenance free
- Magnetic coil technology works on the inductive principle with AC as its driving force
- **Multiple stirring modes:** Clockwise, counterclockwise and oscillate - Adjustable time: 30 seconds to 60 minutes or continuous - Speed: 10 to 2,000 rpm - Capacity: 3L total liquid volume, diameter of beaker not to exceed 165 mm
- **Package includes:** 2 Stir bars 35 mm x 12 mm, silicone anti-skid plate cover and low voltage, double insulated power adapter



Heathrow Scientific
design | function | innovation

CAT. #	DESCRIPTION	PRICE
33-120584	Grey with Purple Accents	520.95
33-120585	Grey with Blue Accents	520.95

MINI MAGNETIC STIRRER

- Capacity: Up to 3L
- Speed: 200 to 2200 rpm (Infinitely adjustable)
- Quiet operation with low vibration - Small but powerful



Heathrow Scientific
design | function | innovation

CAT. #	DESCRIPTION	PRICE
33-120155	Mini Magnetic Stirrer Grey/Purple	155.95

DIGITAL ORBITAL SHAKER

- Maximum load of 4.5 kg with an overload protection cut out
- **Timer range:** 0 to 49 hours, 1 second increments or continuous
- **Platform dimensions:** 27.9 W x 27.9 D cm
- Very low start speed for sensitive sample agitation with high top end speed
- **Speed range:** 20 to 300 rpm 1 rpm increments - Orbit diameter: 19.0 mm (0.75")



Heathrow Scientific
design | function | innovation

CAT. #	DESCRIPTION	PRICE
33-120460	Digital Orbital Shaker	1,343.95



SPROUT® PLUS MINI-CENTRIFUGE

- **Package Includes:** 6 x 1.5/2.0 ml rotor, 16 x 0.2 ml PCR rot, 0.2 ml adapters, 0.5 ml adapters and storage case



CAT. #	DESCRIPTION	PRICE
33-120610	Sprout®, 6,000 rpm - Max RCF: 2,000 Xg, Green	316.95
33-120611	Sprout®, 6,000 rpm - Max RCF: 2,000 Xg, Blue	316.95

MINI VORTEX

- Touch start operation enables mixing with one hand
- Non-slip feet ensures the unit stays in place while mixing
- Speed: 4,500 rpm - Orbit: 6 mm - Capacity: Up to 50 ml tube



CAT. #	DESCRIPTION	PRICE
33-120567	Mini Vortex Gray/Purple	156.95
33-120598	Mini Vortex Gray/Blue	156.95

VORTEXER - VARIABLE SPEED VORTEXING UNIT

- Provides holding slots for six different sizes of tubes and a standard assay plate, while the center pad area can be used for quick hand-held mixing of a vast range of tubes and vessels.
- Variable speed rotary dial provides pre-set speed control (1000, 2000, 3000 rpm) during vortex mixing



CAT. #	DESCRIPTION	PRICE
33-120209	Vortexer® Variable Speed	366.95

MAGFUGE® CENTRIFUGE & MAGNETIC STIRRER

- High-speed Centrifuge and Magnetic Stirrer in ONE
- Digital display, one button programming and tool free rotor exchange
- **MagFuge® Centrifuge:** RCF: 16-9,783 x g (12,500 RPM), Capacity: 12 x 1.5/2.0 ml tubes or 6 x 5 ml tubes
- **MagFuge® Magnetic Stirrer:** Speed Range: 50-2500 rpm, Capacity: 3L maximum
- **Package Includes:** 5 ml Rotor, 1.5/2.0 ml Rotor, Stirrer Rotor; 12 x 0.2 ml adapters, 12 x 0.5 ml adapters, 2 Silicone mats for stirrer lid & 2 Rare-Earth PTFE coated stir bars



CAT. #	DESCRIPTION	PRICE
33-120581	MagFuge® Centrifuge & Magnetic Stirrer	1,435.95

GUSTO® HIGH-SPEED MINI CENTRIFUGE

- ideal for protocols calling for fast spins with a higher RCF requirement. Ideal for PCR®, microfilter cell separation, clinical & diagnostic applications, and HPLC protocols.



CAT. #	DESCRIPTION	PRICE
33-10050	Gusto® High Speed Mini Centrifuge	1,369.95

ULTIDENTBRAND SILICONE LAB MATS

- Protect countertops from hot items up to 200°C (392°F)
- Chemical resistance: chemically inert and does not react with most chemicals
- Mats measure 35 cm x 60 cm, 2 mm thick - Autoclavable at 121° for 15 minutes



CAT. #	DESCRIPTION	PRICE
33-120506	Side One: Solid Yellow - Side Two: Blue design on grey background	71.95
33-120507	Side One: Solid Purple - Side Two: Grey design on purple background	71.95
33-120599	Side One: Solid Blue - Side Two: Solid white with blue rim	71.95

Prices Subject to change without notice.

